

PKanalix 2023R1 single page user guide



1. PKanalix documentation

Version 2023

This documentation is for PKanalix.
©Lixoft

PKanalix performs analysis on PK data set including:

- The **non compartmental analysis** (NCA) – computation of the NCA parameters using the calculation of the λ_z – slope of the terminal elimination phase.
- The **bioequivalence analysis** (BE) – comparison of the NCA parameters of several drug products (usually two, called 'test' and 'reference') using the average bioequivalence.
- The **compartmental analysis** (CA) – estimation of model parameters representing the PK as the dynamics in compartments for each individual. It does not include population analysis performed in **Monolix**.

What else?

- A clear **user-interface** with an **intuitive workflow** to efficiently run the NCA, BE and CA analysis.
- Easily accessible **PK models library** and **auto-initialization** method to improve the convergence of the optimization of CA parameters.
- **Integrated bioequivalence** module to simplify the analysis
- **Automatically generated results and plots** to give an immediate feedback.
- Interconnection with MonolixSuite applications to **export projects to Monolix** for the population analysis.

PKanalix tasks


PKanalix uses the dataset format common for all MonolixSuite applications, see [here](#) for more details. It allows to move your project between applications, for instance export a CA project to Monolix and perform a population analysis with one “click” of a button.

Non Compartmental Analysis	Bioequivalence	Compartmental Analysis
This task consists in defining rules for the calculation of the λ_z (slope of the terminal elimination phase) to compute all the NCA parameters. This definition can be either global via general rules or manual on each individual – with the interactive plots the user selects or removes points in the λ_z calculation.	The average NCA parameters obtained for different groups (e.g a test and a reference formulation) can be compared using the Bioequivalence task. Linear model definition contains one or several fixed effects selected in an integrated module. It allows to obtain a confidence interval compared to the predefined BE limits and automatically displayed in intuitive tables and plots.	This task estimates parameters of a pharmacokinetic model for each individual. The model can be custom or based on one of the user-friendly libraries. Automatic initialization method improves the convergence of parameters for each individual.

All the NCA, BE and/or CA **outputs are automatically displayed** in sortable tables in the Results tab. Moreover, they are **exported in the result folder** in a R-compatible format. **Interactive plots** give an immediate feedback and help to better interpret the results.


The usage of PKanalix is available not only via the user interface, but also via R with a dedicated R-package (detailed [here](#)). All the actions performed in the interface have their equivalent R-functions. It is particularly convenient for reproducibility purpose or batch jobs.

The results of the NCA calculations and bioequivalence calculations have been compared on an extensive number of datasets to the results of WinNonLin and to published results obtained with SAS. All results were identical. See the poster below for more details.



2. Defining a dataset

To start a new PKanalix project, you need to define a dataset by loading a csv file in the Data tab.



The data set format expected in the Data tab is the same as for the entire MonolixSuite, to allow smooth transitions between applications. The columns available in this format and example datasets are detailed [on this page](#). Briefly:

- Each line corresponds to one individual and one time point
- Each line can include a single measurement (also called observation), or a dose amount (or both a measurement and a dose amount)
- Dosing information should be indicated for each individual, even if it is identical for all.

Your dataset may not be originally in this format, and you may want to add information on dose amounts, limits of quantification, units, or filter part of the dataset. To do so, you should proceed in this order:

- 1) **Formatting**: If needed, format your data first by loading the dataset in the Data Formatting tab. Briefly, it allows to:
 - to deal with several header lines
 - merge several observation columns into one
 - add censoring information based on tags in the observation column
 - add treatment information manually or from external sources
 - add more columns based on another file

or


- 1) If the data is already in the right format, **load it directly in the Data tab** (otherwise use the formatted dataset created by data formatting). *If the dataset does not follow a formatting rule, the dataset will not be accepted, but errors will guide you to find what is missing and could be added by data formatting.*

- 2) **Labeling**: Label the columns not recognized automatically to indicate their type and click on ACCEPT.
- 3) **Units**: If you want to work with units, indicate the units in the data and the ones you prefer to use inside the Data tab (if relevant)
- 4) **Filters**: If needed, filter your dataset to use only part of it in the Filters tab
- 5) **Explore**: The interpreted dataset is displayed in Data, and Plots and covariate statistics are generated.

To use a dataset from a Monolix project, or to use simulations from a Simulx project, you can directly **import** or **export** the Monolix/Simulx project which will automatically define the dataset in the data tab.

Labeling


The column type suggested automatically by PKanalix based on the headers in the data can be customized in the preferences. By clicking on Settings>Preferences, the following windows pops up.



In the DATA frame, you can add or remove preferences for each column.

To remove a preference, double-click on the preference you would like to remove. A confirmation window will be proposed.

To add a preference, click on the header type you consider, add a name in the header name and click on "ADD HEADER" as on the following figure.




Notice that all the preferences are shared between Monolix, Datxplorer, and PKanalix.

Resulting plots and tables to explore the data

Once the dataset is accepted:

- **Plots** are automatically generated based on the interpreted dataset to help you proceed with a first data exploration before running any task.



- The **interpreted dataset** appears in Data tab, which incorporates all changes after formatting, setting units, and filtering.

LINE NUMBER	ID	TIME	CONC	AMOUNT	ADMINISTRATION	ROUTE
1	1	0		100000	1	IV
2	1	0.083	3826.27			IV
3	1	0.25	3101.72			IV
4	1	0.5	3955.83			IV
5	1	1	3024.63			IV
6	1	2	3441.82			IV
7	1	8	3869.37			IV
8	1	12	3127.52			IV
9	1	24	2259.83			IV
10	1	32	2092.63			IV
11	1					

- **Covariate Statistics** appear in a section of the data tab.

	MIN	Q1	MEDIAN	Q3	MAX	SD
AGE	21	22	25	28	39	5
HT	176	179	181	183	191	4.26

	ref	test	RT	TR
FORM	1.8	1.8		
SEQ	1.8	1.8		

2.1. Data format for NCA and CA analysis

In PKanalix, a dataset should be loaded in the **Data tab** to create a project. Once the dataset accepted, it is possible to specify **units** and **filter** the dataset, so units and filtering information should not be included in the file loaded in the Data tab.

The data set format used in PKanalix is the same as for the entire MonolixSuite, to allow smooth transitions between applications. In this format:

- Each line corresponds to one individual and one time point.
- Each line can include a single measurement (also called observation), or a dose amount (or both a measurement and a dose amount).
- **Dosing information should be indicated for each individual in a specific column, even if it is the same treatment for all individuals.**
- Headers are free but there can be only one header line.
- Different types of information (dose, observation, covariate, etc) are recorded in different columns, which must be tagged with a column type (see below).

If your dataset is not in this format, in most cases, it is possible to format it in a few steps in the **data formatting** tab, to incorporate the missing information.

- **Overview of most common column types**
- **Example datasets**
 - **Plasma concentration data**
 - **extravascular**
 - **IV infusion**
 - **IV bolus**

- Steady-state data
- BLQ data
- Urine data
- Occasions ("Sort" variables)
- Covariates ("Carry" variables)
- Other useful column-types
- Description of all possible column types

Overview of most common column-types

A data set typically contains only the following columns: **ID**, **TIME**, **OBSERVATION**, **AMOUNT**. The main rules for interpreting the dataset are:

- Cells that do not contain information (e.g AMOUNT column on a line recording a measurement) should have a dot.
- Headers are free but there can be only one header line and the columns must be assigned to one of the available column-types. The full list of column-types is available at the end of this page and a detailed description is given on the [dataset documentation website](#).

The most common column types in addition to ID, TIME, OBSERVATION and AMOUNT are:

- For **IV infusion** data, an additional **INFUSION RATE** or **INFUSION DURATION** is necessary.
- For **steady-state** data, **STEADY-STATE** and **INTERDOSE INTERVAL** column are added.
- If a dose and a measurement occur at the same time, they can be encoded on the same line or on different lines.
- **Sort and carry variables** can be defined using the **OCCASION**, and **CATEGORICAL COVARIATE** and **CONTINUOUS COVARIATE** column-types.
- **BLQ data** are defined using the **CENSORING** and **LIMIT** column-types.

In addition to the typical cases presented above, a few additional column-types may also be convenient:

- **ignoring a line:** with the **IGNORED OBSERVATION** (ignores the measurement only) or **IGNORED LINE** (ignores the entire line, including regressor and dose information). However it is more convenient to filter lines of your dataset without modifying it by using **filters** (available once your dataset is accepted).
- **working with several types of observations:** If several observation types are present in the dataset (for example parent and metabolite), all measurements should still appear in the same **OBSERVATION** column, and another column should be used to distinguish the observation types. If this is not the case in your dataset, **data formatting** enables to merge several observation columns. **To run NCA for several observations** at the same time for the same id, tag the column listing the observation types as **OCCASION**. **To run CA for several observations** with a model including several outputs (for example PK and PD), tag the column listing the observation types as **OBSERVATION ID**. In this case, only one observation type will be available for NCA. It can be selected in the "NCA" section to perform the calculations.
- **working with several types of administrations:** a single data set can contain different types of administrations, for instance IV bolus and extravascular, distinguished using the **ADMINISTRATION ID** column-type. The setting "administration type" in "Tasks>Run" can be chosen separately for each administration id, and the appropriate parameter calculations will be performed.

Example datasets

Below we show how to encode the dataset for most typical situations.

Plasma concentration data

Extravascular

For extravascular administration, the mandatory column-types are **ID** (individual identifiers as integers or strings), **OBSERVATION** (measured concentrations), **AMOUNT** (dose amount administered) and **TIME** (time of dose or measurement).

ID	TIME	AMOUNT	CONC		ID	TIME	AMOUNT	CONC
1	0	100	.		1	0	100	0
1	0	.	0		1	1	.	8.1
1	1	.	8.1	or	1	2	.	7.5
1	2	.	7.5		1	3	.	6.9
1	3	.	6.9		1	4	.	6.5
1	4	.	6.5		1	5	.	5.8
1	5	.	5.8					


To distinguish the extravascular from the IV bolus case, in "Tasks>Run" the administration type must be set to "extravascular".

If no measurement is recorded at the time of the dose, a concentration of zero is added for single dose data, the minimum concentration observed during the dose interval for steady-state data.


Example:

- demo **project_extravascular.pkx**


This data set records the drug concentration measured after single oral administration of 150 mg of drug in 20 patients. For each individual, the first line records the dose (in the "Amount" column tagged as **AMOUNT** column-type) while the following lines record the measured concentrations (in the "Conc" column tagged as **OBSERVATION**). Cells of the "Amount" column on measurement lines contain a dot, and respectively for the concentration column. The column containing the times of measurements or doses is tagged as **TIME** column-type and the subject identifiers, which we will use as sort variable, are tagged as **ID**. Check the **OCCASION** section if more sort variables are needed. After accepting the dataset, the data is automatically assigned as "Plasma".



In the “Tasks/Run” tab, the user must indicate that this is extravascular data. In the “Check lambda_z”, on linear scale for the y-axis, measurements originally present in the data are shown with full circles. Added data points, such as a zero concentration at the dose time, are represented with empty circles. Points included in the λ_z calculation are highlighted in blue.



After running the NCA analysis, PK parameters relevant to extravascular administration are displayed in the “Results” tab.



IV infusion

Intravenous infusions are indicated in the data set via the presence of an **INFUSION RATE** or **INFUSION DURATION** column-type, in addition to the **ID** (individual identifiers as integers or strings), **OBSERVATION** (measured concentrations), **AMOUNT** (dose amount administered) and **TIME** (time of dose or measurement). The infusion duration (or rate) can be identical or different between individuals.

ID	TIME	AMOUNT	CONC	TINF		ID	TIME	AMOUNT	CONC	TINF
1	0	100	.	0.5		1	0	100	0	0.5
1	0	.	0	.		1	1	.	8.1	.
1	1	.	8.1	.	or	1	2	.	7.5	.
1	2	.	7.5	.		1	3	.	6.9	.
1	3	.	6.9	.		1	4	.	6.5	.
1	4	.	6.5	.		1	5	.	5.8	.
1	5	.	5.8	.						


In “Tasks>Run” the administration type must be set to “intravenous”.

If no measurement is recorded at the time of the dose, a concentration of zero is added for single dose data, the minimum concentration observed during the dose interval for steady-state data.


Example:

- demo **project_ivinfusion.pkx**:

In this example, the patients receive an iv infusion over 3 hours. The infusion duration is recorded in the column called “TINF” in this example, and tagged as **INFUSION DURATION**.



In the “Tasks/Run” tab, the user must indicate that this is intravenous data.



IV bolus

For IV bolus administration, the mandatory column-types are **ID** (individual identifiers as integers or strings), **OBSERVATION** (measured concentrations), **AMOUNT** (dose amount administered) and **TIME** (time of dose or measurement).

ID	TIME	AMOUNT	CONC		ID	TIME	AMOUNT	CONC
1	0	100	.		1	0	100	10.4
1	0	.	10.4		1	1	.	8.1
1	1	.	8.1	or	1	2	.	7.5
1	2	.	7.5		1	3	.	6.9
1	3	.	6.9		1	4	.	6.5
1	4	.	6.5		1	5	.	5.8
1	5	.	5.8					


To distinguish the IV bolus from the extravascular case, in “Tasks>Run” the administration type must be set to “intravenous”.

If no measurement is recorded at the time of the dose, the concentration of at time zero is extrapolated using a log-linear regression of the first two data points, or is taken to be the first observed measurement if the regression yields a slope ≥ 0 . See the [calculation details](#) for more information.

Example:

- demo **project_ivbolus.pkx**:

In this data set, 25 individuals have received an iv bolus and their plasma concentration have been recorded over 12 hours. For each individual (indicated in the column “Subj” tagged as **ID** column-type), we record the dose amount in a column “Dose”, tagged as **AMOUNT** column-type. The measured concentrations are tagged as **OBSERVATION** and the times as **TIME**. Check the **OCCASION** section if more sort variables are needed in addition to **ID**. After accepting the dataset, the data is automatically assigned as “plasma”.



In the “Tasks/Run” tab, the user must indicate that this is intravenous data. In the “Check lambda_z”, measurements originally present in the data are shown with full circles. Added data points, such as the C₀ at the dose time, are represented with empty circles. Points included in the λ_z

After running the NCA estimation task, **steady-state specific parameters** are displayed in the "Results" tab.

BLQ data

- "0" indicates that the value in the **OBSERVATION** column is the measurement.
- "1" indicates that the observation is BLQ.

ID	TIME	AMT	CONC	BLQ
27	6	.	4.76	0
27	8	.	3.07	0
27	12	.	1.98	0
27	18	.	2.57	0
27	24	.	1.8	1
27	30	.	1.8	1
28	0	50	LOQ value	mark BLQ data
28	0.5	.	2	0
28	2	.	4.13	0
28	4	.	2.97	0
28	6	.	5.17	0

- zero
- the LOQ value
- the LOQ value divided by 2
- or considered as missing


A **LIMIT** column can be added to record the other limit of the interval (in general zero). This value will not be used by PKanalix but can facilitate the transition from an NCA/CA analysis PKanalix to a population model with Monolix.

Examples:


- In this dataset, the measurements of two different studies (indicated in the STUDY column, tagged as **CATEGORICAL COVARIATE** in order to be carried over) are recorded. For the study 101, the LOQ is 1.8 ug/mL, while it is 1 ug/mL for study 102. The BLQ data are marked with a "1" in the BLQ column, which is tagged as **CENSORING**. The LOQ values are indicated for each BLQ in the CONC column of measurements, tagged as **OBSERVATION**.

LINE NUMBER	ID	TIME	AMT	CONC	OBSERVATION	BLQ	CENSORING	CATEGORICAL COVARIATE
1	1	0	50					
2	1	0.5		3.08				
3	1	2		5.92				
4	1	4		4.7				
5	1	6		4.14				
6	1	8		4.45				
7	1	12		2.44				
8	1	18		2				
9	1	24		1.5				
10	1	30		1.5				
11	1	30		1.5				
12	2	0	50					

In the “Task>NCA>Run” tab, the user can choose how to handle the BLQ data. For the BLQ data before and after the T_{max}, the BLQ data can be considered as missing (as if this data set row would not exist), or replaced by zero (default before T_{max}), the LOQ value or the LOQ value divided by 2 (default after T_{max}). In the “Check lambda_z” tab, the BLQ data are shown in green and displayed according to the replacement value.



For the CA analysis, the replacement value for all BLQ can be chosen in the settings of the “Run” tab (default is Missing). In the “Check init.” tab, the BLQ are again displayed in green, at the LOQ value (irrespective of the chosen method for the calculations).



Urine data

To work with urine data, it is necessary to record the time and amount administered, the volume of urine collected for each time interval, the start and end time of the intervals and the drug concentration in a urine sample of each interval. The time intervals must be continuous (no gaps allowed).

In PKanalix, the start and end times of the intervals are recorded in a single column, tagged as **TIME** column-type. In this way, the end time of an interval automatically acts as start time for the next interval. The concentrations are recorded in the **OBSERVATION** column. The volume column must be tagged as **REGRESSOR** column type. This general column-type of MonolixSuite data sets allows to easily transition to the other applications of the Suite. As several **REGRESSOR** columns are allowed, the user can select which **REGRESSOR** column should be used as volume. The concentration and volume measured for the interval [t₁,t₂] are noted on the t₂ line. The volume value on the dose line is meaningless, but it cannot be a dot. We thus recommend to set it to zero.

A typical urine data set has the following structure. A dose of 150 ng has been administered at time 0. The first sampling interval spans from the dose at time 0 to 4h post-dose. During this time, 410 mL of urine have been collected. In this sample, the drug concentration is 112 ng/mL. The second interval spans from 4h to 8h, the collected urine volume is 280 mL and its concentration is 92 ng/mL. The third interval is marked on the figure: 390mL of uring have been collected from 8h to 12h.

ID	TIME	VOL	CONC	AMOUNT
1	0	0	.	150
1	4	410	112	.
1	8	280	92	.
1	12	390	41	.
1	18	450	29	.
1	24	415	12	.
1	48	1730	2.6	.
1	72	2030	0.1	.


The given data is used to calculate the intervals midpoints, and the excretion rates for each interval. This information is then used to calculate the λ_z and calculate urine-specific parameters. In “Tasks/Check lambda_z”, we display the midpoints and excretion rates. However, in the “Plots>Data viewer”, we display the measured concentrations at the end time of the interval.

Example:

- demo **project_urine.pkx**: urine PK dataset

In this urine PK data set, we record the consecutive time intervals in the “TIME” column tagged as **TIME**. The collected volumes and measured concentration are in the columns “VOL” and “CONC”, respectively tagged as **REGRESSOR** and **OBSERVATION**. Note that the volume and concentration are indicated on the line of the interval end time. The volume on the first line (start time of the first interval, as well as dose line) is set to zero as it must be a double. This value will not be used in the calculations. Once the dataset is accepted, the observation type must be set to “urine” and the regressor column corresponding to the volume defined.

In “Tasks>Check lambda_z”, the excretion rate are plotted on the midpoints time for each individual. The choice of the lambda_z works as usual.



Once the NCA task has run, **urine-specific PK parameters** are displayed in the “Results” tab.

Occasions (“Sort” variables)

The main sort level are the individuals indicated in the **ID** column. **Additional sort levels can be encoded using one or several OCCASION column(s).** **OCCASION** columns contain integer values that permit to distinguish different time periods for a given individual. The time values can restart at zero or continue when switching from one occasion to the next one. The variables differing between periods, such as the treatment for a crossover study, are tagged as **CATEGORICAL** or **CONTINUOUS COVARIATES** (see below). The NCA and CA calculations will be performed on each ID-OCCASION combination. Each occasion is considered independent of the other occasions (i.e a washout is applied between each occasion).

Note: occasions columns encoding the sort variables as integers can easily be added to an existing data set using Excel or R.

With R, the “OCC” column can be added to an existing “data” data frame with a column “TREAT” using `data$OCC <- ifelse(data$TREAT=="ref", 1, 2)`.

With Excel, assuming the sort variable is encoded in the column E with values “ref” and “test”, type `=IF(E2="ref",1,2)` to generate the first value of the “OCC” column and then propagate to the entire column:

A	B	C	D	E	F	G	H
ID	TIME	Y	AMT	TREAT	OCC		
1	0	.	100	ref			
3	1	3.17902	.	ref			
4	1.5	3.85033	.	ref			
5	1	2	4.46786	.	ref		
6	1	2.5	5.22826	.	ref		
7	1	3	4.86716	.	ref		
8	1	4	5.0127	.	ref		
9	1	5	5.55581	.	ref		
10	1	6	3.94428	.	ref		
11	1	8	3.1716	.	ref		
12	1	10	3.19092	.	ref		
13	1	12	3.82381	.	ref		
14	1	18	0.8774	.	ref		
15	1	24	0.39638	.	ref		
16	1	0	100	test			
17	1	1	6.57151	.	test		
18	1	1	7.66164	.	test		


Examples:

- demo **project_occasions1.ppk**: crossover study with two treatments

The subject column is tagged as **ID**, the treatment column as **CATEGORICAL COVARIATE** and an additional column encoding the two periods with integers “1” and “2” as **OCCASION** column.

LINE NUMBER	ID	TIME	OCC	CONC	AMOUNT	TREAT
1	1	0	1		100	ref
2	1	1	1	1.7902		ref
3	1	1.5	1	3.8503		ref
4	1	2	1	4.4678		ref
5	1	2.5	1	5.2283		ref
6	1	3	1	4.8671		ref
7	1	4	1	5.0127		ref
8	1	5	1	5.5558		ref
9	1	6	1	3.9443		ref
10	1	8	1	3.1716		ref
11	1	10	1	3.1909		ref
12	1	12	1	3.8238		ref
13	1	18	1	0.8774		ref
14	1	24	1	0.3964		ref
15	1	0	2		100	test
16	1	1	2	6.5715		test
17	1	1	2	7.6616		test


In the “Check lambda_z” (for the NCA) and the “Check init.” (for CA), each occasion of each individual is displayed. The syntax “1#2” indicates individual 1, occasion 2, according to the values defined in the **ID** and **OCCASION** columns.



In the “Individual estimates” output tables, the first columns indicate the **ID** and **OCCASION** (reusing the data set headers). The covariates are indicated at the end of the table. Note that it is possible to sort the table by any column, including, ID, OCCASION and COVARIATES.

ID	OCC	AUCINF_0_100	AUCINF_0_100_p	AUCINF_0_100_s	AUCINF_0_100_t	AUCINF_0_100_e	AUCINF_0_100_f	AUCINF_0_100_g	AUCINF_0_100_h	AUCINF_0_100_i	AUCINF_0_100_j	AUCINF_0_100_k	AUCINF_0_100_l	AUCINF_0_100_m	AUCINF_0_100_n	AUCINF_0_100_o	AUCINF_0_100_p	AUCINF_0_100_q	AUCINF_0_100_r	AUCINF_0_100_s	AUCINF_0_100_t	AUCINF_0_100_u	AUCINF_0_100_v	AUCINF_0_100_w	AUCINF_0_100_x	AUCINF_0_100_y	AUCINF_0_100_z
1	1	0.888	0.888	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88
1	2	0.888	0.888	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88	88.88	88.88	8.88	8.88
2	1	1.101	1.101	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01
2	2	1.101	1.101	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01	110.10	110.10	11.01	11.01
3	1	1.229	1.229	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29
3	2	1.229	1.229	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29	122.90	122.90	12.29	12.29
4	1	1.464	1.464	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64
4	2	1.464	1.464	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64	146.40	146.40	14.64	14.64
5	1	1.689	1.689	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89
5	2	1.689	1.689	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89	168.90	168.90	16.89	16.89
6	1	0.823	0.823	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23
6	2	0.823	0.823	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23
7	1	0.895	0.895	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95
7	2	0.895	0.895	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95	89.50	89.50	8.95	8.95
8	1	0.823	0.823	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23
8	2	0.823	0.823	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23	82.30	82.30	8.23	8.23

The OCCASION values are available in the plots for stratification, in addition to possible **CATEGORICAL** or **CONTINUOUS COVARIATES** (here “TREAT”).




- demo **project_occasions2.pkx**: study with two treatments and with/without food

In this example, we have three sorting variables: ID, TREAT and FOOD. The TREAT and FOOD columns are duplicated: once with strings to be used as **CATEGORICAL COVARIATE** (TREAT and FOOD) and once with integers to be used as **OCCASION** (occT and occF).

LINE NUMBER	ID	TIME	CONC	AMOUNT	TREAT	FOOD	occT	occF
2	1	0	-	100	ref	yes	1	1
3	1	1	3.01805	-	ref	yes	1	1
4	1	1.5	4.154	-	ref	yes	1	1
5	1	2	5.15407	-	ref	yes	1	1
6	1	2.5	5.47307	-	ref	yes	1	1
7	1	3	5.03303	-	ref	yes	1	1
8	1	4	5.23809	-	ref	yes	1	1
9	1	5	5.47307	-	ref	yes	1	1

In the individual parameters tables and plots, three levels are visible. In the “Individual parameters vs covariates”, we can plot Cmax versus FOOD, split by TREAT for instance (Cmax versus TREAT split by FOOD is also possible).



Covariates (“Carry” variables)

Individual information that need to be carried over to output tables and plots must be tagged as **CATEGORICAL or **CONTINUOUS COVARIATES**.** Categorical covariates define variables with a few categories, such as treatment or sex, and are encoded as strings. Continuous covariates define variables on a continuous scale, such as weight or age, and are encoded as numbers. Covariates will not automatically be used as “Sort” variables. A dedicated **OCCASION** column is necessary (see above).

Covariates will automatically appear in the output tables. Plots of estimated NCA and/or CA parameters versus covariate values will also be generated. In addition, covariates can be used to **stratify** (split, color or filter) any plot. Statistics about the covariates distributions are available in table format in “Results > Cov. stat.” and in graphical format in “Plots > Covariate viewer”.

Note: It is preferable to avoid spaces and special characters (stars, etc) in the strings for the categories of the categorical covariates. Underscores are allowed.

Example:

- demo **project_covariates.pkx**

In this data set, “SEX” is tagged as **CATEGORICAL COVARIATE** and “WEIGHT” as **CONTINUOUS COVARIATE**.

The figure consists of two side-by-side screenshots of the RStudio interface, illustrating the 'Covariate statistics' and 'Covariate visual' windows.

Left Screenshot: Covariate statistics

The 'Covariate statistics' window displays summary statistics for continuous covariates. The table shows the distribution of 'weight' across different categories.

	MIN	Q1	MEDIAN	Q3	MAX	SD
weight	62.34	68.87	69.35	73.44	83.27	4.85

Below the table, there is a section for 'Categorical covariates' showing the distribution across modalities. The table shows the distribution of 'sex'.

	F	M
sex	23	27

Right Screenshot: Covariate visual

The 'Covariate visual' window displays visualizations for the same covariates. It includes a histogram for 'weight', a boxplot for 'sex', and a bar chart for 'sex'.

The histogram for 'weight' shows the frequency distribution of the variable. The boxplot for 'sex' shows the distribution of the variable across the two categories (F and M). The bar chart for 'sex' shows the count of individuals in each category.

Project Settings Export Help

Home Data Tasks **Results** Plots

SUMMARY INDIV ESTIM.

RMA

CON STAT


Individual estimates

Non compartmental analysis individual

	model	Ns_points	lambda_0	Req	Res_adjusted	Open	TO	Tlag	Tlead	Tmean	VLF_40s	VLF_pert	WEIGHT	SEX	ID
	4	0.38	0.939	1.55	0	0	24	1.5	5.7	5.67	70.65	F			1
13	3	1	1	0.65	0	0	24	2	2.29	2.29	67.52	F			2
	4	0.354	0.924	1.15	0	0	24	3	17.37	17.42	69.31	M			3
13	3	0.999	0.992	2.13	0	0	24	1.5	7.52	7.51	66.47	M			4
	3	0.995	0.992	1.8	0	0	24	2	5.76	5.75	69.13	F			5
13	3	0.999	0.996	1.15	0	0	24	5	8.47	8.48	71.19	F			6
	5	0.36	0.947	2.4	0	0	24	3	9.18	9.16	71.81	F			7
13	3	0.973	0.969	3.52	0	0	24	2.5	12.79	12.66	76.76	M			8
	4	0.323	0.913	2.37	0	0	24	5	3.8	3.56	64.49	F			9
13	3	0.995	0.992	1.64	0	0	24	5	3.93	3.95	64.74	F			10
	4	0.399	0.994	2.12	0	0	24	2.5	11.42	11.39	74.94	M			11

The screenshot shows the RStudio interface with the 'Plots' tab selected. The left sidebar contains a 'Data viewer' and a 'Correlation viewer' with a table of correlations. The main plot area displays four plots: two scatter plots (AUC vs. log2FC and AUC vs. log2FC) and two box plots (AUC vs. SEX and AUC vs. SEX). The right sidebar shows the 'Display' tab with a list of individual parameters and a 'Display' button.

Individual parameters	Correlation
AUC	0.0000
AUCNP_2_gred	0.0000
AUCNP_3_gred	0.0000
AUCNP_4_gred	0.0000
AUCNP_5_gred	0.0000
AUCNP_6_gred	0.0000
AUCNP_7_gred	0.0000
AUCNP_8_gred	0.0000
AUCNP_9_gred	0.0000
AUCNP_10_gred	0.0000
AUCNP_11_gred	0.0000
AUCNP_12_gred	0.0000
AUCNP_13_gred	0.0000
AUCNP_14_gred	0.0000
AUCNP_15_gred	0.0000
AUCNP_16_gred	0.0000
AUCNP_17_gred	0.0000
AUCNP_18_gred	0.0000
AUCNP_19_gred	0.0000
AUCNP_20_gred	0.0000
AUCNP_21_gred	0.0000
AUCNP_22_gred	0.0000
AUCNP_23_gred	0.0000
AUCNP_24_gred	0.0000
AUCNP_25_gred	0.0000
AUCNP_26_gred	0.0000
AUCNP_27_gred	0.0000
AUCNP_28_gred	0.0000
AUCNP_29_gred	0.0000
AUCNP_30_gred	0.0000
AUCNP_31_gred	0.0000
AUCNP_32_gred	0.0000
AUCNP_33_gred	0.0000
AUCNP_34_gred	0.0000
AUCNP_35_gred	0.0000
AUCNP_36_gred	0.0000
AUCNP_37_gred	0.0000
AUCNP_38_gred	0.0000
AUCNP_39_gred	0.0000
AUCNP_40_gred	0.0000
AUCNP_41_gred	0.0000
AUCNP_42_gred	0.0000
AUCNP_43_gred	0.0000
AUCNP_44_gred	0.0000
AUCNP_45_gred	0.0000
AUCNP_46_gred	0.0000
AUCNP_47_gred	0.0000
AUCNP_48_gred	0.0000
AUCNP_49_gred	0.0000
AUCNP_50_gred	0.0000
AUCNP_51_gred	0.0000
AUCNP_52_gred	0.0000
AUCNP_53_gred	0.0000
AUCNP_54_gred	0.0000
AUCNP_55_gred	0.0000
AUCNP_56_gred	0.0000
AUCNP_57_gred	0.0000
AUCNP_58_gred	0.0000
AUCNP_59_gred	0.0000
AUCNP_60_gred	0.0000
AUCNP_61_gred	0.0000
AUCNP_62_gred	0.0000
AUCNP_63_gred	0.0000
AUCNP_64_gred	0.0000
AUCNP_65_gred	0.0000
AUCNP_66_gred	0.0000
AUCNP_67_gred	0.0000
AUCNP_68_gred	0.0000
AUCNP_69_gred	0.0000
AUCNP_70_gred	0.0000
AUCNP_71_gred	0.0000
AUCNP_72_gred	0.0000
AUCNP_73_gred	0.0000
AUCNP_74_gred	0.0000
AUCNP_75_gred	0.0000
AUCNP_76_gred	0.0000
AUCNP_77_gred	0.0000
AUCNP_78_gred	0.0000
AUCNP_79_gred	0.0000
AUCNP_80_gred	0.0000
AUCNP_81_gred	0.0000
AUCNP_82_gred	0.0000
AUCNP_83_gred	0.0000
AUCNP_84_gred	0.0000
AUCNP_85_gred	0.0000
AUCNP_86_gred	0.0000
AUCNP_87_gred	0.0000
AUCNP_88_gred	0.0000
AUCNP_89_gred	0.0000
AUCNP_90_gred	0.0000
AUCNP_91_gred	0.0000
AUCNP_92_gred	0.0000
AUCNP_93_gred	0.0000
AUCNP_94_gred	0.0000
AUCNP_95_gred	0.0000
AUCNP_96_gred	0.0000
AUCNP_97_gred	0.0000
AUCNP_98_gred	0.0000
AUCNP_99_gred	0.0000
AUCNP_100_gred	0.0000



Description of all possible column types

Column-types used for all types of lines:

- **ID (mandatory)**: identifier of the individual
- **OCCASION (formerly OCC)**: identifier (index) of the occasion
- **TIME (mandatory)**: time of the dose or observation record
- **DATE/DAT1/DAT2/DAT3**: date of the dose or observation record, to be used in combination with the TIME column
- **EVENT ID (formerly EVID)**: identifier to indicate if the line is a dose-line or a response-line
- **IGNORED OBSERVATION (formerly MDV)**: identifier to ignore the OBSERVATION information of that line
- **IGNORED LINE (from 2019 version)**: identifier to ignore all the informations of that line
- **CONTINUOUS COVARIATE (formerly COV)**: continuous covariates (which can take values on a continuous scale)
- **CATEGORICAL COVARIATE (formerly CAT)**: categorical covariate (which can only take a finite number of values)
- **REGRESSOR (formerly X)**: defines a regression variable, i.e a variable that can be used in the structural model (used e.g for time-varying covariates)
- **IGNORE**: ignores the information of that column for all lines


Column-types used for response-lines:

- **OBSERVATION (mandatory, formerly Y)**: records the measurement/observation for continuous, count, categorical or time-to-event data
- **OBSERVATION ID (formerly YTYPE)**: identifier for the observation type (to distinguish different types of observations, e.g PK and PD)
- **CENSORING (formerly CENS)**: marks censored data, below the lower limit or above the upper limit of quantification
- **LIMIT**: upper or lower boundary for the censoring interval in case of CENSORING column

Column-types used for dose-lines:

- **AMOUNT (mandatory, formerly AMT)**: dose amount
- **ADMINISTRATION ID (formerly ADM)**: identifier for the type of dose (given via different routes for instance)
- **INFUSION RATE (formerly RATE)**: rate of the dose administration (used in particular for infusions)
- **INFUSION DURATION (formerly TINF)**: duration of the dose administration (used in particular for infusions)
- **ADDITIONAL DOSES (formerly ADDL)**: number of doses to add in addition to the defined dose, at intervals INTERDOSE INTERVAL
- **INTERDOSE INTERVAL (formerly II)**: interdose interval for doses added using ADDITIONAL DOSES or STEADY-STATE column types
- **STEADY STATE (formerly SS)**: marks that steady-state has been achieved, and will add a predefined number of doses before the actual dose, at interval INTERDOSE INTERVAL, in order to achieve steady-state

2.2. Data formatting

 PKanalix versions prior to 2023R1 do not include the data formatting module. Instead we provide an [Excel macro](#) to adapt the format of your dataset for PKanalix.

The **dataset format that is used in PKanalix** is the same as for the entire MonolixSuite, to allow smooth transitions between applications. In this format, some rules have to be fulfilled, for example:

- Each line corresponds to one individual and one time point.
- Each line can include a single measurement (also called observation), or a dose amount (or both a measurement and a dose amount).
- Dose amount should be indicated for each individual dose in a column AMOUNT, even if it is identical for all.
- Headers are free but there can be only one header line.

If your dataset is not in this format, in most cases, it is possible to format it in a few steps in the **data formatting** tab, to incorporate the missing information.

In this case, the original dataset should be loaded in the "Format data" box, or directly in the "Data Formatting" tab, instead of the "Data" tab. In the data formatting module, you will be guided to build a dataset in the MonolixSuite format, starting from the loaded csv file. The resulting formatted dataset is then loaded in the Data tab as if you loaded an already-formatted dataset in "Data" directly. Then [as for defining any dataset](#), you can tag columns, accept the dataset, and once accepted, **units** can be specified in the Data tab and the **Filters** tab can be used to select only parts of this dataset for analysis. Note that units and filters are neither information to be included in the data file, nor part of the data formatting process.




The original datafile is NOT modified by PKanalix. Formatting operations are saved in the PKanalix project and applied to the data when the project is loaded. The formatted dataset is not saved by default, but it can be exported by the user as a CSV file.

Jump to:

1. **Data formatting workflow**
2. **Dataset initialization** (mandatory step)
 - **Selecting header lines or lines to exclude** to merge header lines or exclude a line
 - **Tagging mandatory columns** such as ID and TIME
 - **Initialization example** Demo project *CreateOcc_AdmlIdbyCategory.pkx*
3. **Creating occasions from a SORT column** to distinguish different sets of measurements within each subject, (eg formulations). Demo project *CreateOcc_AdmlIdbyCategory.pkx*
4. **Selecting an observation type** (required to add a treatment)
5. **Merging observations from several columns**
 - **As observation ids** to map them to several outputs in CA. Demo project *merge_obsID_ParentMetabolite.pkx*
 - **As occasions** to analyze them simultaneously in NCA. Demo project *merge_occ_ParentMetabolite.pkx*
 - **Option "Duplicate information from undefined columns"**
6. **Specifying censoring from censoring tags** eg "BLQ" instead of a number in an observation column. Demo project *DoseAndLOQ_manual.pkx*
7. **Adding doses in the dataset** Demo project *DoseAndLOQ_manual.pkx*
8. **Reading censoring limits or dosing information from the dataset** "by category" or "from data". Demo projects *DoseAndLOQ_byCategory.pkx* and *DoseAndLOQ_fromData.pkx*
9. **Creating occasions from dosing intervals** to analyze separately the measurements following different doses. Demo project *doseIntervals_as_Occ.mlxtran*
10. **Handling urine data** to merge start and end times in a single column. Demo project *Urine_LOQinObs.pkx*
11. **Adding new columns from an external file**, eg new covariates, or individual parameters estimated in a previous analysis. Demo *warfarin_PKPDseq_project.mlxtran*
12. **Exporting the formatted dataset**

1. Data formatting workflow





When opening a new project, two Browse buttons appear. The first one, under "Data file", can be used to load a dataset already in a MonolixSuite-standard format, while the second one, under "Format data", allows to load a dataset to format in the Data formatting module.

After loading a dataset to format, data formatting operations can be specified in several subtabs: Initialization, Observations, Treatments and Additional columns.

- Initialization is mandatory and must be filled before using the other subtabs.
- Observations is required to enable the Treatments tab.


After Initialization has been validated by clicking on "Next", a button "Preview" is available from any subtab to view in the Data tab the formatted dataset based on the formatting operations currently specified.



 [Back to List of Formatting Steps](#)

2. Dataset initialization

The first tab in Data formatting is named Initialization. This is where the user can select header lines or lines to exclude (in the blue area on the screenshot below) or tag columns (in the yellow area).



LINE NUMBER	ID	TIME_h	CONC_mg_L	SORT	AGE	WT
1	Description: crossover study with two formulations					
2	ID	TIME	CONC		AGE	WT
3		h	mg/L			
4	1	0	BLQ	ref	25	66
5	1	1	NA	ref	25	66

Selecting header lines or lines to exclude

These settings should contain line numbers for lines that should be either handled as column headers or that should be excluded.

- **Header lines:** one or several lines containing column header information. By default, the first line of the dataset is selected as header. If several lines are selected, they are merged by data formatting into a single line, concatenating the cells in each column.
- **Lines to exclude (optional):** lines that should be excluded from the formatted dataset by data formatting.

Tagging mandatory columns

Only the columns corresponding to the following tabs must be tagged in Initialization, while all the other columns should keep the default UNDEFINED tag:

- **ID (mandatory):** subject identifiers
- **TIME (mandatory):** the single time column
- **SORT (optional):** one or several columns containing SORT variables can be tagged as SORT. Occasions based on these columns will be created in the formatted dataset as described in [Section 3](#).
- **START, END and VOLUME (mandatory in case of urine data):** these column tags replace the TIME tag in case of urine data, if the urine collection time intervals are encoded in the dataset with two time columns for the start and end times of the intervals. In that case there should also be a column with the urine volume in each interval. See [Section 10](#) for more details.

UNDEFINED	UNDEFINED
UNDEFINED	TIME
ID	TIME
TIME	TIME
SORT	TIME
Urine headers	
START	h
END	0
VOLUME	
1	1

Initialization example

- demo **CreateOcc_AdmlIdbyCategory.pkx** (the screenshot below focuses on the formatting initialization and excludes other elements present in the demo):

In this demo the first line of the dataset is excluded because it contains a description of the study. The second line contains column headers while the third line contains column units. Since the MonolixSuite-standard format allows only a single header line, lines 2 and 3 are merged together in the formatted dataset.

CreateOcc_AdmlIdbyCategory.pkx [demo] - PKanalix - 2023R1

Project Settings Help

Data Tasks Results Plots

Data formatting Initialization Observations Treatments Additional columns

RESTORE DATA PREVIEW

Data file: C:/Users/pauli/rooft/pkanalix/2023R1/demos/0_data_formatting/data/two_formulations_data.csv

Filters

Covariates statistics

Lines settings

Header lines: 2 3

Lines to exclude: 1

LINE NUMBER	ID	TIME	CONC	FORM	AGE	WT	SEQ
1	Description: crossover study with two formulations						
2	ID	TIME	CONC	FORM	AGE	WT	SEQ
3		h	mg/L				
4	1	0	BLQ	ref	25	66	RT
5	1	1	NA	ref	25	66	RT

Records per page: 10 25 50 100 500 866 867

Page 1 of 18

Showing 1 to 50 of 867 entries

PREVIEW

CreateOcc_AdmlIdbyCategory.pkx [demo] - PKanalix - 2023R1

Project Settings Help

Data Tasks Results Plots

Data formatting Data

Filters

Covariates statistics

Data information

Urine Plasma

Units

TIME column x 1 = hour

CONC column x 1 = milligram / Liter

AMT column x 1 = milligram

LINE NUMBER	ID	TIME	OBSERVATION	CATEGORICAL COVARIATE	CONTINUOUS COVARIATE	CONTINUOUS COVARIATE	CATEGORICAL COVARIATE	OCCASION
2	1	0	BLQ	ref	25	66	RT	1
3	1	1	NA	ref	25	66	RT	1
4	1	2	0.34	ref	25	66	RT	1
5	1	3	0.98	ref	25	66	RT	1

Records per page: 10 25 50 100 500 864

Page 1 of 18

Showing 1 to 50 of 864 entries - 18 Individuals


BACK RESTORE DATA ACCEPT

3. Creating occasions from a SORT column

A SORT variable can be used to distinguish different sets of measurements (usually concentrations) within each subject, that should be analyzed separately by PKanalix (for example: different formulations given to each individual at different periods of time, or multiple doses where concentration profiles are available to be analyzed following several doses).

In PKanalix, these different sets of measurements must be distinguished as OCCASIONS (or periods of time), via the OCCASION column-type. However, a column tagged as OCCASION can only contain integers with occasion indexes. Thus, if a column with a SORT variable contains strings, its format must be adapted by Data formatting, in the following way:


- the user must tag the column as SORT in the Initialization subtab of Data formatting,
- the user validates the Initialization with "Next", then clicks on "Preview" (after optionally defining other data formatting operations),
- the formatted data is shown in Data: the column tagged as SORT is automatically duplicated. The original column is automatically tagged as CATEGORICAL COVARIATE in Data, while the duplicated column, which has the same name appended with "_OCC", is tagged as OCCASION. This column contains occasion indexes instead of strings.



Example:

- demo **CreateOcc_AdmlbyCategory.ppkx** (the screenshot below focuses on the formatting of occasions and excludes other elements present in the demo):


The image below shows lines 25 to 29 from the dataset from the CreateOcc_AdmlbyCategory.ppkx demo, where covariate columns have been removed to simplify the example. This dataset contains two sets of concentration measurements for each individual, corresponding to two different drug formulations administered on different periods. The sets of concentrations are distinguished with the FORM column, which contains "ref" and "test" categories (reference/test formulations). The column is tagged as SORT in Data formatting Initialization. After clicking on "Preview", we can see in the Data tab that a new column named FORM_OCC has been created with occasion indexes for each individual: for subject 1, FORM_OCC=1 corresponds to the reference formulation because it appears first in the dataset, and FORM_OCC=2 corresponds to the test formulation because it appears second in the dataset.



[Back to List of Formatting Steps](#)

4. Selecting an observation type

The second subtab in Data formatting allows to select one or several observation types. An observation type corresponds to a column of the dataset, that contains a type of measurements (usually drug concentrations, but it can also be PD measurements for example). Only columns that have not been tagged as ID, TIME or SORT are available as observation type.



This action is optional and can have several purposes:

- If doses must be added by Data formatting (see [Section 7](#)), specifying the column containing observations is mandatory, to avoid duplicating observations on new dose lines.
- If several observation types exist in different columns (for example: concentrations for different analytes, or measurements for PK and PD), they

must be specified in Data formatting to be merged into a single observation column (see [Section 5](#)).


- In the MonolixSuite-standard format, the column containing observations can only contain numbers, and no string except “.” for a missing observation. Thus if this column contains strings in the original dataset, it must be adapted by Data formatting, with two different cases:
 - if the strings are tags for censored observations (usually BLQ: below the limit of quantification), they can be specified in Data formatting to adapt the encoding of the censored observations (see [Section 6](#)),
 - any other string in the column is automatically replaced by “.” by Data formatting.

 [Back to List of Formatting Steps](#)

5. Merging observations from several columns

The MonolixSuite-standard format allows a single column containing all observations (such as concentrations or PD measurements). Thus if a dataset contains several observation types in different columns (for example: concentrations for different analytes, or measurements for PK and PD), they must be specified in Data formatting to be merged into a single observation column.

In that case, different settings can be chosen in the area marked in orange in the screenshot below:



- The user must choose between distinguishing observation types with observation ids or occasions.
- The user can unselect the option “Duplicate information from undefined columns”.

As observation ids

After selecting the “Distinguish observation types with: observation ids” option and clicking “Preview,” the columns for different observation types are combined into a single column called “OBS.” Each row of the dataset is duplicated for each observation type, with one value per observation type. Additionally, an “OBSID” column is created, with the name of the observation type corresponding to the measurement on each row.

This option is recommended for joint modeling of observation types, such as CA in PKanalix or population modeling in Monolix. It is important to note that NCA cannot be performed on two different observation ids simultaneously, so it is necessary to choose one observation id for the analysis.

Example:

- demo **merge_obsID_ParentMetabolite.pkx** (the screenshot below focuses on the formatting of observations and excludes other elements present in the demo):

This demo involves two columns that contain drug parent and metabolite concentrations. When merging both observation types with observation ids, a new column called OBSID is generated with categories labeled as “PARENT” and “METABOLITE.”

The top screenshot shows the 'Data file' tab with the file path `...\\pkanalix\\2023R1\\demos\\0_data_formatting\\data\\parent_metabolite_data.csv`. The 'Lines settings' tab shows 'Header lines' set to 1 and 'Lines to exclude' set to an empty field. The data table below has columns: LINE NUMBER, ID, TIME, PARENT, METABOLITE, DOSE, and SEX. The bottom screenshot shows the 'Data information' tab with 'METABOLITE' set to 'Plasma' and 'PARENT' set to 'Plasma'. The 'Units' tab shows 'TIME' as 'hour', 'CONC' as 'milligram', and 'AMT' as 'milligram'. The data table below has columns: LINE NUMBER, ID, TIME, DOSE, SEX, OBS, and OBSID. A settings dialog box on the right shows 'Distinguish observation types with:' with 'observation ids' selected and 'Duplicate information from undefined columns' checked. A 'PREVIEW' button is also visible.

As occasions

After selecting the "Distinguish observation types with: occasions" option and clicking "Preview," the columns for different observation types are combined into a single column called "OBS." Each row of the dataset is duplicated for each observation type, with one value per observation type. Additionally, two columns are created: an "OBSID_OCC" column with the index of the observation type corresponding to the measurement on each row, and an "OBSID_COV" with the name of the observation type.

This option is recommended for NCA, which can be run on different occasions for each individual. However, joint modeling of the observation types with CA or population modeling with Monolix cannot be performed with this option.

Example:

- demo **merge_occ_ParentMetabolite.pkx**:

This demo involves two columns that contain drug parent and metabolite concentrations. When merging both observation types with occasions, two new columns called OBSID_OCC and OBSID_COV are generated with OBSID_OCC=1 corresponding to OBSID_COV="PARENT" and OBSID_OCC=2 corresponding to OBSID_COV="METABOLITE."

merge_obsID_ParentMetabolite.plx [demo] - PKanalix - 2023R1

Project Settings Help

Data Tasks Plots

Data formatting Initialization Observations Treatments Additional columns RESTORE DATA PREVIEW

Data file

...\\pkanalix\2023R1\demos\0_data_formatting\data\parent_metabolite_data.csv

BROWSE...

Lines settings

Header lines: 1

Lines to exclude:

LINE NUMBER	ID	TIME	PARENT	METABOLITE	DOSE	SEX
1	ID	TIME	PARENT	METABOLITE	DOSE	SEX
2	1	1	13.4806	1.00136	1000	M
3	1	2	24.7161	2.66854	1000	M
4	1	4	34.6926	9.95048	1000	M

Records per page: 10 25 50 100 500 721

<< Page 1 of 15 >>

Showing 1 to 50 of 721 entries

merge_occ_ParentMetabolite.plx [demo] - PKanalix - 2023R1

Project Settings Help

Data Tasks Plots

Data formatting Data Filters Covariates statistics

Data file

Data information

Units

TIME column x 1 = hour

CONC column x 1 = milligram / milliliter

AMT column x 1 = milligram /

Disp

LINE NUMBER	ID	TIME	DOSE	SEX	OBS	OBSID_OCC	OBSID_COV
2	1	1	1000	M	13.4806	1	PARENT
3	1	1	1000	M	1.00136	2	METABOLITE
4	1	2	1000	M	24.7161	1	PARENT
5	1	2	1000	M	2.66854	2	METABOLITE

Records per page: 10 25 50 100 500 867 900 1440

<< Page 1 of 29 >>

Showing 1 to 50 of 1,440 entries - 40 Individuals

BACK RESTORE DATA ACCEPT

Settings

Distinguish observation types with:

☐ observation ids (recommended for joint CA and population modeling)

☒ occasions (recommended for NCA)

☒ Duplicate information from undefined columns

PREVIEW

Duplicate information from undefined columns

When merging two observation columns into a single column, all other columns will see their lines duplicated. The data formatting will know how to treat columns which have been tagged in the Initialization tab, but not the other columns (header "UNDEFINED") which are not used for data formatting. A checkbox enables to decide if the information from these columns should be duplicated on the new lines, or if "." should be used instead. The default option is to duplicate information, because in general, the undefined columns correspond to covariates with one value per individual, so this value is the same for the two lines that correspond to the same id.

merge_obsID_ParentMetabolite.plx [demo] - PKanalix - 2023R1

Project Settings Help

Data Tasks Plots

Data formatting Initialization Observations Treatments Additional columns RESTORE DATA PREVIEW

Data file

...\\pkanalix\2023R1\demos\0_data_formatting\data\parent_metabolite_data.csv

BROWSE...

Lines settings

Header lines: 1

Lines to exclude:

DOSE and SEX are undefined

LINE NUMBER	ID	TIME	PARENT	METABOLITE	DOSE	SEX
1	ID	TIME	PARENT	METABOLITE	DOSE	SEX
2	1	1	13.4806	1.00136	1000	M
3	1	2	24.7161	2.66854	1000	M
4	1	4	34.6926	9.95048	1000	M
5	1	6	38.251	17.8245	1000	M

Records per page: 10 25 50 100 500 721

<< Page 1 of 15 >>

Showing 1 to 50 of 721 entries

merge_obsID_ParentMetabolite.plx [demo] - PKanalix - 2023R1

Project Settings Help

Data Tasks Plots

Data formatting Data Filters Covariates statistics

Data file

Data information

Units

TIME column x 1 = hour

CONC column x 1 = milligram / milliliter

AMT column x 1 = milligram /

Add covariates

None

LINE NUMBER	ID	TIME	DOSE	SEX	AMT	OBS	CENS	LIMIT	OBSID
1	1	0	1000	M	1000				
2	1	1	1000	M		13.4806	notCensored		PARENT
3	1	1	1000	M		1.00136	notCensored		METABOLITE
4	1	1	1000	M		24.7161	notCensored		PARENT
5	1	2	1000	M		2.66854	notCensored		METABOLITE

Records per page: 10 25 50 100 500 1440

<< Page 1 of 29 >>

Showing 1 to 50 of 1,440 entries - 40 Individuals

Settings

Distinguish observation types with:

☒ observation ids (recommended for joint CA and population modeling)

☐ occasions (recommended for NCA)

☒ Duplicate information from undefined columns

PREVIEW

It is rare that you need to uncheck this box. An example where you should not duplicate the information is if you already have a column Amount in

the MonolixSuite format, so with a dose amount only at the dosing time, and "." everywhere else. If you do not want to specify amount again in data formatting, and simply want to merge observation columns as observation ids, you should not duplicate the lines of the Amount column which is undefined. Indeed, the dose amounts have been administered only once.

[Back to List of Formatting Steps](#)

6. Specifying censoring from censoring tags

In the MonolixSuite-standard format, censored observations are encoded with a 1 or -1 flag in a column tagged as CENSORING in the Data tab, while exact observations have a 0 flag in that column. In addition, on rows for censored observations, the LOQ is indicated in the observation column: it is the LLOQ (lower limit of quantification) if CENSORING=1 or the ULOQ (upper limit of quantification) if CENSORING=-1. Finally, to specify a censoring interval, an additional column tagged as LIMIT in the Data tab must exist in the dataset, with the other censoring bound.

The Data Formatting module can take as input a dataset with censoring tags directly in the observation column, and adapt the dataset format as described above. After selecting one or several observation types in the Observations subtab (see [Section 4](#)), all strings found in the corresponding columns are displayed in the "Censoring tags" on the right of the observation types. If at least one string is found, the user can then define some censoring associated with an observation type and with one or several censoring tags with the button "Add censoring". 3 types of censoring can be defined:

- **LLOQ:** this corresponds to left-censoring, where the censored observation is below a lower limit of quantification (LLOQ), that must be specified by the user. In that case Data Formatting replaces the censoring tags in the observation column by the LLOQ, and creates a new CENS column tagged as CENSORING in the Data tab, with 1 on rows that had censoring tags before formatting, and 0 on other rows.
- **ULOQ:** this corresponds to right-censoring, where the censored observation is above an upper limit of quantification (ULOQ), that must be specified by the user. Here Data Formatting replaces the censoring tags in the observation column by the ULOQ, and creates a new CENS column tagged as CENSORING in the Data tab, with -1 on rows that had censoring tags before formatting, and 0 on other rows.
- **Interval:** this is for interval-censoring, where the user must specify two bound of a censoring interval, to which the censored observation belong. Data Formatting replaces the censoring tags in the observation column by the upper bound of the interval, and creates two new columns: a CENS column tagged as CENSORING in the Data tab, with 1 on rows that had censoring tags before formatting, and 0 on other rows, and a LIMIT column with the lower bound of the censoring interval on rows that had censoring tags before formatting, and "." on other rows.

For each type of censoring, available options to define the limits are:

- **“Manual”**: limits are defined manually, by entering the limit values for all censored observations.
- **“By category”**: limits are defined manually for different categories read from the dataset.
- **“From data”**: limits are directly read from the dataset.

The options “by category” and “from data” are described in detail in [Section 8](#).

Example:

- demo **DoseAndLOQ_manual.pkx** (the screenshot below focuses on the formatting of censored observations and excludes other elements present in the demo):

In this demo there are two censoring tags in the CONC column: BLQ1 (from Study 1) and BLQ2 (from Study 2), that correspond to different LLOQs. An interval censoring is defined for each censoring tag, with manual limits, where LLOQ=0.06 for BLQ1 and LLOQ=0.1 for BLQ2, and the lower limit of the censoring interval being 0 in both cases.

The screenshot shows the PKAnalix software interface. The top window is titled 'DoseAndLOQ_manual.pkx [demo] - PKAnalix - 2023R1'. It has tabs for 'Data', 'Tasks', 'Plots', 'Initialization', 'Observations', 'Treatments', and 'Additional columns'. The 'Data' tab is active, showing 'Data file' and 'Lines settings'. The 'Lines settings' section has 'Header lines: 1-5' and 'Lines to exclude:'. Below this is a table with columns: LINE NUMBER, ID, TIME_h, CONC_mg_L, and STUDY. The table shows data for lines 1 to 5, with line 4 having a value of 0.34 in the CONC_mg_L column and line 5 having a value of 0.34. The bottom window is also titled 'DoseAndLOQ_manual.pkx [demo] - PKAnalix - 2023R1'. It has tabs for 'Data', 'Tasks', 'Plots', 'Data information', 'Units', 'TIME', 'CONC', 'AMT', 'CATEGORICAL COVARIATE', 'CENSORING', and 'LIMIT'. The 'Data' tab is active, showing a table with columns: LINE NUMBER, ID, TIME_h, CONC_mg_L, STUDY, CENS, and LIMIT. The table shows data for lines 2 to 6, with line 2 having a value of 0.06 in the CONC_mg_L column and line 3 having a value of 0.06. The 'CENS' column contains values like 1, 1, 0, 0, 0. The 'LIMIT' column contains values like 0, 0, ., ., ..

[Back to List of Formatting Steps](#)

7. Adding doses in the dataset

Datasets in MonolixSuite-standard format should contain all information on doses, as dose lines. An AMOUNT column records the amount of the administrated doses on dose-lines, with “.” on response-lines. In case of infusion, an INFUSION DURATION or INFUSION RATE column records the infusion duration or rate. If there are several types of administration, an ADMINISTRATION ID column can distinguish the different types of doses with integers.

If doses are missing from a dataset, the Data Formatting module can be used to add dose lines and dose-related columns: after initializing the dataset, the user can specify one or several treatments in the Treatments subtab. The following operations are then performed by Data Formatting:

- a new dose line is inserted in the dataset for each defined dose, with the dataset sorted by subject and times. On such a dose line, the values from the next line are duplicated for all columns, except for the observation column in which “.” is used for the dose line.
- A new column AMT is created with “.” on all lines except on dose lines, on which dose amounts are used. The AMT column is automatically tagged as AMOUNT in the Data tab.
- If administration ids have been defined in the treatment, an ADMID column is created, with “.” on all lines except on dose lines, on which administration ids are used. The ADMID column is automatically tagged as ADMINISTRATION ID in the Data tab.
- If an infusion duration or rate has been defined, a new INFDUR (for infusion duration) or INFRATE (for infusion rate) is created, with “.” on all lines

except on dose lines. The INF DUR column is automatically tagged as INFUSION DURATION in the Data tab, and the INFRATE column is automatically tagged as INFUSION RATE.

For each treatment, the dosing schedule can be defined as:

- **regular**: for regularly spaced dosing times, defined with the start time, inter-dose interval, and number of doses. A “repeat cycle” option allows to repeat the regular dosing schedule to generate a more complex regimen.
- **manual**: a vector of one or several dosing times, each defined manually. A “repeat cycle” option allows to repeat the manual dosing schedule to generate a more complex regimen.
- **external**: an external text file with columns id (optional), occasions (optional), time (mandatory), amount (mandatory), admid (administration id, optional), tinf or rate (optional), that allows to define individual doses.

While dose amounts, administration ids and infusion durations or rates are defined in the external file for external treatments, available options to define them for treatments of type “manual” or “regular” are:

- **“Manual”**: this applies the same amount (or administration id or infusion duration or rate) to all doses.
- **“By category”**: dose amounts (or administration id or infusion duration or rate) are defined manually for different categories read from the dataset.
- **“From data”**: dose amounts (or administration id or infusion duration or rate) are directly read from the dataset.

The options “by category” and “from data” are described in detail in [Section 8](#).

There is a “common settings” panel on the right:

- **dose intervals as occasions**: this creates a column to distinguish the dose intervals as different occasions (see [Section 9](#)).
- **infusion type**: If several treatments correspond to infusion administration, they need to share the same type of encoding for infusion information: as infusion duration or as infusion rate.

Example:

- demo **DoseAndLOQ_manual.pkx** (the screenshot below focuses on the formatting of doses and excludes other elements present in the demo):

In this demo, doses are initially not included in the dataset to format. A single dose at time 0 with an amount of 600 is added for each individual by Data Formatting. This creates a new AMT column in the formatted dataset, tagged as AMOUNT.

The top screenshot shows the 'Data file' tab with the file path 'C:/Users/jan/Desktop/multiple_BQ_tags_data_nocov.csv'. The 'Lines settings' tab shows 'Header lines: 1-2' and 'Lines to exclude:'. The bottom screenshot shows the 'Data' tab with a table of data. A callout box on the right shows the 'Treatment 1' settings for 'Time' and 'Amount'.

LINE NUMBER	ID	TIME_h	CONC_mg_L	STUDY
1	ID	TIME	CONC	STUDY
2		h	mg/L	
3	1	0	BLQ1	1
4	1	1	BLQ1	1
5	1	2	0.34	1
6	1	3	0.98	1

Records per page: 10 25 50 100 434
Page 1 of 9
Showing 1 to 50 of 434 entries

The bottom screenshot shows the 'Data' tab with a table of data. A callout box on the right shows the 'Treatment 1' settings for 'Time' and 'Amount'.

LINE NUMBER	ID	TIME_h	CONC_mg_L	STUDY	AMT
2	1	0	.	1	600
3	1	0	.	1	.
4	1	1	.	1	.
5	1	2	0.34	1	.
6	1	3	0.98	1	.

Records per page: 10 25 50 100 450
Page 1 of 9
Showing 1 to 50 of 450 entries - 18 individuals

The callout box on the right shows the 'Treatment 1' settings for 'Time' and 'Amount'.

Treatment 1

REGULAR MANUAL EXTERNAL

Time

0

Repeat cycle

Amount MANUAL FROM DATA 600

By category No category

Infusion

Administration id

PREVIEW

↑ Back to List of Formatting Steps

8. Reading censoring limits or dosing information from the dataset

When defining censoring limits for observations (see [Section 6](#)) or dose amounts, administration ids, infusion duration or rate for treatments (see [Section 7](#)), two options allow to define different values for different rows, based on information already present in the dataset: "by category" and "from data".

By category

It is possible to define manually different censoring limits, dose amounts, administration ids, infusion durations, or rates for different categories within a dataset's column. After selecting this column in the "By category" drop-down menu, the different modalities in the column are displayed and a value must be manually assigned each modality.

The left screenshot shows the 'Observation type based on CONC_mg_L' settings. The 'Interval censoring' section shows 'Upper limit' set to 'MANUAL' and 'BLQ'. The 'By category' dropdown is set to 'STUDY'. The table shows modalities and their upper limits.

By category	Modalities	Upper limit
STUDY	SD_400mg	0.01
	SD_500mg	0.1
	SD_600mg	0.1

The right screenshot shows the 'Treatment 1' settings for 'Time' and 'Amount'.

Treatment 1

REGULAR MANUAL EXTERNAL

Time

0

Repeat cycle

Amount MANUAL FROM DATA

By category

Modalities	AMOUNT
SD_400mg	400
SD_500mg	500
SD_600mg	600

- For censoring limits, the censoring limit used to replace each censoring tag depends on the modality on the same row.
- For doses, the value chosen for the newly created column (AMT for amount, ADMID for administration id, INF DUR for infusion duration, INFRATE for infusion rate) on each new dose line depends on the modality on the first row found in the dataset for the same individual and the same time as the dose, or the next time if there is no line in the initial dataset at that time, or the previous time if no time is found after the dose.

Example:

- demo **DoseAndLOQ_byCategory.pkx** (the screenshot below focuses on the formatting of doses and excludes other elements present in the demo):

In this demo there are three studies distinguished in the STUDY column with the categories "SD_400mg", "SD_500mg" and "SD_600mg". In Data Formatting, a single dose is manually defined at time 0 for all individuals, with different amounts depending the STUDY category. In addition, censoring interval is defined for the censoring tags BLQ, with an upper limit of the censoring interval (lower limit of quantification) that also depends

on the STUDY category. Three new columns – AMT for dose amounts, CENS for censoring tags (0 or 1), and LIMIT for the lower limit of the censoring intervals – are created by Data Formatting. A new dose line is then inserted at time 0 for each individual.

The screenshot displays the 'DoseAndLOQ_byCategory.pkx [demo]' interface. The 'Data formatting' tab is active, showing a 'Data file' section with a file path and a 'Lines settings' section with 'Header lines' (1-3) and 'Lines to exclude' (4-5). The 'Data' table shows columns: LINE NUMBER, ID, TIME_h, CONC_mg_L, STUDY. The 'Data information' section shows 'Units' set to 'Plasma' and 'Time' column set to 'hour'. The 'Data' table shows columns: LINE NUMBER, ID, TIME_h, CONC_mg_L, STUDY, AMT, CENS, LIMIT. The 'Observation type based on CONC_mg_L' section shows 'Interval censoring' settings. The 'Treatment 1' section shows 'Time' set to 0 and 'Amount' set to 400, 500, or 600. The 'PREVIEW' button is visible.

From data

The option "From data" is used to directly read censoring limits, dose amounts, administration ids, infusion durations, or rates from a dataset's column. The column must contain either numbers or numbers inside strings. In that case, the first number found in the string is extracted (including decimals with .).

- For censoring limits, the censoring limit used to replace each censoring tag is read from the selected column on the same row.
- For doses, the value chosen for the newly created column (AMT for amount, ADMID for administration id, INFDUR for infusion duration, INFRATE for infusion rate) on each new dose line is read from the selected column on the first row found in the dataset for the same individual and the same time as the dose, or the next time if there is no line in the initial dataset at that time, or the previous time if no time is found after the dose.

Example:

- demo **DoseAndLOQ_fromData.pkx** (the screenshot below focuses on the formatting of doses and censoring and excludes other elements present in the demo):

In this demo there are three studies distinguished in the STUDY column with the categories "SD_400mg", "SD_500mg" and "SD_600mg". In Data Formatting, a single dose is manually defined at time 0 for all individuals, with the amount read the STUDY column. In addition, censoring interval is defined for the censoring tags BLQ, with an upper limit of the censoring interval (lower limit of quantification) read from the LLOQ_mg_L column. Three new columns – AMT for dose amounts, CENS for censoring tags (0 or 1), and LIMIT for the lower limit of the censoring intervals – are created by Data Formatting. A new dose line is then inserted at time 0 for each individual, with amount 400, 500 or 600 for studies SD_400mg, SD_500mg and SD_600mg respectively.

Back to List of Formatting Steps






9. Creating occasions from dosing intervals

The option "Dose intervals as occasions" in the Treatments subtab of Data Formatting allows to create an occasion column to distinguish dose intervals. This is useful if the sets of measurements following different doses should be analyzed independently for a same individual.

Example:

- demo **doseIntervals_as_Occ.mlxtran** (Monolix demo in the folder 0.data_formatting, here imported into PKanalix):

This demo imported from a Monolix demo has an initial dataset in Monolix-standard format, with multiple doses encoded as dose lines with dose amounts in the AMT column. When using this dataset directly into Monolix or PKanalix, a single analysis is done on each individual concentration profile considering all doses, which means that NCA would be done on the concentrations after the last dose only, and modeling (CA in PKanalix or population modeling in Monolix) would be estimated with a single set of parameter values for each individual. If instead we want to run separate analyses on the sets of concentrations following each dose, we need to distinguish them as occasions with a new column added with the Data Formatting module. To this end, we define the same treatment as in the initial dataset with Data Formatting (here as regular multiple doses) with the option "Dose intervals as occasions" selected. After clicking Preview, Data Formatting adds two new columns: an AMT1 column with the new doses, to be tagged as AMOUNT instead of the AMT column that will now be ignored, and a DOSE_OCC column to be tagged as OCCASION.





Back to List of Formatting Steps

10. Handling urine data

In PKanalix-standard format, the start and end times of urine collection intervals must be recorded in a single column, tagged as **TIME** column-type, where the end time of an interval automatically acts as start time for the next interval (see [here](#) for more details). If a dataset contains start and end times in two different columns, they can be merged into a single column by Data Formatting. This is done automatically by tagging these two columns as **START** and **END** in the Initialization subtab of Data Formatting (see [Section 2](#)). In addition the column containing urine collection volume must be tagged as **VOLUME**.

Example:

- demo **Urine_LOQinObs.plx**:



Back to List of Formatting Steps

11. Adding new columns from an external file

The last subtab is used to insert additional columns in the dataset from a separate file. The external file must contain a table with a column named **ID** or **id** with the same subject identifiers as in the dataset to format, and other columns with a header name and individual values (numbers or strings).

There can be only one value per individual, which means that the additional columns inserted in the formatted dataset can contain only a constant value within each individual, and not time-varying values.

Examples of additional columns that can be added with this option are:

- individual parameters estimated in a previous analysis, to be read as regressors to avoid estimating them. **Time-varying regressors are not handled.**
- new covariates.

If occasions are defined in the formatted dataset, it is possible to have an occasion column in the external file and values defined per subject-occasion.

Example:

- demo **warfarin_PKPDseq_project.mlxtran** (Monolix demo in the folder 0.data_formatting, here imported into PKanalix):

This demo imported from a Monolix demo has an initial PKPD dataset in Monolix-standard format. The option “Additional columns” is used to add the PK parameters estimated on the PK part of the data in another Monolix project.

The top screenshot shows the PKanalix interface with the 'Additional columns' tab selected. A file named 'warfarinPK_regressors.txt' is loaded, and a table of PK parameters is displayed:

id	Tlag_mode	ka_mode	V_mode	Cl_mode
1	0.881068	1.27612	7.68831	0.111825
2	0.739306	0.906915	8.35584	0.128526
3	0.909005	0.955194	8.46397	0.112645
4	0.893931	1.42947	5.08504	0.072557
5	0.851363	1.36498	11.5598	0.176061

The bottom screenshot shows the 'Data' tab with the formatted dataset loaded. The table displays columns for ID, TIME, AMOUNT, OBSERVATION, and various covariates, including the PK parameters from the previous table:

LINE NUMBER	id	time	amt	dv	dvid	wt	sex	age	Tlag_mode	ka_mode	V_mode	Cl_mode
1	1	0	100			66.7	1	50	0.881068	1.27612	7.68831	0.111825
2	1	0		98	2	66.7	1	50	0.881068	1.27612	7.68831	0.111825
3	1	24		9.2	1	66.7	1	50	0.881068	1.27612	7.68831	0.111825


[Back to List of Formatting Steps](#)

12. Exporting the formatted dataset

Once data formatting is done and the new dataset is accepted, the project can be saved and it is possible to export the formatted dataset as a csv file from the main menu Export > Export formatted data.



Note that if you did some data formatting directly in Monolix (instead of PKanalix), the possibility to save the project and export the formatted data is enabled only after loading a model in the structural model tab.



↑ Back to List of Formatting Steps

2.3. Units: display and scaling

- It **displays units of the NCA and CA parameters** calculated automatically from information provided by a user about units of measurements in a dataset. Units are shown in the results tables, on plots and are added to the saved files.
- It allows for **scaling** values of a dataset **to represent outputs in preferred units**, which facilitates the interpretation of the results. Scaling is done in the PKanalix data manager, and does not change the original data file.

- Units definition
- Units of the NCA parameters
- Units of the CA parameters
- Units display
- Units preferences

LINE NUMBER	ID	TIME_h	PERIOD	AMT_mg	ADM	Conc_mg_L	CENS	LIMIT	FORM	AGE	WT
1#1									ref	25	66
2	1	0	1	600	1				ref	25	66
3	1	0	1			0.06	intervalCensored	0	ref	25	66
4	1	1	1			0.06	intervalCensored	0	ref	25	66
5	1	2	1			0.34	notCensored		ref	25	66
6	1	3	1			0.98	notCensored		ref	25	66
7	1	4	1			1.44	notCensored		ref	25	66
8	1	5	1			1.81	notCensored		ref	25	66

Units definition

Units of the NCA and CA parameters are considered as combinations of units of: time, amount and volume. For instance, unit of AUC is [concentration unit * time unit] = [amount unit / volume unit * time unit]. These quantities correspond to columns in a typical PK dataset: time of measurements, observed concentration, dose amount (and volume of collected urine when relevant).

The "Units" block allows to define the preferred units for the output NCA parameters (purple frame below), which are related to the units of the data set columns (green frame below) via scaling factors (in the middle). The output units are displayed in results and plots after running the NCA and CA tasks.

Input data	Scaling factor	Output units
TIME column × 1	=	hour
CONC column × 1	=	milligram / liter
AMT column × 1	=	milligram / kilogram

Input data

Output units

In PKanalix version 2023R1, the amount can now also be specified in mass of the administered dose per body mass or body surface area. This eliminates the need to edit NCA or CA parameters by post-processing unit conversion.

Concentration is a quantity defined as "amount per volume". It has two separate units, which are linked (and equal) to AMT and VOLUME respectively. Changing the amount unit of CONC will automatically update the amount unit of AMT. This constraint allows to apply simplifications to the output parameters units, for instance have Cmax_D as [mL] and not [mg/mL/mg].

Input data	Scaling factor	Output units
TIME column × 1	=	hour
CONC column × 1	=	nanogram / milliliter
AMT column × 1	=	nanogram
VOLUME column × 1	=	milliliter

Units without data conversion: Output units correspond to units of the input data.

In other words, desired units of the NCA and CA parameters correspond to units of measurements in the dataset. In this case, select from the drop-down menus units of the input data and keep the default scaling factor (equal to one). All calculations are performed using values in the original dataset and selected units are displayed in the results tables and on plots.

Units conversion: output units are different from the units of the input data.

NCA and CA parameters can be calculated and displayed in any units, not necessarily the same as used in the dataset. The scaling factors (by default equal to 1) multiply the corresponding columns in a dataset and transform them to a dataset in new units. PKanalix shows data with new, scaled values after having clicked on the "accept" button. **This conversion occurs only internally and the original dataset file remains unchanged.** So, in this case, select desired output units from the list and, knowing units of the input data, scale the original values to the new units.

Example 1: input time units in [hours] and desired output units in [days]

For instance, let measurement times in a dataset be in hours. To obtain the outputs in days, set the output time unit as "days" and the scaling factor to 1/24, as shown below. It reads as follows:

(values of time in **hours** from the dataset) * (1/24) = time in **days**.

Input data	Scaling factor	Output units
TIME column × 1	= 1/24	day
CONC column × 1	= 24	milliliter
AMT column × 1	= 1/24	kilogram

After accepting the scaling, a dataset shown in the Data tab is converted internally to a new data set. It contains original values multiplied by scaling factors. Then, all computations in the NCA and CA tasks are performed using new values, so that the results correspond to selected units.

Example 2: input data set concentration units in [ng/mL] and amount in [mg]

	ID	TIME	AMOUNT	OBSERVATION
LINE NUMBER	ID	time	Amount	DV
2	271	0	100	.
3	271	0	.	2950
4	271	0.17	.	2164
5	271	0.25	.	1884
6	271	0.5	.	1366
7	271	1	.	844
8	271	1.5	.	625
9	271	2	.	478
10	271	4	.	256
11	271	6	.	147

Let's assume that we have a data set where the

concentration units are [ng/mL]=[ug/L] and the dose amount units are [mg], as presented above. It is not possible to indicate these units directly, as the amount unit of CONC and AMT must be the same. One option would be to indicate the CONC as [mg/kL] and the AMT as [mg] but having the volume unit as [kL] is not very convenient. We will thus use the scaling factors to harmonize the units of the concentration and the amount.

If we would like to have the output NCA parameters using [ug] and [L], we can define the CONC units as [μg/L] (as the data set input units) and scale the AMT to convert the amount column from [mg] to [ug] unit with a scaling factor of 1000. After clicking "accept" on the bottom right, the values in the amount column of the (internal) data set have been multiplied by 1000 by PKanalix such that they are now in [μg].

Data file

C:\Users\koolman\OneDrive\pkanalix\pkanalix\2023\1\demo\2\case_study\data\N2000_1\bolus_singledose.csv

BROWSE...

Data information

☐ Urine ☒ Plasma

Units

TIME column * 1 = hour

CONC column * 1 = microgram / liter

AMT column * 1000 = microgram /

Add covariates

NONE

LINE NUMBER	ID	time	Amount	DV	GENDER	AGE	HT	WT
1	271		100000		F	69	157	63.7
2	271	0			F	69	157	63.7
3	271	0		2950	F	69	157	63.7
4	271	0.17		2164	F	69	157	63.7
5	271	0.25		1884	F	69	157	63.7
6	271	0.5		1366	F	69	157	63.7
7	271	1		844	F	69	157	63.7

Amount in data set has been multiplied by 1000 by PKanalix to obtain μg

If we would like to have the output NCA parameters using [ng] and [mL], we can define the CONC units as [ng/mL] (as the data set input units) and scale the AMT to convert the amount column from [mg] to [ng] unit with a scaling factor of 1000000.

Data file

C:\Users\koolman\OneDrive\pkanalix\pkanalix\2023\1\demo\2\case_study\data\N2000_1\bolus_singledose.csv

BROWSE...

Data information

☐ Urine ☒ Plasma

Units

TIME column * 1 = hour

CONC column * 1 = nanogram / milliliter

AMT column * 1000000 = nanogram /

Add covariates

NONE

LINE NUMBER	ID	time	Amount	DV	GENDER	AGE	HT	WT
1	271		1e+08		F	69	157	63.7
2	271	0			F	69	157	63.7
3	271	0		2950	F	69	157	63.7
4	271	0.17		2164	F	69	157	63.7
5	271	0.25		1884	F	69	157	63.7
6	271	0.5		1366	F	69	157	63.7
7	271	1		844	F	69	157	63.7

Units for the NCA parameters

Rs _q	no unit
Rs _q _adjusted	no unit
Corr_XY	no unit
No_points_lambda_z	no unit
Lambda_z	time-1
Lambda_z_lower	time
Lambda_z_upper	time
HL_lambda_z	time
Span	no unit

Lambda_z_intercept	no unit
T0	time
Tlag	time
Tmax_Rate	time
Max_Rate	amount.time-1
Mid_Pt_last	time
Rate_last	amount.time-1
Rate_last_pred	amount.time-1
AURC_last	amount
AURC_last_D	grading
Vol_UR	volume
Amount_recovered	amount
Percent_recovered	% [not calculable when grading -> set as NaN]
AURC_all	amount
AURC_INF_obs	amount
AURC_PerCentExtrap_obs	%
AURC_INF_pred	amount
AURC_PerCentExtrap_pred	%
C0	amount.volume-1
Tmin	time
Cmin	amount.volume-1
Tmax	time
Cmax	amount.volume-1
Cmax_D	grading.volume-1
Tlast	time
Clast	amount.volume-1
AUClast	time.amount.volume-1
AUClast_D	grading.time.volume-1
AUMClast	time2.amount.volume-1
AUCall	time.amount.volume-1
AUCINF_obs	time.amount.volume-1
AUCINF_D_obs	grading.time.volume-1
AUCINF_pred	time.amount.volume-1
AUCINF_D_pred	grading.time.volume-1
AUC_PerCentExtrap_obs	%
AUC_PerCentBack_Ext_obs	%
AUMCINF_obs	time2.amount.volume-1
AUMC_PerCentExtrap_obs	%

Vz_F_obs	volume.grading-1
Cl_F_obs	volume.time-1.grading-1
Cl_obs	volume.time-1.grading-1
Cl_pred	volume.time-1.grading-1
Vss_obs	volume.grading-1
Clast_pred	amount.volume-1
AUC_PerCentExtrap_pred	%
AUC_PerCentBack_Ext_pred	%
AUMCINF_pred	time2.amount.volume-1
AUMC_PerCentExtrap_pred	%
Vz_F_pred	volume.grading-1
Cl_F_pred	volume.time-1.grading-1
Vss_pred	volume.grading-1
Tau	time
Ctau	amount.volume-1
Ctrough	amount.volume-1
AUC_TAU	time.amount.volume-1
AUC_TAU_D	grading.time.volume-1
AUC_TAU_PerCentExtrap	%
AUMC_TAU	time2.amount.volume-1
Vz	volume.grading-1
Vz_obs	volume.grading-1
Vz_pred	volume.grading-1
Vz_F	volume.grading-1
CLss_F	volume.time-1.grading-1
CLss	volume.time-1.grading-1
CAvg	amount.volume-1
FluctuationPerCent	%
FluctuationPerCent_Tau	%
Accumulation_index	no unit
Swing	no unit
Swing_Tau	no unit
Dose	amount.grading-1
N_Samples	no unit
MRTlast	time
MRTINF_obs	time
MRTINF_pred	time
AUC_lower_upper	time.amount/volume

AUC_lower_upper_D	grading.time.volume-1
CAVG_lower_upper	amount/volume
AURC_lower_upper	amount
AURC_lower_upper_D	grading


Units for the PK parameters


Units are available only in PKanalix. When exporting a project to Monolix, values of PK parameters are re-converted to the original dataset unit. Below, volumeFactor is defined implicitly as: $\text{volumeFactor} = \text{amtFactor} / \text{concFacto}$, where “factor” is the scaling factor used in the “Units” block in the Data tab.

PARAMETER	UNIT	INVERSE of UNITS
V (V1, V2, V3)	volume.grading-1	value/volumeFactor
k (ka, k12, k21, k31, k13, Ktr)	1/time	value*timeFactor
Cl, Q (Q2, Q3)	volume/time.grading-1	value*timeFactor/volumeFactor
Vm	amount/time.grading-1	value*timeFactor/amountFactor
Km	concentration	value/concFactor
T (Tk0, Tlag, Mtt)	time	value/timeFactor
alpha, beta, gamma	1/time	value*timeFactor
A, B, C	grading/volume	value*volumeFactor

Units display

To visualise units, switch on the “units toggle” and accept a dataset. Then, after running the NCA and CA tasks, units are displayed:

- In the results table next to the parameter name (violet frame), in the table copied with “copy” button  (blue frame) and in the saved .txt files in the result folder (green frame).



The screenshot shows the PKanalix interface. The 'Results' tab is active, displaying a table of individual estimates. A blue box highlights the 'Units' column, and a green box highlights the 'Copy' button. A Notepad window shows the copied data with units.

id	AUCINF_obs (h·g·mL ⁻¹)	AUC_%Extrap_obs (%)	AUClast (h·g·mL ⁻¹)	CL_F_obs (mL·h ⁻¹)	Clast (g·mL ⁻¹)	Cmax (g·mL ⁻¹)	Lambda_z (h ⁻¹)	Tlast (h)	Tmax (h)	Vz_F_obs (mL)
101	249.91	8.31	229.14	0.6	2.29	20.48	0.11	24	2	5.45
102	160.81	10.26	144.31	0.93	1.6	14.14	0.097	24	2.5	9.65
103	251.54	13.43	217.76	0.6	3.12	19.87	0.092	24	2.5	6.46
104	213.71	11.29	189.58	0.7	2.28	22.38	0.095	24	3	7.42
105	135.83	16.51	113.41	1.1	1.52	13.14	0.068	24	2.5	16.33
106	115.87	6.11	108.79	1.29	0.81	12.28	0.11	24	3	11.27
107	239.48	10.19	195.93	0.63	3.37	17.51	0.077	24	2.5	8.09
108	175.63	9.05	159.74	0.85	1.8	14.39	0.11	24	2	7.54
109	136.09	19.24	118.07	1.1	1.53	12.13	0.085	24	3	13
110	229.15	8.35	210.02	0.65	2.22	16.67	0.12	24	1	5.65

The Notepad window shows the copied data with units:


```
id,AUCINF_obs,AUC_PercentExtrap_obs,AUClast,CL_F_obs,Clast,Cmax,Lambda_z,Tlast,Tmax,Vz_F_obs
id,AUCINF_obs,AUClast,CL_F_obs,Clast,Cmax,Lambda_z,Tlast,Tmax,Vz_F_obs
h·g·mL-1,%,h·g·mL-1,mL·h-1,g·mL-1,g·mL-1,h-1,h,h,mL
101,249.9078,8.30835,229.143,0.600224,2.286,20.4776,0.110099,24,2,5.45167
102,160.811,10.2615,144.309,0.932775,1.59526,14.136,0.096673,24,2.5,9.64876
103,251.542,13.4281,217.764,0.596323,3.11717,19.8727,0.0922858,24,2.5,6.4617
104,213.707,11.2908,189.578,0.701897,2.28098,22.383,0.0945321,24,3,7.42496
105,135.834,16.5097,113.408,1.10429,1.51674,13.1413,0.067634,24,2.5,16.3275
106,115.87,6.10671,108.795,1.29455,0.81243,12.2776,0.114817,24,3,11.2749
107,239.481,10.1871,195.926,0.626355,3.3704,17.5126,0.0773837,24,2.5,8.09427
```

- On plots if the “units” display is switched on in the General plots settings



Units preferences

The units selected by default when starting a new PKnaalix project can be chosen in the menu Settings > Preferences.



2.4. Filtering a data set

Starting on the 2020 version, filtering a data set to only take a subpart into account in your modelization is possible. It allows to make filters on some specific IDs, times, measurement values,... It is also possible to define complementary filters and also filters of filters. It is accessible through the filters item on the data tab.

- Creation of a filter
- Filtering actions
- Filters with several actions
- Other filters: filter of filter and complementary filters

LINE NUMBER	ID	TIME	DV	AMT	SEX	WEIGHT
1	1	0		100	F	70.6466
2	1	0		100	F	70.6466
3	1	1	16.6314		F	70.6466
4	1	1.5	21.5208		F	70.6466
5	1	2	11.746		F	70.6466
6	1	2.5	20.0199		F	70.6466
7	1	3	11.9624		F	70.6466
8	1	4	10.7658		F	70.6466

Creation of a filter

To create a filter, you need to click on the data set name. You can then create a “child”. It corresponds to a subpart of the data set where you will define your filtering actions.

action	header	operator	value
select ids	ID	=	1

You can see on the top (in the green rectangle) the action that you will complete and you can CANCEL, ACCEPT, or ACCEPT & APPLY with the bottoms on the bottom.

Filtering a data set: actions

In all the filtering actions, you need to define

- An **action**: it corresponds to one of the following possibilities: select ids, remove ids, select lines, remove lines.
- A **header**: it corresponds to the column of the data set you wish to have an action on. Notice that it corresponds to a column of the data set that was tagged with a header.
- An **operator**: it corresponds to the operator of choice (=, ≠, <, >, or ≥).
- A **value**. When the header contains numerical values, the user can define it. When the header contains strings, a list is proposed.

For example, you can

- Remove the ID 1 from your study:

action	header	operator	value
remove ids	ID	=	1

In that case, all the IDs except ID = 1 will be used for the study.

- Select all the lines where the time is less or equal 24:

action	header	operator	value
select lines	TIME	≤	24

In that case, all lines with time strictly greater than 24 will be removed. If a subject has no measurement anymore, it will be removed from the study.

- Select all the ids where SEX equals F:

action	header	operator	value
select ids	SEX	=	F

In that case, all the male will be removed of the study.

- Remove all Ids where WEIGHT less or equal 65:

action	header	operator	value
remove ids +	WEIGHT	≤	65

In that case, only the subjects with a weight over 65 will be kept for the study.

In any case, the interpreted filter data set will be displayed in the data tab.

Filters with several actions

In the previous examples, we only did one action. It is also possible to do several actions to define a filter. We have the possibility to define UNION and/or INTERSECTION of actions.

INTERSECTION

By clicking by the + and – button on the right, you can define an intersection of actions. For example, by clicking on the +, you can define a filter corresponding to intersection of

- The IDs that are different to 1.
- The lines with the time values less than 24.

action	header	operator	value
select ids +	ID	≠	1
select lines +	TIME	≤	24

Thus in that case, all the lines with a time less than 24 and corresponding to an ID different than 1 will be used in the study. If we look at the following data set as an example

ID	AMT	TIME	CONC	WEIGHT	SEX
1	4.02	0		79.6	M
1		0.25	2.84	79.6	M
1		0.57	6.57	79.6	M
1		1.12	10.5	79.6	M
1		2.02	9.66	79.6	M
1		3.82	8.58	79.6	M
1		5.1	8.36	79.6	M
1		7.03	7.47	79.6	M
1		9.05	6.89	79.6	M
1		12.12	5.94	79.6	M
1		24.37	3.28	79.6	M
2	4.4	0		72.4	M
2		0.27	1.72	72.4	M
2		0.52	7.91	72.4	M
2		1	8.31	72.4	M
2		1.92	8.33	72.4	M
2		3.5	6.85	72.4	M
2		5.02	6.08	72.4	M
2		7.03	5.4	72.4	M
2		9	4.55	72.4	M
2		12	3.01	72.4	M
2		24.3	0.9	72.4	M
3	4.53	0		70.5	M
3		0.27	4.4	70.5	M
3		0.58	6.9	70.5	M
3		1.02	8.2	70.5	M
3		2.02	7.8	70.5	M
3		3.62	7.5	70.5	M
3		5.08	6.2	70.5	M
3		7.07	5.3	70.5	M
3		9	4.9	70.5	M
3		12.15	3.7	70.5	M
3		24.17	1.05	70.5	M

Initial data set

Resulting data set after action: select IDs ≠ 1

ID	AMT	TIME	CONC	WEIGHT	SEX
2	4.4	0		72.4	M
2		0.27	1.72	72.4	M
2		0.52	7.91	72.4	M
2		1	8.31	72.4	M
2		1.92	8.33	72.4	M
2		3.5	6.85	72.4	M
2		5.02	6.08	72.4	M
2		7.03	5.4	72.4	M
2		9	4.55	72.4	M
2		12	3.01	72.4	M
2		24.3	0.9	72.4	M
3	4.53	0		70.5	M
3		0.27	4.4	70.5	M
3		0.58	6.9	70.5	M
3		1.02	8.2	70.5	M
3		2.02	7.8	70.5	M
3		3.62	7.5	70.5	M
3		5.08	6.2	70.5	M
3		7.07	5.3	70.5	M
3		9	4.9	70.5	M
3		12.15	3.7	70.5	M
3		24.17	1.05	70.5	M

Resulting data set after action: select lines with time ≤ 24

Resulting data set for the study as the intersection of the two actions

ID	AMT	TIME	CONC	WEIGHT	SEX
1	4.02	0		79.6	M
1		0.25	2.84	79.6	M
1		0.57	6.57	79.6	M
1		1.12	10.5	79.6	M
1		2.02	9.66	79.6	M
1		3.82	8.58	79.6	M
1		5.1	8.36	79.6	M
1		7.03	7.47	79.6	M
1		9.05	6.89	79.6	M
1		12.12	5.94	79.6	M
1		24.37	3.28	79.6	M
2	4.4	0		72.4	M
2		0.27	1.72	72.4	M
2		0.52	7.91	72.4	M
2		1	8.31	72.4	M
2		1.92	8.33	72.4	M
2		3.5	6.85	72.4	M
2		5.02	6.08	72.4	M
2		7.03	5.4	72.4	M
2		9	4.55	72.4	M
2		12	3.01	72.4	M
2		24.3	0.9	72.4	M
3	4.53	0		70.5	M
3		0.27	4.4	70.5	M
3		0.58	6.9	70.5	M
3		1.02	8.2	70.5	M
3		2.02	7.8	70.5	M
3		3.62	7.5	70.5	M
3		5.08	6.2	70.5	M
3		7.07	5.3	70.5	M
3		9	4.9	70.5	M
3		12.15	3.7	70.5	M
3		24.17	1.05	70.5	M

UNION

By clicking by the + and – button on the bottom, you can define an union of actions. For example, in a data set with a multi dose, I can focus on the first and the last dose. Thus, by clicking on the +, you can define a filter corresponding to union of

- The lines where the time is strictly less than 12.
- The lines where the time is greater than 72.

select lines +	time	<	12
UNION			
select lines +	time	>	72
UNION			
select lines +	amt	=	40

id	time	amt	y
1	0	40	.
1	0.5	.	3.365
1	1	.	2.4718
1	3	.	2.5617
1	5	.	1.8819
1	7	.	2.1349
1	9	.	1.6741
1	11	.	2.0317
1	12	40	.
1	12.5	.	4.219
1	13	.	3.3258
1	15	.	3.4157
1	17	.	2.7359
1	19	.	2.9889
1	21	.	2.5281
1	23	.	2.8857
1	24	40	.
...
1	72	40	.
1	72.5	.	6.0911
1	73	.	5.874
1	75	.	6.0025
1	77	.	4.7898
1	79	.	4.1839
1	81	.	4.2116
1	83	.	3.2958

Initial data set

Resulting data set after action:
select lines where the time is strictly less than 12

id	time	amt	y
1	0	40	.
1	0.5	.	3.365
1	1	.	2.4718
1	3	.	2.5617
1	5	.	1.8819
1	7	.	2.1349
1	9	.	1.6741
1	11	.	2.0317
1	12	40	.
1	12.5	.	4.219
1	13	.	3.3258
1	15	.	3.4157
1	17	.	2.7359
1	19	.	2.9889
1	21	.	2.5281
1	23	.	2.8857
1	24	40	.
...
1	72	40	.
1	72.5	.	6.0911
1	73	.	5.874
1	75	.	6.0025
1	77	.	4.7898
1	79	.	4.1839
1	81	.	4.2116
1	83	.	3.2958

Resulting data set after action:
select lines where the time is greater than 72

id	time	amt	y
1	0	40	.
1	0.5	.	3.365
1	1	.	2.4718
1	3	.	2.5617
1	5	.	1.8819
1	7	.	2.1349
1	9	.	1.6741
1	11	.	2.0317
1	12	40	.
1	12.5	.	4.219
1	13	.	3.3258
1	15	.	3.4157
1	17	.	2.7359
1	19	.	2.9889
1	21	.	2.5281
1	23	.	2.8857
1	24	40	.
...
1	72	40	.
1	72.5	.	6.0911
1	73	.	5.874
1	75	.	6.0025
1	77	.	4.7898
1	79	.	4.1839
1	81	.	4.2116
1	83	.	3.2958

Resulting data set after action:
select lines where amt equals 40

id	time	amt	y
1	0	40	.
1	12	40	.
1	24	40	.
1	72	40	.

id	time	amt	y
1	0	40	.
1	0.5	.	3.365
1	1	.	2.4718
1	3	.	2.5617
1	5	.	1.8819
1	7	.	2.1349
1	9	.	1.6741
1	11	.	2.0317
1	12	40	.
1	12.5	.	4.219
1	13	.	3.3258
1	15	.	3.4157
1	17	.	2.7359
1	19	.	2.9889
1	21	.	2.5281
1	23	.	2.8857
1	24	40	.
...
1	72	40	.
1	72.5	.	6.0911
1	73	.	5.874
1	75	.	6.0025
1	77	.	4.7898
1	79	.	4.1839
1	81	.	4.2116
1	83	.	3.2958

Considered data set for the study
as the union of the three actions

Notice that, if just define the first two actions, all the dose lines at a time in]12, 72[will also be removed. Thus, to keep having all the doses, we need to add the condition of selecting the lines where the dose is defined.

In addition, it is possible to do any combination of INTERSECTION and UNION.

Other filters: filter of filter and complementary filters

Filtering a data set can be nested. Based on the definition of a filter, by clicking on the filter, it is possible to create:

- A **child**: it corresponds to a new filter with the initial filter as the source data set.
- A **complement**: corresponds to the complement of the filter. For example, if you defined a filter with only the IDs where the SEX is F, then the complement corresponds to the IDs where the SEX is not F.



2.5. Import a Monolix/Simulx project

MonolixSuite applications are interconnected and projects can be exported/imported between different applications. This interconnection is guaranteed by using the same model syntax (mlxtran language) and the same dataset format.


There are two options to create a PKanalix project from a Monolix or Simulx project*

- **Import a Monolix project to PKanalix** to run NCA or CA on the original dataset in the Monolix project
- **Export from Monolix to PKanalix** to run NCA or CA on the original dataset, or on the VPC simulations or individual fits
 - Export with individual fits dataset.
 - Export with VPC simulations
- **Export from Simulx to PKanalix** to run NCA or CA on simulations

*Note: Import and export to PKanalix is available starting from the 2023 version. In the previous versions, you can only export a PKanalix project to Monolix.

Import a Monolix project to PKanalix

To import a Monolix project, open PKanalix application and in the home page select IMPORT FROM:



In PKanalix, it is possible to import only a Monolix project with its original dataset. When you select in the home screen IMPORT FROM > Monolix, then PKanalix will create a new, untitled, project with:


- Dataset tagging and formatting steps if present as in the Monolix project.
- Dataset filters if applied as in the Monolix project.
- NCA settings: "administration type" set to intravenous or extravascular depending on the administration in a model selected in the Monolix project.
- NCA settings: "observation ID to use" set to the first one alphabetically (if obsid are strings) or numerically (if obsid are integers).
- NCA settings: default PKanalix settings for the integral method, treatment of BLQ values, parameters.
- Acceptance criteria: not selected.
- Bioequivalence: default PKanalix settings
- CA model: set to the same structural model and the same mapping between observation ids and model outputs as in the Monolix project
- CA initial parameter values: equal to the initial estimates of the Monolix project
- CA parameters constraints: none for normally distributed parameters, positive for log-normally distributed parameters; bounded with limits imported from the Monolix project for logit-normally and probit-normally distributed parameters.
- CA calculations settings: default PKanalix settings

After the import, you can save the PKanalix project as usual, edit it and run the analysis.

Remark: Importing a Monolix project to PKanalix is only with the original dataset used to create a Monolix project. Export from Monolix to PKanalix has more options. You can choose to import a Monolix project with its original dataset, vpc dataset or with individual fits dataset, see the next section. Import from Simulx is currently unavailable.

Export from Monolix to PKanalix

To export a Monolix project to PKanalix click EXPORT PROJECT TO in the top menu "Export":



When you export a Monolix project, in the export pop-up window you can choose:

- to which application you want to export your current project: PKanalix or Simulx: select "PKanalix"
- which dataset you want to use in the export: original dataset, vpc dataset, individual fits dataset

By default, Monolix will copy all files, e.g. dataset and model, next to the new PKanalix project. To keep current location of these files, switch the toggle “Generated files next to project” off. Click the “EXPORT” button at the bottom to confirm. PKanalix application will open automatically with a predefined project called “untitled”. It will contain the same pre-defined elements as during the import, see the description above. If you selected export with vpc dataset or individual fits dataset, then Monolix creates a .csv files with vpc simulations or individual model predictions in the MonolixSuite data format. These dataset will be automatically loaded in the new PKanalix project.

Export with individual fits dataset

The .csv generated file contains model predictions on a fine time grid (specified in the Plots tasks settings in Monolix) obtained with different individual parameters:

	ID	TIME	IGNORE	IGNORE	IGNORE	OBSERVATION	AMOUNT	CONTINUOUS COVARIATE	CONTINUOUS COVARIATE	CATEGORICAL COVARIATE
LINE NUMBER	ID	TIME	popPred_medianCOV	popPred	indivPred_mean	indivPred_mode	amt	age	wt	sex
1	100							50	66.7	1
2	100	0					100	50	66.7	1
3	100	0	0	0	0	0		50	66.7	1
4	100	1.5e-15	0	0	0	0		50	66.7	1

- **popPred**: model predictions with individual parameters corresponding to the population parameters estimated by Monolix and the impact of the individual covariates, but random effects set to zero. For instance, $V_i = V_{pop} \times \exp^{\beta_{V,WT} \times WT_i}$
- **popPred_medianCOV**: model predictions with individual parameters corresponding to the population parameters and the impact of the covariates using the population median for continuous covariates and reference category for categorical covariates, and random effects set to zero. For instance, $V_i = V_{pop} \times \exp^{\beta_{V,WT} \times WT_{med}}$
- **indivPred_mode**: [when EBEs task has run] model predictions with individual parameters corresponding to the EBEs (conditional mode) estimated by Monolix (as displayed in Monolix Results > Indiv. param > Cond. mode.). Tagged as OBSERVATION by default.
- **indivPred_mean**: [when COND. DISTRIB. task has run] model predictions with individual parameters corresponding to the conditional mean estimated by Monolix (as displayed in Monolix Results > Indiv. param > Cond. mean.)

The dataset contains also the individual design: doses, occasions, regressors, covariates.

Export with vpc dataset



The .csv generated dataset in the Monolix format constructed from the VPC simulations and individual design (doses, occasions, regressors, covariates)

	STRATIFICATION CATEGORICAL COVARIATE	IGNORE	ID	TIME	IGNORE	IGNORE	OBSERVATION	OBSERVATION ID	AMOUNT	CONTINUOUS COVARIATE	CONTINUOUS COVARIATE	CATEGORICAL COVARIATE
LINE NUMBER	rep	ID	uid	TIME	timeAfterLastDose	doseid	OBSERVATION	obsid	amt	age	wt	sex
1	1		100_rep1							50	66.7	1
2	1	100	100_rep1	0	0	1			100	50	66.7	1
3	1	100	100_rep1	0.5	0.5	1	0.0582826	y1		50	66.7	1
4	1	100	100_rep1	1	1	1	-0.296608	y1		50	66.7	1

- **rep**: simulation replicate. Tagged as stratification categorical covariate
- **uid**: unique id obtained by concatenation of “rep” and “id”. Tagged as ID by default.
- **time/timeAfterLastDose**: observation times or times after last dose
- **doseid**: administration id
- **y1/observation**: simulated observation values. Default header is “y1” if a model has only one output, and “observation” in case of several outputs (obsid column is created automatically). Tagged automatically as OBSERVATION


3. Non Compartmental Analysis

One of the main feature of PKanalix is the calculation of the parameters in the Non Compartmental Analysis framework.

- [NCA task](#)
- [NCA results](#) 
- [NCA plots](#) 
- [NCA outputs](#)

NCA task

Main function of PKanalix is to calculate PK metrics (non-compartmental parameters). In the Tasks tab, you find a dedicated section “NCA” (left panel). There are two subtabs: RUN and CHECK λ_z .



- RUN tab contains section with calculation task buttons (for non-compartmental analysis and [bioequivalence](#)) and the settings: calculation, acceptance criteria, parameters and bioequivalence design. The meaning of all the settings and their default is in a separate page [here](#).
- CHECK λ_z allows to define the calculation of λ_z – it is a graphical preview of each individual data and the regression line, more details are [here](#).

To calculate NCA parameters click on the button NON-COMPARTMENTAL ANALYSIS. If you click on RUN button, then PKanalix runs all selected tasks (“in green”). Results and plots are automatically generated in separate tabs.

NCA results

When you run the NCA task, results are automatically generated and displayed in the RESULTS tab. There are three tables: individual estimates, summary and points for λ_z .

Non compartmental analysis results per individual

Individual estimates of the NCA parameters are displayed in the table in the tab “INDIV. ESTIM.”

project_aPCS9_SAD.plx [demo] - PKanalix - 2023R1

Project Settings Export Help

Data Tasks Results Plots

NCA

POINTS FOR λ_z SUMMARY INDIV. ESTIM.

Individual estimates

Non compartmental analysis results per individual

Records per page: 10 18

Page 1 of 2

Showing 1 to 10 of 18 entries

id	AUC _{0-∞} /Dose (d·L ⁻¹)	AUC _{0-∞} (d·mg·L ⁻¹)	AUC _{0-tlast} (d·mg·L ⁻¹)	CL/F (L·d ⁻¹)	C _{last} (mg·L ⁻¹)	C _{max} (mg·L ⁻¹)	C _{max} /Dose (L ⁻¹)	K _{el} (d ⁻¹)	Intercept (d ⁻¹)	V _d /F (L)	DOSE_mg
1	2.26	338.76	332.66	0.44	0.42	19.72	0.13	0.069	3.23	6.43	150mg
2	1.8	270.59	269.76	0.55	0.12	12.52	0.063	0.15	6.13	3.8	150mg
3	1.32	197.78	194.52	0.76	0.39	11.03	0.074	0.12	4.09	6.32	150mg
4	0.91	136.34	133.84	1.1	0.27	8.04	0.054	0.11	3.25	10.16	150mg
5	0.96	144.35	141.83	1.04	0.29	8.35	0.056	0.11	3.6	9.04	150mg
6	1.16	174.44	170.51	0.86	0.46	8.68	0.058	0.12	4.15	7.34	150mg
7	0.95	285.36	283.54	1.05	0.19	15.31	0.051	0.1	4.25	10.02	300mg
8	2.1	600.4	621.81	0.48	0.77	26.68	0.089	0.09	4.81	5.31	300mg
9	2.7	810.86	809.78	0.37	0.11	31	0.1	0.098	5.98	3.79	300mg
10	2.39	717.26	715.79	0.42	0.12	33.47	0.11	0.081	4.68	5.14	300mg

All the computed parameters depend on the type of observation and administration. Description of all the parameters is [here](#). You can copy this table to a word or an excel file using an icon on the top right (purple frame).

Summary statistics on NCA results

Statistics on non-compartmental parameters are displayed in the summary table in the sub-tab "SUMMARY". Description of all the summary calculation is [here](#).

project_covariates.plx [demo] - PKanalix - 2023R1

Project Settings Export Help

Data Tasks Results Plots

NCA

POINTS FOR λ_z SUMMARY INDIV. ESTIM.

Summary

Statistics on non compartmental analysis results

Records per page: 10 13

Page 1 of 2

Showing 1 to 10 of 25 entries

Elapsed time (seconds): 0.141

	MIN	Q1	MEDIAN	Q3	MAX	MEAN	SD	SE	CV	NTOT	NOBS	NMISS	GEOMEAN	GEOSD	GEOCV	HARMMEAN
AUC _{0-∞} (mg·h·L ⁻¹)	18893.42	42220.7	49741.47	61394.08	86402.34	62049.43	86324.49	7011.02	85.1	63	63	0	52825.99	1.4	43.78	41660.07
AUC _{0-tlast} (mg·h·L ⁻¹)	25669.35	42220.49	49307.8	61349.58	106193.24	61021.74	98427.48	6899.96	81.27	63	63	0	52607.04	1.4	42.51	41452.14
CL/F (L·d ⁻¹)	0.02	0.038	0.038	0.045	0.065	0.038	0.009	0.0011	22.8	63	63	0	0.038	1.24	33.37	0.037
C _{last} (mg·L ⁻¹)	0.1	0.12	0.13	0.14	0.32	0.14	0.17	0.021	83.42	63	63	0	0.12	1.57	47.76	0.14
C _{max} (mg·L ⁻¹)	21.0	27.27	46.42	83.04	245.4	70.42	36.44	4.34	49.01	63	63	0	66.98	2.83	44.4	61.25
K _{el} (d ⁻¹)	0.00089	0.00096	0.0011	0.0013	0.0018	0.0011	0.00094	0.000093	31.89	63	63	0	0.00099	1.47	34.45	0.00079
Intercept	-0.51	2.47	0.35	3.72	4.75	0.1	1.05	0.15	35.81	63	63	0	NaN	NaN	NaN	NaN
Number of points	3	4	5	14	34	10.74	7.8	0.88	78.39	63	63	0	8.26	2.11	56.07	4.38
K required	0.09	0.05	0.09	0.09	0.09	0.09	0.023	0.0028	2.31	63	63	0	0.09	1.02	2.35	0.09
Adjusted R squared	0.89	0.88	0.89	0.89	0.89	0.89	0.024	0.009	2.43	63	63	0	0.89	1.03	2.5	0.89

The summary table can be

- split and filtered according to the categorical and continuous covariates tagged in the data set
- filtered according to the acceptance criteria defined in the [NCA settings](#)

Summary

Statistics on non compartmental analysis results

Records per page: 10 13

Page 1 of 2

Showing 1 to 10 of 25 entries

Elapsed time (seconds): 0.004

	AGE	SEX	MIN	Q1	MEDIAN	Q3	MAX	MEAN	SD	SE	CV	NTOT	NOBS	NMISS
AUC _{0-∞} obs	20-38.5	0	391.472997	577.139521	663.464203	731.876034	805.511991	839.521525	143.40919	58.428206	22.435724	4	4	0
AUC _{0-∞} pred	20-38.5	1	429.032774	740.902488	802.402501	1285.049323	5403.045484	1441.47845	1543.229423	521.074474	108.444257	9	9	0
AUC _{0-∞} obs	38.5-66.5	0	402.777208	539.923232	694.246896	813.085525	957.089739	837.542583	194.124685	80.468384	28.491884	4	4	0
AUC _{0-∞} pred	38.5-66.5	1	886.349804	666.894931	835.844676	949.093241	1140.070335	127.146894	200.352483	40.408547	24.222116	11	11	0
AUC _{0-∞} obs	66.5-85	0	467.747093	605.877438	474.414744	779.494331	1039.796156	707.128564	171.412442	54.205374	24.240433	10	10	0
AUC _{0-∞} pred	66.5-85	1	607.144244	407.443564	688.135459	1028.417599	1163.190774	797.092017	238.759314	97.473082	28.956172	4	4	0
C _{max}	20-38.5	0	32.36	45.82	55.12	59.07	93.6	54.891467	20.409878	8.332286	35.891295	4	4	0
C _{max}	20-38.5	1	21.7	34.77	64.28	89.955	237.7	83.443333	62.740358	20.920119	75.213148	9	9	0
C _{max}	38.5-66.5	0	58.64	59.78	72.43	85.72	87.4	72.544467	14.748978	6.021082	20.324177	4	4	0
C _{max}	38.5-66.5	1	44.22	67.935	80.15	88.6275	100.72	74.930909	18.408218	5.410989	24.344487	11	11	0
C _{max}	66.5-85	0	59.64	45.48	76.54	78.39	94.93	74.944	10.707247	3.389329	14.283185	10	10	0
C _{max}	66.5-85	1	59.05	42.39	70.995	94.40	94.65	70.421667	15.793252	6.447168	20.938941	4	4	0
Flag_AUC_%Extrap_pred	20-38.5	0	1	1	1	1	1	1	1	0	0	4	4	0
Flag_AUC_%Extrap_pred	20-38.5	1	1	1	1	1	1	1	1	0	0	9	9	0
Flag_AUC_%Extrap_pred	38.5-66.5	0	1	1	1	1	1	1	1	0	0	4	4	0

Stratification Groups

AGE

HT

WT

SEX

Split Filter

% extrapolated AUC

Adjusted R2

Span

20 x 38.5

AGE

38.5 x 66.5

66.5 x 85

HT

WT

SEX

- Split by covariates:** When you apply "split" of a table, the values of the splitting covariates are displayed in the first columns of the summary table (blue highlight). The order of these columns corresponds to the order of the clicks to setup the splitting covariates (orange highlight). You can discretized continuous covariates into groups by defining the group limits (yellow highlight), and create new categories for categorical covariates (drug and drop to merge existing categories).

It is currently not possible to split the table in several subtables (instead of splitting the rows), nor to choose the orientation of the table (NCA parameters as columns for instance). See the [reporting documentation](#) for details how these options can be applied in a report.

- **Filter by acceptance criteria:** you can remove from the calculations in the summary table parameter values from subjects that do not meet the acceptance criteria

[Read details on filtering rules](#)

When you save a PKanalix project, the stratification settings of the results are saved in the result folder and are reloaded when reloading a PKanalix project. **The table in the <result folder>/PKanalix/IndividualParameters/nca/ncaIndividualParametersSummary.txt takes into account the split definition.** This table is generated when clicking “run” in the task tab (usually without splits, as the result tab is not yet available to define them) and also upon saving the project (with splits if defined).

Points for lambda_z

In the sub-tab “POINTS FOR λ_z ” you find table of concentration points used for the terminal slope (λ_z) calculation.

ID	time	concentration	BLQ	included for λ_z
1	0	0	0	0
1	0.13	2.48	0	0
1	1	14.93	0	0
1	3	19.72	0	0
1	5	17.61	0	1
1	7	14.21	0	1
1	9	13.01	0	1
1	14	9.86	0	1
1	21	5.66	0	1
1	28	4.08	0	1

This table is displayed only if number of rows <5000. Above this, the tab appears as usual but with incomplete data. You can find the full table in the result folder: <result folder>/PKanalix/IndividualParameters/nca/pointsIncludedForLambdaZ.txt.

NCA plots

The PLOTS tab shows plot associated to the individual parameters:

- **Correlation between NCA parameters:** displays scatter plots for each pair of parameters. It allows to identify correlations between parameters, which can be used to see the results of your analysis and see the coherence of the parameters for each individuals.
- **Distribution of the NCA parameters:** shows the empirical distribution of the parameters. it allows to verify distribution of parameters over the individuals.
- **NCA parameters w.r.t. covariates:** displays the individual parameters as a function of the covariates. It allows to identify correlation effects between the individual parameters and the covariates.
- **NCA individual fits:** shows the lambdaZ regression line for each individual.

NCA outputs

After running the NCA task, the following files are available in the result folder: <resultFolder>/PKanalix/IndividualParameters/nca

- **Summary.txt** contains the summary of the **NCA parameters** calculation, in a format easily readable by a human (but not easy to parse for a computer)
- **ncaIndividualParametersSummary.txt** contains the summary of the **NCA parameters** in a computer-friendly format.
 - The first column corresponds to the name of the parameters

- The second column corresponds to the CDISC name of the parameters
- The other columns correspond to the several elements describing the summary of the parameters (as explained [here](#))
- **ncaIndividualParameters.txt** contains the **NCA parameters** for each subject-occasion along with the covariates.
 - The first line corresponds to the name of the parameters
 - The second line corresponds to the CDISC name of the parameters
 - The other lines correspond to the value of the parameters
- **pointsIncludedForLambdaZ.txt** contains for each individual the concentration points used for the lambda_z calculation.
 - id: individual identifiers
 - occ: occasions (if present). The column header corresponds the data set header of the column(s) tagged as occasion(s).
 - time: time of the measurements
 - concentration: concentration measurements as displayed in the NCA individual fits plot (i.e after taking into the BLQ rules)
 - BLQ: if the data point is a BLQ (1) or not (0)
 - includedForLambdaZ: if this data point has been used to calculate the lambda_Z (1) or not (0)

The files **ncaIndividualParametersSummary.txt** and **ncaIndividualParameters.txt** can be exported in R for example using the following command

```
read.table("/path/to/file.txt", sep = ",", header = T)
```

Remarks

- To load the individual parameters using PKanalix name as headers, your just need to skip the second line

```
ncaParameters = read.table("/path/to/file.txt", sep = ",", header = T);  
ncaParameters[-1,] # to remove the CDISC name line
```

- To load the individual parameters using CDISC as headers, your just need to skip the second line

```
ncaParameters = read.table("/path/to/file.txt", sep = ",", header = T, skip = 1)
```

- The separator is the one defined in the user preferences. We set “,” in this example as it is the one by default.


3.1. Check lambda_z

The λ_z is the slope of the terminal elimination phase. If estimated, the **NCA parameters** will be extrapolated to infinity. Check lambda_z tab visualizes the linear regression model for each individual and allows to control the points included in the calculation of the λ_z parameter.

General rule

The “Check lambda_z” tab shows the data of each individual plotted by default on y-axis log-scale.

- Points used in the λ_z calculation are in blue.
- Points not used in the λ_z calculation are in grey.
- The λ_z curve is in green.



A main rule to select which points to include in the λ_z calculation applies to all individuals. You find it on top in the right panel – red frame in the figure above. There are four methods in a drop-down menu:

- “R2” chooses the number of points that optimize the correlation coefficient.

- “Adjusted R2” (called “Best Fit” in Winonlin) optimizes the adjusted correlation coefficient, taking into account the number of points. The “maximum number of points” and “minimum time” allow to further restrict the points selection.
- “Interval” selects all points within a given time interval.
- “Points” corresponds to the number of terminal points taken into account in the calculation.

The calculation of the λ_z is via the linear regression. The selection of the weighting used in the cost function is in the right panel (blue highlight). Three options are available: uniform (default), $1/Y$, $1/Y^2$.


See the [Settings](#) page and the [Calculations](#) page for more details. In addition to the main rule, it is possible to define more precisely (manually) which points to include or exclude for each individual.

Manual choice of measurements for the λ_z calculation


Charts in the check lambda_z tab are interactive. For each individual specific points can be selected manually to be included or excluded from the λ_z calculation. There are two methods to do it: define a range (time interval) with a slider or select each point individually.

Change of the range


For a given individual, the grey arrow below each individual chart defines a range – time interval when measured concentrations are used in the lambda-z calculation. New points to include are colored in light blue, and points to remove are colored in dark grey, as in the following figure:



The button APPLY INCLUDE/EXCLUDE in the top – right corner of the chart area applies the manual modifications of the points in the calculation of the λ_z . If modifications concern several individuals, then clicking on the “Apply include/exclude” button updates all of them automatically.



The plots are updated with new points and new λ_z curve. For modified plots, the “return arrow” appears in the top – left corner, red frame in the figure below. It resets the individual calculation to the general rule. The drop-down menu in the “Apply include/exclude” (blue frame) allows to reset all individuals to the general rule. Clear selection the current selection made before clicking on the “Apply include/exclude” button.




Include/exclude single points

It is possible to manually modify each point to include it in the calculation of λ_z or to remove it from the calculation:


- If you click on an unused point (light grey), it will change its color into light blue and will be selected for inclusion (figure on the left below).
- If you click on a used point (blue), it will change its color into dark grey and will be selected for exclusion (figure on the right below).

Clicking on the "Apply include/exclude" button applies these changes. The "return arrow" in the top - left corner of each individual plot resets the individual selection to the general rule one by one, while the "Reset to general rule" button in the "Apply include/exclude" drop-down menu resets all of them at once.





Remarks

- If the general rule is modified, only the individuals for which the general rule applies will be updated. Individual for which the points have been manually modified will not be updated.
- If the λ_z can not be computed, the background of the plot will be in red as in the following figure.



3.2. Data processing and calculation rules

This page presents the rules applied to pre-process the data set and the calculation rules applied for the NCA analysis.

- [Data processing](#)
 - [Ignored data](#)
 - [Forbidden situations](#)
 - [Additional points](#)
 - [BLQ data](#)
 - [Steady-state data](#)
 - [Urine data](#)
- [Calculation rules](#)
 - [Lambda_z](#) 
 - [AUC calculation](#) 
 - [Interpolation for partial AUC](#)

Data processing

Ignored data

All observation points occurring before the last dose recorded for each individual are excluded. Observation points occurring at the same time as the last dose are kept, irrespective of their position in the data set file.

Note that for plasma data, negative or zero concentrations are *not* excluded.

Forbidden situations

For plasma data, mandatory columns are ID, TIME, OBSERVATION, and AMOUNT. For urine data, mandatory columns are ID, TIME, OBSERVATION, AMOUNT and one REGRESSOR (to define the volume).

Two observations at the same time point will generate an error.

For urine data, negative or null volumes and negative observations generate an error.

Additional points

For plasma data, if an individual has no observation at dose time, a value is added:

- **Extravascular and Infusion data:** For single dose data, a concentration of zero. For steady-state, the minimum value observed during the dosing interval.
- **IV Bolus data:** the concentration at dose time (C0) is extrapolated using a log-linear regression (i.e log(concentration) versus time) with uniform weight of first two data points. In the following cases, C0 is taken to be the first observed measurement instead (can be zero or negative):
 - one of the two observations is zero
 - the regression yields a slope ≥ 0

BLQ data

Measurements marked as BLQ data with a "1" in the CENSORING column will be replaced by zero, the LOQ value or the LOQ value divided by 2, or considered as missing (i.e excluded) depending on the setting chosen. They are then handled as any other measurement. The LOQ value is indicated in the OBSERVATION column of the data set.

Steady-state

Steady-state is indicated using the [STEADY-STATE](#) and [INTERDOSE INTERVAL](#) column-types. Equal dosing intervals are assumed. Observation points occurring after the dose time + interdose interval are excluded for Cmin and Cmax, but not for lambda_z. Dedicated parameters are computed such as the AUC in the interdose interval, and some specific formula should be considered for the clearance and the volume for example. More details can be found [here](#).

Urine

Urine data is assumed to be single-dose, irrespective of the presence of a STEADY-STATE column. For the NCA analysis, the data is not used directly. Instead the intervals midpoints and the excretion rate for each interval (amount eliminated per unit of time) are calculated and used:

$$\text{midpoint} = \frac{\text{start time} + \text{end time}}{2}$$
$$\text{excretion rate} = \frac{\text{concentration} \times \text{volume}}{\text{end time} - \text{start time}}$$

Calculation rules

Lambda_z

PKanalix tries to estimate the slope of the terminal elimination phase, called λ_z , as well as the intercept called Lambda_z_intercept. λ_z is calculated via a linear regression between $Y=\log(\text{concentrations})$ and the $X=\text{time}$. Several weightings are available for the regression: uniform, $1/Y$ and $1/Y^2$.

Zero and negative concentrations are excluded from the regression (but not from the NCA parameter calculations). The number of points included in the linear regression can be chosen via the "[Main rule](#)" setting. In addition, the user can define specific points to include or exclude for

each individual (see [Check lambda_z page](#) for details). When one of the automatic “main rules” is used, points prior to Cmax, and the point at Cmax for non-bolus models are never included. Those points can however be included manually by the user. If λ_z can be estimated, [NCA parameters](#) will be extrapolated to infinity.

R2 rule: the regression is done with the three last points, then four last points, then five last points, etc. If the R2 for n points is larger than or equal to the R2 for (n-1) points – 0.0001, then the R2 value for n points is used. Additional constraints of the measurements included in the λ_z calculation can be set using the “maximum number of points” and “minimum time” settings. If strictly less than 3 points are available for the regression or if the calculated slope is positive, the λ_z calculation fails.

Adjusted R2 rule: the regression is done with the three last points, then four last points, then five last points, etc. For each regression the adjusted R2 is calculated as:

$$\text{Adjusted R2} = 1 - \frac{(1-R^2) \times (n-1)}{(n-2)}$$

with (n) the number of data points included and (R^2) the square of the correlation coefficient.

If the adjusted R2 for n points is larger than or equal to the adjusted R2 for (n-1) points – 0.0001, then the adjusted R2 value for n points is used. Additional constraints of the measurements included in the λ_z calculation can be set using the “maximum number of points” and “minimum time” settings. If strictly less than 3 points are available for the regression or if the calculated slope is positive, the λ_z calculation fails.

Interval: strictly positive concentrations within the given time interval are used to calculate λ_z . Points on the interval bounds are included. Semi-open intervals can be defined using +/- infinity.

Points: the n last points are used to calculate λ_z . Negative and zero concentrations are excluded *after* the selection of the n last points. As a consequence, some individuals may have less than n points used.

AUC calculation

The following linear and logarithmic rule apply to calculate the AUC and AUMC over an interval [t1, t2] where the measured concentrations are C1 and C2. The total AUC is the sum of the AUC calculated on each interval. If the logarithmic AUC rule fails in an interval because C1 or C2 are null or negative, then the linear interpolation rule will apply for that interval.

Linear formula:

$$\begin{aligned} AUC|_{t_1}^{t_2} &= (t_2 - t_1) \times \frac{C_1 + C_2}{2} \\ AUMC|_{t_1}^{t_2} &= (t_2 - t_1) \times \frac{t_1 \times C_1 + t_2 \times C_2}{2} \end{aligned}$$

Logarithmic formula:

$$\begin{aligned} AUC|_{t_1}^{t_2} &= (t_2 - t_1) \times \frac{C_2 - C_1}{\ln(\frac{C_2}{C_1})} \\ AUMC|_{t_1}^{t_2} &= (t_2 - t_1) \times \frac{t_2 \times C_2 - t_1 \times C_1}{\ln(\frac{C_2}{C_1})} - (t_2 - t_1)^2 \times \frac{C_2 - C_1}{\ln(\frac{C_2}{C_1})^2} \end{aligned}$$

Interpolation formula for partial AUC

When a partial AUC is requested at time points not included in the original data set, it is necessary to add an additional measurement point. Those additional time points can be before or after the last observed data point.

Note that the partial AUC is not computed if a bound of the interval falls before the dosing time.

Additional point before last observed data point

Depending on the choice of the “Integral method” setting, this can be done using a linear or log formula to find the added concentration C* at requested time t*, given that the previous and following measurements are C1 at t1 and C2 at t2.

Linear interpolation formula:

$$C^* = C_1 + \left| \frac{t^* - t_1}{t_2 - t_1} \right| \times (C_2 - C_1)$$

Logarithmic interpolation formula:

$$C^* = \exp\left(\ln(C_1) + \left| \frac{t^* - t_1}{t_2 - t_1} \right| \times (\ln(C_2) - \ln(C_1))\right)$$

If the logarithmic interpolation rule fails in an interval because C1 or C2 are null or negative, then the linear interpolation rule will apply for that interval.

Additional point after last observed data point

If λ_z is not estimable, the partial area will not be calculated. Otherwise, λ_z is used to calculate the additional concentration C*:

$$C^* = \exp(\text{Lambda_z_intercept} - \lambda_z \times t)$$

3.3. NCA parameters

The following page describes all the parameters computed by the non compartmental analysis. Parameter names are fixed and cannot be changed.

- Parameters related to λ_z
- Parameters related to plasma/blood measurements
- Parameters related to plasma/blood measurements specific to steady state dosing regimen
- Parameters related to urine measurements
- Parameter renamings

Parameters related to λ_z

Name	PKPARAMCD CDISC	PKPARAM CDISC	UNITS	DESCRIPTION
Rsq	R2	R Squared	no unit	Goodness of fit statistic for the terminal (log-linear) phase between the linear regression and the data
Rsq_adjusted	R2ADJ	R Squared Adjusted	no unit	Goodness of fit statistic for the terminal elimination phase, adjusted for the number of points used in the estimation of λ_z
Corr_XY	CORRXY	Correlation Between TimeX and Log ConcY	no unit	Correlation between time (X) and log concentration (Y) for the points used in the estimation of λ_z
No_points_lambda_z	LAMZNPT	Number of Points for Lambda z	no unit	Number of points considered in the λ_z regression
Lambda_z	LAMZ	Lambda z	1/time	First order rate constant associated with the terminal (log-linear) portion of the curve. Estimated by linear regression of time vs. log concentration
Lambda_z_lower	LAMZLL	Lambda z Lower Limit	time	Lower limit on time for values to be included in the λ_z calculation
Lambda_z_upper	LAMZUL	Lambda z Upper Limit	time	Upperlimit on time for values to be included in the λ_z calculation
HL_Lambda_z	LAMZHL	Half-Life Lambda z	time	Terminal half-life = $\ln(2)/\text{Lambda_z}$
Lambda_z_intercept	-	-	no unit	Intercept found during the regression for (λ_z), i.e. value of the regression (in log-scale) at time 0, i.e. the regression writes $\log(\text{Concentration}) = -\text{Lambda_z} * t + \text{Lambda_z_intercept}$
Span	-	-	no unit	Ratio between the sampling interval of the measurements used for the λ_z and the terminal half-life = $(\text{Lambda_z_upper} - \text{Lambda_z_lower}) * \text{Lambda_z} / \ln(2)$

Parameters related to plasma/blood measurements

Name	PKPARAMCD CDISC	PKPARAM CDISC	UNITS	DESCRIPTION
Tlag	TLAG	Time Until First Nonzero Conc	time	Tlag is the time prior to the first measurable (non-zero) concentration. Tlag is 0 if the first observation after the last dose is not 0 or LOQ. The value is set to 0 for non extravascular input.
T0			time	Time of the dose
Dose			amount	Amount of the dose
N_Samples	-	-	no unit	Number of samples in the individuals.
C0	C0	Initial Conc	amount/volume	If a PK profile does not contain an observation at dose time (C0), the following value is added <i>Extravascular and Infusion data</i> . For single dose data, a concentration of zero. For steady-state, the minimum observed during the dose interval. <i>IV Bolus data</i> . Log-linear regression of first two data points to back-extrapolate C0.
Tmax	TMAX	Time of CMAX	time	Time of maximum observed concentration. - For non-steady-state data, the entire curve is

				considered. – For steady-state data, Tmax corresponds to points collected during a dosing interval. If the maximum observed concentration is not unique, then the first maximum is used.
Cmax	CMAx	Max Conc	amount/volume	Maximum observed concentration, occurring at Tmax. If not unique, then the first maximum is used.
Cmax_D	CMAxD	Max Conc Norm by Dose	1/volume	Maximum observed concentration divided by dose. $Cmax_D = Cmax/Dose$
Tlast	TLST	Time of Last Nonzero Conc	time	Last time point with measurable concentration
Clast	CLST	Last Nonzero Conc	amount/volume	Concentration of last time point with measurable concentration
AUClast	AUCLST	AUC to Last Nonzero Conc	time.amount/volume	Area under the curve from the time of dosing to the last measurable positive concentration. The calculation depends on the Integral method setting.
AUClast_D	AUCLSTD	AUC to Last Nonzero Conc Norm by Dose	time/volume	Area under the curve from the time of dosing to the last measurable concentration divided by the dose. $AUClast_D = AUClast/Dose$
AUMClast	AUMCLST	AUMC to Last Nonzero Conc	time ² .amount/volume	Area under the moment curve (area under a plot of the product of concentration and time versus time) from the time of dosing to the last measurable concentration.
MRTlast	MRTIVLST	MRT Intravasc to Last Nonzero Conc	time	[if intravascular] Mean residence time from the time of dosing to the time of the last measurable concentration, for a substance administered by intravascular dosing. $MRTlast_iv = AUMClast/AUClast - TI/2$, where TI represents infusion duration.
MRTlast	MRTEVLST	MRT Extravasc to Last Nonzero Conc	time	[if extravascular] Mean residence time from the time of dosing to the time of the last measurable concentration for a substance administered by extravascular dosing. $MRTlast_ev = AUMClast/AUClast - TI/2$, where TI represents infusion duration.
AUCall	AUCALL	AUC All	time.amount/volume	Area under the curve from the time of dosing to the time of the last observation. If the last concentration is positive AUClast=AUCall. Otherwise, AUCall will not be equal to AUClast as it includes the additional area from the last measurable concentration down to zero or negative observations.
AUCINF_obs	AUCIFO	AUC Infinity Obs	time.amount/volume	AUC from Dosing_time extrapolated to infinity, based on the last observed concentration. $AUCINF_obs = AUClast + Clast/Lambda_z$
AUCINF_D_obs	AUCIFOD	AUC Infinity Obs Norm by Dose	time/volume	AUCINF_obs divided by dose $AUCINF_D_obs = AUCINF_obs/Dose$
AUC_PerCentExtrap_obs	AUCPEO	AUC %Extrapolation Obs	%	Percentage of AUCINF_obs due to extrapolation from Tlast to infinity. $AUC_ \%Extrap_obs = 100*(1 - AUClast / AUCINF_obs)$
AUC_PerCentBack_Ext_obs	AUCPBEO	AUC %Back Extrapolation Obs	%	Applies only for intravascular bolus dosing. the percentage of AUCINF_obs that was due to back extrapolation to estimate C(0).
AUMCINF_obs	AUMCIFO	AUMC Infinity Obs	time ² .amount/volume	Area under the first moment curve to infinity using observed Clast $AUMCINF_obs = AUMClast + (Clast/Lambda_z)*(Tlast + 1.0/Lambda_z)$

AUMC_PerCentExtrap_obs	AUMCPEO	AUMC % Extrapolation Obs	%	Extrapolated (% or total) area under the first moment curve to infinity using observed Clast $AUMC_ \%Extrap_obs = 100*(1 - AUMClast / AUMCINF_obs)$
MRTINF_obs	MRTIVFO	MRT Intravasc Infinity Obs	time	[if intravascular] Mean Residence Time extrapolated to infinity for a substance administered by intravascular dosing using observed Clast $MRTINF_obs_iv = AUMCINF_obs/AUCINF_obs - TI/2$, where TI represents infusion duration.
MRTINF_obs	MRTEVIFO	MRT Extravasc Infinity Obs	time	[if extravascular] Mean Residence Time extrapolated to infinity for a substance administered by extravascular dosing using observed Clast $MRTINF_obs_ev = AUMCINF_obs/AUCINF_obs$
Vz_F_obs	VZFO	Vz Obs by F	volume	[if extravascular] Volume of distribution associated with the terminal phase divided by F (bioavailability) $Vz_F_obs = Dose/Lambda_z/AUCINF_obs$
Cl_F_obs	CLFO	Total CL Obs by F	volume/time	[if extravascular] Clearance over F (based on observed Clast) $Cl_F_obs = Dose/AUCINF_obs$
Vz_obs	VZO	Vz Obs	volume	[if intravascular] Volume of distribution associated with the terminal phase $Vz_obs = Dose/Lambda_z/AUCINF_obs$
Cl_obs	CLO	Total CL Obs	volume/time	[if intravascular] Clearance (based on observed Clast) $Cl_obs = Dose/AUCINF_obs$
Vss_obs	VSSO	Vol Dist Steady State Obs	volume	[if intravascular] An estimate of the volume of distribution at steady state based last observed concentration. $Vss_obs = MRTINF_obs*Cl_obs$
Clast_pred	-	-	amount/volume	$Clast_pred = exp(Lambda_z_intercept - Lambda_z* Tlast)$ The values alpha (corresponding to the y-intercept obtained when calculating λ_z) and lambda_z are those values found during the regression for λ_z
AUCINF_pred	AUCIFP	AUC Infinity Pred	time.amount/volume	Area under the curve from the dose time extrapolated to infinity, based on the last predicted concentration, i.e., concentration at the final observation time estimated using the linear regression performed to estimate λ_z . $AUCINF_pred = AUClast + Clast_pred/Lambda_z$
AUCINF_D_pred	AUCIFPD	AUC Infinity Pred Norm by Dose	time/volume	AUCINF_pred divided by dose $= AUCINF_pred/Dose$
AUC_PerCentExtrap_pred	AUCPEP	AUC %Extrapolation Pred	%	Percentage of AUCINF_pred due to extrapolation from Tlast to infinity $AUC_ \%Extrap_pred = 100*(1 - AUClast / AUCINF_pred)$
AUC_PerCentBack_Ext_pred	AUCPBEP	AUC %Back Extrapolation Pred	%	Applies only for intravascular bolus dosing. The percentage of AUCINF_pred that was due to back extrapolation to estimate C(0).
AUMCINF_pred	AUMCIFP	AUMC Infinity Pred	time2.amount/volume	Area under the first moment curve to infinity using predicted Clast $AUMCINF_pred = AUMClast + (Clast_pred/Lambda_z)*(Tlast+1/Lambda_z)$
AUMC_PerCentExtrap_pred	AUMCPEP	AUMC % Extrapolation Pred	%	Extrapolated (% or total) area under the first moment curve to infinity using predicted Clast $AUMC_ \%Extrap_pred = 100*(1 - AUMClast / AUMCINF_pred)$
MRTINF_pred	MRTIVIFP	MRT Intravasc Infinity Pred	time	[if intravascular] Mean Residence Time extrapolated to infinity for a substance

				administered by intravascular dosing using predicted Clast
MRTINF_pred	MRTEVIFP	MRT Extravasc Infinity Pred	time	[if extravascular] Mean Residence Time extrapolated to infinity for a substance administered by extravascular dosing using predicted Clast
Vz_F_pred	VZFP	Vz Pred by F	volume	[if extravascular] Volume of distribution associated with the terminal phase divided by F (bioavailability) $= \text{Dose}/\text{Lambda}_z/\text{AUCINF_pred}$
Cl_F_pred	CLFP	Total CL Pred by F	volume/time	[if extravascular] Clearance over F (using predicted Clast) $\text{Cl_F_pred} = \text{Dose}/\text{AUCINF_pred}$
Vz_pred	VZP	Vz Pred	volume	[if intravascular] Volume of distribution associated with the terminal phas $\text{Vz_pred} = \text{Dose}/\text{Lambda}_z/\text{AUCINF_pred}$
Cl_pred	CLP	Total CL Pred	volume/time	[if intravascular] Clearance (using predicted Clast) $= \text{Dose}/\text{AUCINF_pred}$
Vss_pred	VSSP	Vol Dist Steady State Pred	volume	[if intravascular] An estimate of the volume of distribution at steady state based on the last predicted concentration. $\text{Vss_pred} = \text{MRTINF_pred} * \text{Cl_pred}$
AUC_lower_upper	AUCINT	AUC from T1 to T2	time.amount/volume	AUC from T1 to T2 (partial AUC)
AUC_lower_upper_D	AUCINTD	AUC from T1 to T2 Norm by Dose	time/volume	AUC from T1 to T2 (partial AUC) divided by Dose
CAVG_lower_upper	CAVGINT	Average Conc from T1 to T2	amount/volume	Average concentration from T1 to T2

Parameters related to plasma/blood measurements specific to steady state dosing regimen

In the case of repeated doses, dedicated parameters are used to define the steady state parameters and some specific formula should be considered for the clearance and the volume for example. Notice that all the calculation dedicated to the area under the first moment curve (AUMC) are not relevant.

Name	PKPARMCD CDISC	PKPARM CDISC	UNITS	DESCRIPTION
Tau	-		time	The (assumed equal) dosing interval for steady-state data.
Ctau	CTAU	Conc Trough	amount/volume	Concentration at end of dosing interval. <i>If the observed concentration does not exist, the value is interpolated. It it cannot be interpolated, it is extrapolated using lambda_z. If lambda_z has not been computed, it is extrapolated as the last observed value.</i>
Ctrough	CTROUGH	Conc Trough	amount/volume	Concentration at end of dosing interval. <i>If the observed concentration does not exist, the value is NaN.</i>
AUC_TAU	AUCTAU	AUC Over Dosing Interval	time.amount/volume	The area under the curve (AUC) for the defined interval between doses (TAU). The calculation depends on the Integral method setting.
AUC_TAU_D	AUCTAUD	AUC Over Dosing Interval Norm by Dose	time/volume	The area under the curve (AUC) for the defined interval between doses (TAU) divided by the dose. $\text{AUC_TAU_D} = \text{AUC_TAU}/\text{Dose}$
AUC_TAU_PerCentExtrap	-	-	%	Percentage of AUC due to extrapolation in steady state. $\text{AUC_TAU_}\% \text{Extrap} = 100 * (\text{AUC} [\text{Tlast}, \text{tau}] \text{ if } \text{Tlast} \leq \text{tau}) / \text{AUC_TAU};$

AUMC_TAU	AUMCTAU	AUMC Over Dosing Interval	time2.amount/volume	The area under the first moment curve (AUMC) for the defined interval between doses (TAU).
Vz_F	VZFTAU	Vz for Dose Int by F	volume	[if extravascular] The volume of distribution associated with the terminal slope following extravascular administration divided by the fraction of dose absorbed, calculated using AUC_TAU. $Vz_F = Dose/Lambda_z/AUC_TAU$
Vz	VZTAU	Vz for Dose Int	volume	[if intravascular] The volume of distribution associated with the terminal slope following intravascular administration, calculated using AUC_TAU. $Vz = Dose/Lambda_z/AUC_TAU$
CLss_F	CLFTAU	Total CL by F for Dose Int	volume/time	[if extravascular] The total body clearance for extravascular administration divided by the fraction of dose absorbed, calculated using AUC_TAU . $CLss_F = Dose/AUC_TAU$
CLss	CLTAU	Total CL for Dose Int	volume/time	[if intravascular] The total body clearance for intravascular administration, calculated using AUC_TAU. $CLss = Dose/AUC_TAU$
Cavg	CAVG	Average Concentration	amount/volume	AUCTAU divided by Tau. $Cavg = AUC_TAU /Tau$
FluctuationPerCent	FLUCP	Fluctuation%	%	The difference between Cmin and Cmax standardized to Cavg, between dose time and Tau. $Fluctuation\% = 100.0 * (Cmax -Cmin)/Cavg$
FluctuationPerCent_Tau	-	-	%	The difference between Ctau and Cmax standardized to Cavg, between dose time and Tau. $Fluctuation\%_Tau = 100.0 * (Cmax -Ctau)/Cavg$
Accumulation_Index	AILAMZ	Accumulation Index using Lambda z	no unit	Theoretical accumulation ratio: Predicted accumulation ratio for area under the curve (AUC) calculated using the Lambda z estimated from single dose data. $Accumulation_Index = 1.0/(1.0 -exp(-Lambda_z *Tau))$
Swing	-	-	no unit	The degree of fluctuation over one dosing interval at steady state $Swing = (Cmax -Cmin)/Cmin$
Swing_Tau	-	-	no unit	$Swing_Tau = (Cmax -Ctau)/Ctau$
Tmin	TMIN	Time of CMIN observation.	time	Time of minimum concentration sampled during a dosing interval.
Cmin	CMIN	Min Conc	amount/volume	Minimum observed concentration between dose time and dose time + Tau.
Cmax	CMAX	Max Conc	amount/volume	Maximum observed concentration between dose time and dose time + Tau.
MRTINF_obs	MRTIVIFO or MRTEVIFO	MRT Intravasc Infinity Obs or MRT Extravasc Infinity Obs	time	Mean Residence Time extrapolated to infinity using predicted Clast, calculated using AUC_TAU.

Parameters related to urine measurements


Name	PKPARMCD CDISC	PKPARM CDISC	UNITS	DESCRIPTION
T0			time	Time of the last administered dose (assumed to be zero unless otherwise specified).
Dose	-	-	amount	Amount of the dose

N_Samples	–	–	no unit	Number of samples in the individuals.
Tlag	TLAG	Time Until First Nonzero Conc	time	Midpoint prior to the first measurable (non-zero) rate.
Tmax_Rate	ERTMAX	Midpoint of Interval of Maximum ER	time	Midpoint of collection interval associated with the maximum observed excretion rate.
Max_Rate	ERMAX	Max Excretion Rate	amount/time	Maximum observed excretion rate.
Mid_Pt_Last	ERTLST	Midpoint of Interval of Last Nonzero ER	time	Midpoint of collection interval associated with Rate_last.
Rate_last	ERLST	Last Meas Excretion Rate	amount/time	Last measurable (positive) rate.
AURC_last	AURCLST	AURC to Last Nonzero Rate	amount	Area under the urinary excretion rate curve from time 0 to the last measurable rate.
AURC_last_D	AURCLSTD	AURC to Last Nonzero Rate Norm by Dose	no unit	The area under the urinary excretion rate curve (AURC) from time zero to the last measurable rate, divided by the dose.
Vol_UR	VOLPK	Sum of Urine Vol	volume	Sum of Urine Volumes
Amount_Recovered			amount	Cumulative amount eliminated.
Percent_Recovered			%	100*Amount_Recovered/Dose
AURC_all	AURCALL	AURC All	amount	Area under the urinary excretion rate curve from time 0 to the last rate. This equals AURC_last if the last rate is measurable.
AURC_INF_obs	AURCIFO	AURC Infinity Obs	amount	Area under the urinary excretion rate curve extrapolated to infinity, based on the last observed excretion rate.
AURC_PerCentExtrap_obs	AURCPEO	AURC % Extrapolation Obs	%	Percent of AURC_INF_obs that is extrapolated
AURC_INF_pred	AURCIFP	AURC Infinity Pred	amount	Area under the urinary excretion rate curve extrapolated to infinity, based on the last predicted excretion rate.
AURC_PerCentExtrap_pred	AURCPEP	AURC % Extrapolation Pred	%	Percent of AURC_INF_pred that is extrapolated
AURC_lower_upper	AURCINT	AURC from T1 to T2 (partial AUC)	amount	The area under the urinary excretion rate curve (AURC) over the interval from T1 to T2.
AURC_lower_upper_D	AURCINTD	AURC from T1 to T2 Norm by Dose	no unit	The area under the urinary excretion rate curve (AURC) over the interval from T1 to T2 divided by Dose
Rate_last_pred	–	–	amount/time	The values alpha and Lambda_z are those values found during the regression for lambda_z

Parameter renamings

Starting from the 2023 version of MonolixSuite, default parameter names used in results tables and plots in the graphical user interface and reports do not match the names present in the Name column of the tables above. The names that are used can be customized in the NCA parameters renamings section of Preferences. The section contains two columns:

- **Parameter** – this column contains parameter names from the tables above,
- **Alias** – this column contains parameter names used in the interface and the values in cells can be changed.



The Alias column is prefilled with default aliases. By clicking on an alias, PKanalix allows users to input characters. There are several options present, when in the input mode:

- Undo – reverse the last action,
- Redo – redo an undone action,
- Subscript x_2 – formats selected text and subsequently entered characters as subscript,
- Superscript x^2 – formats selected text and subsequently entered characters as superscript,
- CF – clears formatting of a selected text,
- Cancel – cancels the changes and restores the alias,
- Accept – accepts the changes (changes can be accepted by clicking anywhere in the interface as well).

3.4. NCA settings

The following page describes all the settings for the parameters calculations.

- [Calculations settings](#)
- [\$\lambda_z\$ calculation settings](#)
- [Parameters to compute](#)
- [Acceptance criteria settings](#)
- [Global: Obs id to use](#)

Calculations

These settings corresponds to the settings impacting the calculation of the NCA parameters.

- **Administration type:** *Intravenous* or *Extravascular*. This defines the drug type of administration. IV bolus and IV infusions must be set as “intravenous”.
 - **Integral method:** Method for AUC and AUMC calculation and interpolation. Trapezoidal refers to formula to calculate the AUC and Interpolation refers to the formula to add an additional point in case of partial AUC with a time point not originally present in the data set. See the [Calculation rules](#) for details on the formulas.
 - *Linear Trapezoidal Linear* (Linear trapezoidal, linear interpolation): For AUC calculation linear formula. For interpolation linear formula.
 - *Linear Log Trapezoidal* (Linear/log trapezoidal, linear/log interpolation): For AUC, linear before Cmax and log after Cmax. For interpolation, linear before Cmax and log after.
 - *Linear Up Log Down* (LinUpLogDown trapezoidal, LinUpLogDown interpolation): For AUC, linear if concentration is going up or is stable, log if going down. For interpolation, linear if concentration is going up or is stable, log if going down.
 - *Linear Trapezoidal linear/log* (Linear trapezoidal, linear/log interpolation): For AUC, linear formula. For interpolation, linear before Cmax, log after Cmax. If several Cmax, the first one is used.
- For all methods, if an observation value is less than or equal to zero, the program defaults to the linear trapezoidal or interpolation rule for that point. Similarly, if adjacent observations values are equal to each other, the program defaults to the linear trapezoidal or interpolation rule.
- **Partial AUC time:** Define if the user would like to compute a partial AUC in a specific time interval. It is possible to define several partial AUC intervals, applicable to all individuals.
 - **BLQ before Tmax:** In case of BLQ measurements, this corresponds to the method by which the BLQ data before Tmax should be replaced. Possible methods are “missing”, “LOQ”, “LOQ/2” or “0” (default value is 0).
 - **BLQ after Tmax:** In case of BLQ measurements, this corresponds to the method by which the BLQ data after Tmax should be replaced. Possible methods are “missing”, “LOQ”, “LOQ/2” or “0” (default value is LOQ/2).

Settings impacting the λ_z calculation

This settings are impacting the calculation of λ_z . See also the [Check lambda_z page](#).

- **Main rule** for the λ_z estimation. It corresponds to the rule to define the measurements used for the λ_z calculation. Possible rules are "R2", "interval", "points" or "adjustedR2" (called Best Fit in Winonlin) (default value is "adjustedR2"). See the [calculation rules](#) for more details.
 - In case of "R2" or "adjustedR2" rule, the user has the possibility to define the *maximum number of points* and/or *minimum time for the λ_z estimation*. It allows to constrain the points tested for inclusion.
 - In case of "Interval" rule, the user has the possibility to define the *time interval* to consider.
 - In case of "Points" rule, the user has the possibility to define the *number of points* to consider.
- **Weighting method** used for the regression that estimates λ_z . Possible methods are "Y", "Y2" or "uniform" (default value is "uniform").

Selecting NCA parameters to compute

The NCA parameters to compute can be selected from the "Parameter" column to be placed in the "Computed parameters" column, either by clicking on one or several parameters and using the arrows, or by double clicking on the parameters.

It is recommended to always compute Rsq_adjusted, No_points_lambda_z, Lambda_z and Lambda_z_intercept as these parameters are necessary to display the terminal slope and information box in the [NCA fit plot](#).

The list of parameters computed by default when starting a new PKanalix project can be defined in the Settings > Preferences.

Acceptance criteria

These settings corresponds to the settings impacting the acceptance criteria. When acceptance criteria are defined, flags indicating if the condition is met or not are added to the output result table. Statistics in the summary table can be calculated using only individuals who satisfy one or several acceptance criteria. The filtering option is available directly in the Summary sub-tab of the NCA results.

- **Adjusted R2**: It corresponds to the threshold of the adjusted R2 for the estimation of λ_z . If activated, it will fill the value of Flag_Rsq_adjusted: if Rsq_adjusted > threshold, then Flag_Rsq_adjusted=1. The default value is .98.
- **% extrapolated AUC**: It corresponds to the threshold of the percentage of the total predicted AUC (or AUC for urine models) due to the extrapolation to infinity after the last measurement. If activated, it will fill the value of Flag_AUC_%Extrap_pred: if AUC_%Extrap_pred < threshold, then Flag_AUC_%Extrap_pred=1. The default value is 20%.
- **Span**: It corresponds to the threshold of the span. The span corresponds to the ratio of the sampling interval length for the λ_z calculation and the half-life. If activated, it will fill the value of Flag_Span: if Span > threshold, then Flag_Span=1. The default value is 3.

Global

- **Obs. ID to use**: when a column OBSERVATION ID has defined several types of measurements in the data set, this setting permit to choose which one to use for the NCA analysis

3.5. Parameters summary

This pages provides the information on the statistics used for the summary of the individual parameters for both NCA tasks and CA tasks. All the statistics are made with respect to the NOBS individuals where the parameter has value. For example, in the NCA task, the parameter λ_z might not be possible to compute. Thus, the value associated to these subjects will be missing and not taken into account in the statistics.

- **MIN**: Minimum of the parameter over the NOBS individuals
- **Q1**: First quartile of the parameter over the NOBS individuals
- **MEDIAN**: Median of the parameter over the NOBS individuals
- **Q3**: Third quartile of the parameter over the NOBS individuals
- **MAX**: Maximum of the parameter over the NOBS individuals
- **MEAN**: Arithmetic mean of the parameter over the NOBS individuals
- **SD**: Standard deviation of the parameter over the NOBS individuals
- **SE**: Standard error of the parameter over the NOBS individuals
- **CV**: Coefficient of variation of the parameter over the NOBS individuals
- **NTOT**: Total number of individuals
- **NOBS**: Number of individuals with a valid value
- **NMISS**: Number of individuals with no valid value
- **GEOMEAN**: Geometric mean of the parameter over the NOBS individuals
- **GEOSD**: Geometric standard deviation of the parameter over the NOBS individuals
- **GEOCV**: Geometric coefficient of variation over the NOBS individuals
- **HARMMEAN**: Harmonic mean of the parameter over the NOBS individuals

4. Compartmental Analysis

Overview of the workflow

One of the main features of PKanalix is the calculation of the parameters in the compartmental analysis framework. It finds parameters of a compartmental model representing the PK dynamics for each individual using the Nelder-Mead algorithm.


NEW: PKanalix 2023 version performs individual model fit using any model with continuous outputs – loaded from the [built-in library](#) or a [custom model](#).

Compartmental analysis workflow includes:

- [CA task](#)
- [CA model](#)
- [CA check init.](#)
- [CA results](#)
- [CA plots](#)
- [CA outputs](#)
- [CA export to Monolix/Simulx](#)

CA task

To access the compartmental analysis task in PKanalix use a dedicated section (left panel) in the “Tasks” tab, as in the following figure.




This task contains three sub-tabs:

- The first one called “MODEL” allows to select a model, e.g. from a PKanalix models library or your directory, edit it and match the observations from a dataset with model outputs. Complete guide to model definition is in the [CA model page](#).
- The second one called “CHECK INIT.” shows model predictions obtained with the initial parameter values for each individual together with the data points. It helps to initialize model parameters by changing values manually or using the Auto-init function, as explained in the [Initialization page](#).
- The third one called “Run” contains a button to run the calculations and the settings for the model and the calculations. The meaning of all the settings and their default values is defined the [Settings page](#).

CA model


In this section you specify a model. You can edit the selected model directly in the interface. The right hand side section is to match observations from a dataset with model outputs, see [here](#) for a complete guide.



PKanalix accepts custom models – select BROWSE and load a .txt file with the model. There are also built-in libraries (PK, PK double absorption, TMDD, parent – metabolite, etc.), to facilitate the use of the software. Note that models can include multiple outputs. Models used in all MonolixSuite applications have to be in the **mlxtran** language. Note: For PKanalix versions before 2023, only the classical PK models, from the PK model library, are available.

CA check init.

The “Check init” tab helps to initialize model parameters. It shows the model predictions obtained with the initial model parameter values and the individual designs (doses and regressors). Each sub plot shows one individual with overlaid individual data points. This feature is very useful to find “good” initial values, as explained [here](#).



CA results

When you run the CA task, Pkanalix shows automatically the results in the “Results” tab, see the [CA Results](#) page for a complete guide. There are three tables:

project_extravascular.pkl [demo] * - PKanalix - 2023R1

Project Settings Export Help

Data Tasks Results Plots

SUMMARY INDIV. ESTIM. COST

NCA

CA

Individual estimates

Compartmental analysis results per individual

Records per page: 10 20

<< < Page 1 of 2 > >>

Showing 1 to 10 of 20 entries

id	Cost	ka	V	Ct
101	0.32	0.098	0.38	0.6
102	0.48	0.091	0.66	0.94
103	0.22	0.076	0.27	0.58
104	0.5	0.085	0.72	0.72
105	0.23	0.1	0.72	1.16
106	0.24	0.12	1.79	1.31
107	0.37	0.065	0.94	0.6
108	0.43	0.092	0.95	0.86
109	0.43	0.086	0.96	1.12
1010	0.23	0.095	0.44	0.65

- INDIV. ESTIM contains compartmental analysis results for each individual.
- SUMMARY contains statistics on compartmental analysis results and summary calculation (described [here](#)). You can split or filter this table by covariates present in a dataset.
- COST contains cost function information.

All the computed parameters depend on the chosen model. To copy a table and paste in a word or an excel document click on the icon on the top left corner of a table.

CA plots

After running the CA task, several plots associated to the individual parameters are automatically generated and displayed in the “Plots” tab.

- **Individual fits**: displays model fit for each individual on a separate plot.
- **Correlation between CA parameters**: displays scatter plots for each pair of parameters and allows to identify correlations between parameters.
- **Distribution of the CA parameters**: displays the empirical distribution of the parameters and allows to analyse their distribution over the individuals.
- **CA parameters w.r.t. covariates**: displays the individual parameters as a function of the covariates and allows to identify correlation effects between the individual parameters and the covariates.
- **Observations vs predictions** displays observations from a dataset versus model predictions computed using individual parameters and is useful to detect misspecifications in the model.

CA outputs

After running the CA task, the following files are automatically generated in the project Result folder/IndividualParameters/ca :

- **summary.txt** contains the summary of the CA parameters calculation, in a format easily readable by a human (but not easy to parse for a computer).
- **cost.txt** contains the information about total cost and information criteria for a current model.
- **caIndividualParametersSummary.txt** contains the summary of the CA parameters in a computer friendly format useful for post-processing:
 - The first column corresponds to the name of the parameters.
 - The other columns correspond to the several elements describing the summary of the parameters (as explained [here](#)).
- **caIndividualParameters.txt** contains the CA parameters for each subject-occasion along with the covariates.
 - The first line corresponds to the name of the parameters.
 - The other lines correspond to the value of the parameters.


You can export the files **caIndividualParametersSummary.txt** and **caIndividualParameters.txt** in R for example using the following command

```
read.table("/path/to/file.txt", sep = ",", header = T)
```

Remark: the separator is the one defined in the user preferences. We set “,” in this example as it is the one by default.

CA export to Monolix/Simulx

Using the top menu “Export”, a PKanalix project (dataset information, model, estimated parameters) can be exported to Monolix or Simulx for further analysis, see [Export to Monolix/Simulx page](#) for a complete guide.



- **Monolix** for population modeling and further model development and diagnosis.
- **Simulx** to use a CA model and estimated individual parameters for simulations of new scenarios, eg. new dosing regimens.

Note: For PKanalix versions before 2023, only export to Monolix is available.


Next section: **CA model**

4.1. CA model


The first step to run the compartmental analysis task in PKanalix is selecting a model. Models in PKanalix are .txt files and are loaded in the Model sub tab of the Task tab. You can:

- Load a **custom model** from your directory by clicking the BROWSE button.
- Load a **model from PKanalix library** by clicking the LOAD FROM LIBRARY button.
- Write a new model from scratch in the **mlxtran language** by clicking the NEW MODEL button.

Note: In PKanalix versions before 2023, only models from the PK library are available.




After loading a model, PKanalix switches automatically the MODEL sub-tab to the CHECK INIT. sub-tab for **parameters initialization**. To view the model file go back to the MODEL sub-tab. The left section shows the model file content. Clicking on the EDIT MODEL button allows to edit the model, save and apply it to the current project. For more details, see the **custom model** documentation page. The right hand side shows the mapping between the data set observations ids and the model outputs. For more details, see the **mapping in Monolix** documentation page. This feature is the same in PKanalix and Monolix.



Saving model information in the PKanalix project file

In PKanalix, a model file is a .txt file and its name is displayed on top of the MODEL sub-tab. If a model is from the PKanalix library, then the file name contains a prefix "lib:". For custom models there is the full path to the file.


Remark on using custom models. Be careful when you change a location of a model file. If it was used in any PKanalix project, then you will not be able to re-open this project. You will see an error "Missing model file at 'original_path_to_modelFile/model_filename.txt'". To avoid such situation, use a project settings (In the top menu "Settings") and switch on the toggle "Save the dataset and the model beside the project".



Next section: [library of models](#)

4.1.1. Libraries of models


PKanalix provides libraries of ready-to-use models for compartmental analysis. This simplifies the use of the software by decreasing the necessity of coding. You can pick models from the libraries, use them directly or edit. They are accessible in the Model sub-tab of the CA task by clicking the button "LOAD FROM LIBRARY".




Libraries are organized into several categories:

- **PK model library** includes standard compartmental models. It is possible to select different administration routes (bolus, infusion, first-order absorption, zero-order absorption, with or without Tlag), different number of compartments (1, 2 or 3 compartments), and different types of eliminations (linear or Michaelis-Menten). The PK library models can be used with single or multiple doses data, and with two different types of administration in the same data set (oral and bolus for instance).
- **PD model library** includes direct response models such as Emax and Imax with various baseline models, and turnover response models. These models are PD models only and the drug concentration over time must be defined in the data set and passed as a **regressor**.
- **PK/PD model library** provides all standard combinations of pharmacokinetic and pharmacodynamic models.
- **PK double absorption model library** includes double absorption models. They take into account all the combinations of absorption types and delays for two absorptions. The absorptions can be specified as simultaneous or sequential, and with a pre-defined or independent order.
- **Parent-Metabolite model library** includes models describing parent drug and one metabolite, with or without the first pass effect, uni and bidirectional transformation and up to three compartments for parent and metabolite.
- **Target-mediated drug disposition (TMDD) model library** – includes a large number of TMDD models corresponding to different approximations, different administration routes, different parametrizations, and different outputs.
- **Tumor growth inhibition (TGI) library** – includes a wide range of models for tumour growth (TG) and tumour growth inhibition (TGI) that are available in the literature. Models correspond to different hypotheses on the tumor or treatment dynamics.

In each library, available models are displayed as a list. With the top filtering panel you can easily selected a model based on its characteristics such as administration type, number of compartments or elimination mechanism.



To select a model click on its corresponding line – when all filters are used, only one model is shown. The model will be automatically loaded and displayed in the Model tab. Model filename is displayed on top of the "Model file" section, and has a prefix "lib" to distinguish it from custom models.




Next section: [Custom model](#)


4.1.2. Custom model

In the compartmental analysis task in PKanalix you can use any custom model* – any model written in the [Mlxtran language](#). To load an existing model, click on "BROWSE" in the Model sub tab, and to write a new model from scratch click "NEW MODEL" button.

*Note: In PKanalix versions before 2023, only models from the PK library are available.



You can also choose a model from the library and then edit it directly in the Model tab by clicking on the button "EDIT MODEL":



In the edition mode, there are special features that help to write models in mlxtran. It includes functions to color syntax and to detect mlxtran syntax errors, to zoom editing window. In addition, it gives functions and keyword suggestions for the mlxtran language.

You can save a new model and automatically apply it to a current project – click **SAVE & APPLY**. To save a model with a different name, use “**SAVE AS & APPLY**” or “**SAVE AS**”. Note that models from the library cannot be overwritten.

Merging two models from libraries

In MonolixSuite application there is a dedicated application to write and edit models – [mlxEditor](#). Its functionalities are integrated in the interface, as described above, but you can also use [mlxEditor](#) independently of PKanalix. The following video shows an example how to merge two models using [mlxEditor](#).

Writing a model from scratch

The following video describes how to write a custom model. Note that this video uses the old version of [mlxEditor](#), but the functionalities are the same.

Understanding the error messages


The error messages generated when syntax errors occur in the model are very informative and quickly help to get the model right. This video explains the most common error messages..

Next section: [Initialization of model parameters](#)

4.2. CA Task #


The following page describes how compartmental analysis calculates individual parameters of a model selected in the **MODEL** tab. Three elements impact this calculation:

- [Initialization](#) of model parameters and their constrains on their values
- [Cost function](#) used in the optimization algorithm
- [Estimation](#) algorithm



4.2.1. Parameters initialization #

Parameters of a model you selected in the **MODEL** tab are listed in the **RUN** tab > **Settings**, see image below. Each parameter has two attributes: an initial value and a constraint.



Constraint sets a range of allowed values for selected parameter during estimation (and as initial values). There are three types of constraints:

- none: a parameter is not constraint and its value can be from minus to plus "infinity" (in practice its minus or plus 10e16),
- positive: (default option) a parameter is strictly greater then zero,
- bounded: a parameter is strictly inside a bounded interval which limits you can specify manually.


You can set an initial value of a parameter manually or using the auto-initialization algorithm. Selection of "good" initial values – capturing the characteristic features of a selected model (e.g two slopes on log scale for a 2-compartment model) – is important during the estimation of parameters with the optimization algorithm. It can help to avoid local minima during the optimization process and decrease the runtime of the algorithm.

CHECK INIT.

In the Tasks tab of the CA section there is a dedicated sub-tab CHECK INIT. to help you in the initialization of parameters. The goal of the CHECK INIT. sub tab is to:

- Visually check parameters initialization. It displays the model predictions obtained with the initial model parameters values and the individual designs (doses and regressors) for each individual together with the data points.
- Help finding "good" initial values: manually or with the auto-init function.

CHECK INIT. tab in PKanalix uses the same auto-init algorithm and has the same features as CHECK INIT. in Monolix.




- You can modify the initial values of the parameters, shown on the bottom of the screen, manually by typing new values or using the auto-init algorithm – click the button RUN in the right panel.
- Switch toggle buttons in the X-axis and Y-axis sections to set log-scale and to apply the same limits for all individual plots. It gives a better comparison between the individuals.

- If there are not enough points for the prediction (e.g. there are a lot of doses), change the grid size by increasing the number of points (last section in the right panel).

IMPORTANT: Click on the “**SET AS INITIAL VALUES**” button on the top of the plot to save the initial values for the optimization step.

If there are several observation types which obs-ids have been mapped to model outputs (for example a parent and a metabolite, or a PK and a PD observation), then switch the plots between them in the “Output” section on the top-right corner, see figure below.




Reference in the “check initial estimates”

The following video describes how to use the Reference option and when it can be helpful. It uses an example in Monolix, but this feature works in the same way in PKanalix.

[/vc_column_text][vc_column][vc_row]

Adding a reference helps to see better how changing parameter values impacts the prediction. To add a reference click on the icon next to “Set as initial values” button. It saves the current curve (i.e. current parameters values) as a reference and adds it to the list on top of the right panel, see below. The solid red curve corresponds to the current set of parameters (displayed at the bottom), while the dashed one corresponds to the reference. At any time, you can restore the reference as the current curve (green arrow next to the reference curve name), delete the reference (red cross) or delete all references by clicking on the trash icon above the list.



Auto – init: Automatic initialization of the parameters

PKanalix has an automatic algorithm for the initialization. It is in the right panel in the Auto-init section and you can launch it by clicking on the button RUN. The following video describes how it works using an example in Monolix.

<https://youtu.be/IULGI9WqahgI>

After clicking on the button RUN, PKanalix computes initial model parameters that best fit the data points using a pooled fit approach, starting from parameter values currently used in the panel at the bottom. By default, the algorithm uses all individual data from all observations mapped to a model output. You can change the set of individuals in the “Ids” selection panel just below the RUN button.

The algorithm is a custom optimization method, on the pooled data. The purpose is not to find a perfect match but rather to have all the parameters in the good range for starting the Nelder-Mead optimization on each individual.


- While auto-init is running, the pop-up shows the evolution of the cost of the optimization algorithm over the iterations. Stop the algorithm at any time, e.g. if you see that the cost has decreased sufficiently and you want to check the parameter values.
- Note that the more individuals you select, the longer the run will take.
- Selecting one or few individual that show clearly model characteristics (eg a third compartment, a complex absorption) can help the auto-init algorithm to find set of parameters that is sensitive to specific model features.

After running the auto-init, PKanalix interface updates automatically parameter values and plots with new predictions. To use these parameters as initial values in the optimization process, you click on the button “SET AS INITIAL VALUES”.

Note that the auto-init procedure takes into account the current initial values. Sometimes the auto-init might give poor results. To improve it change manually the parameter values before running the auto-init again.

4.2.2. Estimation

The estimation of model parameters is the key function of the compartmental analysis task. In PKanalix it is performed in the framework of generalized least squares estimation. The optimization algorithm consists in minimizing a **cost function** with an improved version of the Nelder-Mead simplex algorithm, starting from multiple guesses in addition to the initial values provided. To run the algorithm click on the button COMPARTMENTAL ANALYSIS in the RUN tab.



The goal of the optimization is to find model parameters that give the best model fit to the observed data – minimize the objective (cost) function.

- The optimization algorithm performs well in general. However, even if it tries to explore the parameter space, the core algorithm (Nelder-Mead) is local. Therefore **proper initialization of parameters** is crucial.
- The cost function is described in [this](#) documentation page.
- In case of **censored data**, you can choose a method to treat BLQ data during the optimization: missing, equal to zero, equal to LOQ or equal to LOQ/2.
- By default, the optimization is performed for each individual independently. Switching on the toggle “**Pooled fit**” in the calculation settings will fit the same parameters to all individuals (data from all individuals is analysed together).
- If occasions are present and pooled fit is not selected, parameters are all considered occasion-dependent so the optimization for every individual includes $N_{param} \times N_{occ}$ parameters (with N_{param} the number of parameters in the structural model and N_{occ} the number of occasions for this individual).

4.2.3. Cost function, AIC and BIC

Cost function and likelihood

Cost function

The cost function that is minimized in the optimization performed during CA is described here. The cost is reported in the [result tables](#) after running CA, as individual cost in the table of individual parameters, and as total cost in a separate table, together with the likelihood and additional criteria.

Formulas for the cost function used in PKanalix are derived in the framework of generalized least squares described in the publication [Banks, H. T., & Joyner, M. L. \(2017\). AIC under the framework of least squares estimation. Applied Mathematics Letters, 74, 33-45](#), with an additional scaling to put the same weight on different observation types.

The total cost function, shown in the COST table in the results, is a sum of individual cost per each occasion (sum over id-occ):

$$Cost = \sum_{id-occs} X_{scaling} \sum_{\substack{obstypes \\ \text{with data for} \\ id-occ}} \frac{1}{N_{obs}(id\#occ, obstype)} \sum_{j \text{ in time}} (res_j / w_j)^2$$

Individual cost shown in case of indiv fit and minimized in indiv fit

where in the above formula:

- $N_{obs}(id\#occ, obsType)$ denotes the number of observations for one individual, one occasion and one observation type (ie one model output mapped to data, eg parent and metabolite, or PK and PD, see [here](#)).
- $res_j = Y_{pred}(t_j) - Y_{obs}(t_j)$ are residuals at measurement times t_j , where Y_{pred} are the predictions for every time point t_j , and Y_{obs} the observations read from the dataset.
- w_j are weights $1, Y_{obs}, Y_{pred}, \sqrt{|Y_{obs}|}, \sqrt{|Y_{pred}|}, \sqrt{|Y_{obs}| \cdot |Y_{pred}|}$ depending on the choice in the CA task > Settings > Calculations, as explained in [this tutorial video](#).

Calculations

Cost function

$(Y_{pred} - Y_{obs})^2 / |Y_{pred} * Y_{obs}|$ ▼

$(Y_{pred} - Y_{obs})^2$

$(Y_{pred} - Y_{obs})^2 / Y_{obs}$

$(Y_{pred} - Y_{obs})^2 / Y_{pred}$

$(Y_{pred} - Y_{obs})^2 / Y_{obs}^2$

$(Y_{pred} - Y_{obs})^2 / Y_{pred}^2$

$(Y_{pred} - Y_{obs})^2 / |Y_{pred} * Y_{obs}|$

with the same parameters for all

$X_{scaling}$ is a scaling parameter given by:

$$X_{scaling} = \frac{\sum_{id-occs} n_{obstypes}(id\#occ)}{\sum_{id-occs} \sum_{obstypes} \left(\frac{1}{N_{obs}(id\#occ, obstype)} \right)}$$

with $n_{obstypes}(id\#occ)$ number of observation IDs mapped to a model output for one individual and one occasion. Note that if there is only one observation type (model output mapped to an observation ID), $X_{scaling} = 1$ and the cost function is a simple sum of squares over all observations, weighted by the w_j .

Likelihood

The likelihood is not directly used for the optimization in CA, but it is possible to derive it based on the optimal cost obtained with the following formula (details in [Banks et al](#)):

$$-2LL = (1 + \ln(2\pi)) N_{tot} + 2 \sum_{j \text{ in data}} \ln(w_j) + N_{tot} \ln\left(\frac{Cost}{N_{tot}}\right)$$

where

- -2LL stands for -2 x ln(Likelihood)
- N_{tot} is the total number of observations for all individuals, all occasions and all mapped model outputs
- w_j are the weights as defined above

Additional Criteria

Additional information criteria are the Akaike Information Criteria (AIC) and the Bayesian Information Criteria (BIC). They are used to compare models. They penalize -2LL by the number of model parameters used for the optimization N_{param} , and the number of observations: N_{tot} , total number of observations, and $N_{obs}(id\#occ)$, number of observations for every subject-occasion. This penalization enables to compare models with different levels of complexity. Indeed, with more parameters, a model can easily give a better fit (lower -2LL), but with too many parameters, a model loses predictive power and confidence in the parameter estimates. If not only the -2LL, but also AIC and BIC criteria decrease after changing a model, it means that the increased likelihood is worth increasing the model complexity.

$$\begin{aligned} \text{Indiv fit} \\ AIC &= -2LL + 2(N_{param}N_{idocc} + 1) \\ BIC &= -2LL + \ln(N_{tot}) + N_{param} \sum_{id-occs} \ln(N_{obs}(id\#occ)) \end{aligned}$$

This formula simplifies if the "pooled fit" option is selected in the calculation settings of the CA task:

$$\begin{aligned} \text{Pooled fit} \\ AIC &= -2LL + 2(N_{param} + 1) \\ BIC &= -2LL + \ln(N_{tot})(N_{param} + 1) \end{aligned}$$

AIC derivation in the case of individual fits

The total AIC for all individuals penalizes -2LL with 2x the number of parameters used in the optimization. For the structural model, we use $N_{param}N_{id\#occ}$, where N_{param} is the number of parameters in the structural model and $N_{id\#occ}$ is the number of subject-occasions (since there is a set of parameters of each id#occ). For the statistical model, we also calculate implicitly one parameter σ^2 , which is the variance of the residual error model shared by all the individuals (it is derived using the sum of squares of residuals as described in [Banks et al](#)). Therefore, AIC penalizes -2LL by $2(N_{param}N_{idocc} + 1)$.

BIC derivation in the case of individual fits

The total BIC should penalize -2LL by (nb of parameters) x ln(nb of data points). In this case, we should distinguish between the data for each id#occ, which is used to optimize N_{param} parameters, and the data for all individuals, which is used all together to calculate the residual variance σ^2 . To penalize each individual's data in the optimization with N_{param} , we add the term $N_{param} \times \ln(N_{obs}(id\#occ))$ for each id#occ. To penalize the calculation of the error model parameter with all the observations (substitution of σ as described in [Banks et al](#)), we add the term $1 \times \ln(N_{tot})$.

4.3. CA results

PKanalix automatically generates results after running the CA task calculations. This page contains a complete list of output tables and files.

- [Output tables](#)
- [Output files](#)

Output tables

Compartmental analysis results are displayed in tables in the "RESULTS" tab, in the CA section.

project_covariates.plx [demo] * - PKanalix - 2023R1

Project Settings Export Help

Data Tasks Results Plots

SUMMARY INDIV. ESTIM. COST

Individual estimates

Compartmental analysis results per individual

id	Cost	ka	V	Cl	WEIGHT	SEX
1	0.82	4.91	5.42	0.44	70.65	F
2	0.68	1.07	3.02	0.5	67.52	F
3	0.37	2.59	20.63	0.91	65.31	M
4	0.32	1.32	9.11	0.97	66.67	M
5	0.22	0.6	4.63	0.61	68.13	F
6	0.53	0.93	5.74	0.36	71.19	F
7	0.64	0.49	5.01	0.54	71.81	F
8	0.25	0.13	2.29	1.75	76.76	M
9	0.27	0.68	3.96	0.36	64.49	F
10	0.29	0.67	4.52	0.32	64.74	F

Records per page: 10 25 50

<< Page 1 of 5 >>

Showing 1 to 10 of 50 entries

INDIV. ESTIM.

In this table (shown above) you find results of the compartmental analysis for each individual: model parameters values, cost function values, covariates (from the dataset). By default, the table shows 10 records (ids) per page. To increase it and to see following pages, use settings below the table.

SUMMARY

This table contains statistics on compartmental analysis results and summary calculation. The complete list with description is [here](#).

project_covariates.plx [demo] * - PKanalix - 2023R1

Project Settings Export Help

Data Tasks Results Plots

SUMMARY INDIV. ESTIM. COST

Summary

Statistics on compartmental analysis results

	MIN	Q1	MEDIAN	Q3	MAX	MEAN	SD	SE	CV	GEOMEAN	GEOSD	GEOCV	HARMMEAN
Cost	0.15	0.27	0.37	0.5	1.18	0.41	0.19	0.027	46.93	0.37	1.53	44.51	0.34
ka	0.13	0.56	0.8	1.32	45016.23	2326.96	9462.37	1338.18	406.64	1.58	15.26	4100.49	0.62
V	2.29	4.63	6.79	9.11	20.63	7.27	3.58	0.51	49.22	6.52	1.6	49.92	5.86
Cl	0.25	0.46	0.81	0.97	2.12	0.8	0.41	0.058	51.53	0.7	1.67	55.09	0.62


Elapsed time (seconds): 0.15

Elapsed time (seconds): 0.15

For better comparison you can split or filter this table using covariates present in a dataset. When you select split/filter PKanalix automatically recalculates the values and displays them. The following Feature of the week video explain how this feature works on the NCA individual estimates table example.

COST

This table contains information about model comparison indices calculated for the model selected in the CA task. It includes: total cost (sum of individual costs values shown in the INDIV. ESTIM. sub tab), -2LL, AIC and BIC, see [here](#) for details about the formulas. This information criteria allow to compare statistically the fit obtained with different models. The rule of thumb is: the lower the value of the information criteria, the better the model fit. There is no absolute value to compare with, important is only the relative difference between two models.



Output Files

You can find the compartmental analysis results also in the project results folder in the "IndividualParameters > ca" directory. There are four .txt files:

summary.txt

Description: human readable summary file.

Outputs:

- **Header:** project file name, date and time of run, PKanalix version
- **Table:** Summary statistics of the individual cost function and individual model parameters

caIndividualParameters.txt

Description: Individual parameters estimates

Outputs:

- **ID:** subject name and occasion (if applicable). If there is one type of occasion, there will be an additional(s) column(s) defining the occasions.
- **Cost:** value of the individual cost function
- **parameterName:** individual model parameter values estimated during the compartmental analysis.
- **COVname:** continuous covariates values corresponding to all data set columns tagged as "Continuous covariate".
- **CATname:** modalities associated to the categorical covariates.

caIndividualParametersSummary.txt

Description: Statistics summary of individual parameters estimates

Outputs:

- **ID:** subject name and occasion (if applicable). If there is one type of occasion, there will be an additional(s) column(s) defining the occasions.
- **Parameter:** cost and model parameters
- **Statistics:** min, Q1, median, Q3, max, mean, SD, SE, CV, geoMean, geoSD, geoCV, harmMean

cost.txt

Description: Total cost and information criteria of a model

Outputs:

- **Cost:** Total cost, see formulas above
- **Log-likelihood estimation:** -2LL, AIC, BIC

4.4. Export to Monolix/Simulx


MonolixSuite applications are interconnected and you can export & import a PKanalix project between different applications. This interconnection is guaranteed by using the same model syntax (mlxtran language) and the same dataset format.

There are two options to export a PKanalix project*

- **Export to Monolix** for population model development, estimation and diagnosis. See [Monolix](#) documentation.
- **Export to Simulx** for simulations of new scenarios, e.g. new dosing regimens. It requires loading a model in CA. See [Simulx](#) documentation.

**Note: Export of a PKanalix project to Simulx is available starting from the 2023 version. In the previous versions, you can only export a PKanalix project only to Monolix.*

To export a PKanalix project to Monolix or Simulx click **EXPORT PROJECT TO** in the top menu "Export", select the target software (Monolix or Simulx), and confirm the process by clicking on the button EXPORT at the bottom. By default, generated files, e.g. dataset and model file, will be saved next to the target project.



Export to Monolix

Exporting a PKanalix project to Monolix opens automatically a Monolix application and creates a new Monolix project with predefined elements:

- Dataset tagging and formatting steps (if present) as in the PKanalix project.
- Dataset filters (if applied) as in the PKanalix project.
- Structural model and initial parameters values as selected in the CA task in the PKanalix project. If there is no model in the PKanalix project, then Monolix will use a default model from the library with oral or iv administration (depending on the NCA settings), one compartment and linear elimination.
- Statistical model and tasks: default Monolix settings

Monolix will open in the "Statistical model and tasks" tab, and you can immediately click the button RUN to execute estimation of parameters and generate diagnostic plots.

Export to Simulx

Exporting a PKanalix project to Simulx creates a Simulx project with all elements that could be useful for simulating the model automatically. In the Definition tab:

- Model, individual parameters estimates and output variables – imported from the CA task in PKanalix.
- Occasions, covariates, treatments and regressors – imported from the dataset used in the PKanalix project.

Moreover, Simulx automatically sets exploration and simulation scenarios. The Exploration (tab) contains one group to explore effects of model parameters and treatment on a single individual. The Simulation (tab) contains one group to re-simulate individuals from the PKanalix project.

You will find more information about imported elements and how to create new scenarios in the [Simulx documentation webpage](#).


5. Bioequivalence

After having run the "NCA" task to calculate the NCA PK parameters, the "Bioequivalence" task can be launched to compare the calculated PK parameters between two groups, usually one receiving the "test" formulation and one the "reference" formulation.

- [Bioequivalence task](#)
- [Bioequivalence results](#)
- [Bioequivalence plots](#)
- [Bioequivalence outputs](#)

Bioequivalence task

The task "Bioequivalence" is available if at least one categorical covariate column has been tagged in the "Data" tab. The "Bioequivalence" task can only be run after having run the "NCA" task. It can be run by clicking on the "Bioequivalence" task button (pink highlight below) or by selecting it as part of the scenario (checkbox in the upper right corner of the "Bioequivalence" button) and clicking the "Run" button. The parameters to be included in the bioequivalence analysis are selected in the box "Parameters to compute / BE" (orange highlight), where the user can also indicate if the parameter should be log-transformed or not. The fixed effect included in the linear model can be chosen in the box "Bioequivalence design" (blue highlight).



Bioequivalence plots

In the "Plots" tab, several plots are displayed:

- **Sequence by period:** This plot allows to visualize, for each parameter, the mean and standard deviation for each period, sequence and formulation.
- **Subject by formulation:** This plot allows to visualize, for each parameter, the subject-by-formulation interaction and the inter-subject variability.
- **Confidence intervals:** This plot gives a visual representation of the confidence intervals of the ratio of each parameter.
- **NCA individual fits:** This plot can also be used to display the individual concentration data with one subplot per individual with all its periods.

Bioequivalence outputs

After running the Bioequivalence task, the following files are available in the result folder <result folder>/PKanalix/IndividualParameters/be:

- **anova_XXX.txt:** these is one such file for each NCA parameter included in the bioequivalence analysis. It contains the ANOVA table with columns 'Factor' (factors included in the linear model), 'Df' (degrees of freedom), 'SumSq' (sum of squares), 'MeanSq' (mean squares), 'FValue' (F-value), 'Pr(>F)' (p-value). For the residuals (last line), the two last columns are empty.
- **confidenceInterval_XXX.txt:** these is one such file per non-reference formulation. It contains the confidence interval table with columns 'Parameter' (parameter name), 'AdjustedMeanTest' (adjusted mean for the test formulation), 'NTest' (number of individuals for the test formulation), 'AdjustedMeanRef' (adjusted mean for the reference formulation), 'NRef' (number of individuals for the ref formulation), 'Difference' (formulation difference – see [calculation rules](#)), 'CIRawLower' (lower confidence interval bound for the difference), 'CIRawUpper' (upper confidence interval bound for the difference), 'Ratio' (formulation ratio – see [calculation rules](#)), 'CILower' (lower confidence interval bound for the ratio), 'CIUpper' (upper confidence interval bound for the ratio), 'Bioequivalence' (1 if the CI for the ratio falls within the BE limits, 0 otherwise)
- **estimatedCoefficients_XXX.txt:** these is one such file for each NCA parameter included in the bioequivalence analysis. It contains the estimated coefficient for each category of each factor included in the linear model, as well as the intercept. This information is only available as an output table and is not displayed in the GUI. The table columns are 'name' (factor followed by the category), 'estimate' (estimated coefficient), 'se' (standard error), 'tValue' (estimate divided by the SE), and 'Pr(>|t|)' (p-value).
- **variationCoef.txt:** This file contains the standard deviation and coefficient of variation from each NCA parameter. Columns are 'Parameter', 'SD' and 'CV(%)'

5.1. Bioequivalence calculation rules

Detection of "crossover" versus "parallel" design

In PKanalix, the design is detected based on the data set columns. If no OCCASION column is present, the design is detected as "parallel". On the opposite, if one or several columns have been tagged as OCCASION, the design is detected as "crossover" (repeated or non-repeated). The detected design is displayed in the bioequivalence settings for the users information and cannot be changed. Depending on the detected design, the factors selected by default for the linear model are different.

Linear model

The linear model can only include fixed effects, as recommended by the FDA for parallel and non-repeated crossover design (see [Statistical Approaches to Establishing Bioequivalence](#), page 10) and by the EMA for parallel, non-repeated crossover and repeated crossover designs (see [Guideline On The Investigation Of Bioequivalence](#), page 15). In addition, "id" is automatically considered as nested in "sequence" and no additional nesting can be defined. Interaction terms and random effects are not supported.

According to the regulatory guidelines, the default models are:

- **parallel design:**

$$\log(\text{NCAparam}) \sim \text{FORM}$$

which can also be written as

$$\log(\text{NCAparam}) = \text{intercept} + \beta_{\text{test}}[\text{ if FORM=test}] + \varepsilon$$

in case of two formulations "ref" and "test", with ε a normal random variable (i.e the residuals).

- **crossover design:**

$$\log(\text{NCAparam}) \sim \text{SEQ} + \text{ID} + \text{PERIOD} + \text{FORM}$$

which can also be written as

$$\log(\text{NCAparam}) = \text{intercept} + \beta_{\text{test}}[\text{ if FORM=test}] + \beta_{\text{period2}}[\text{ if PERIOD=2}] + \beta_{\text{seqTR}}[\text{ if SEQ=TR}] + \beta_{\text{id2}}[\text{ if ID=id2}] + \beta_{\text{id3}}[\text{ if ID=id3}] + \dots$$

in case of several individuals (id1, id2, etc), two sequences ("RT" and "TR"), two periods ("1" and "2") and two formulations ("ref" and "test").

The model parameters are calculated using a QR factorization.

The β parameters are called the coefficients. They are saved in the output file **estimatedCoefficients_XXX.txt**. We are interested in the coefficient representing the formulation effect β_{test} , which is also called the point estimate.

Difference and ratio

In the Results > Confidence intervals table, the "difference" corresponds to the point estimate β_{test} (see above). The "ratio" is calculated depending on the log-transformation choice:

- without log-transformation: $\text{Ratio} = \text{LSM}_{\text{test}} / \text{LSM}_{\text{ref}} \times 100$ with LSM the least square mean (also called adjusted mean, see below).
- with log-transformation: $\text{Ratio} = \exp(\beta_{\text{test}}) \times 100$

Confidence intervals

The confidence interval is first calculated for the difference (CI_{diff}) and then for the ratio (CI_{ratio}).

$$CI_{\text{diff}} = \beta_{\text{test}} \pm t(1 - \alpha, df) \times SE$$

with β_{test} the point estimate ("difference", see above), $t(1 - \alpha, df)$ the quantiles of a Student t-distribution at the α level and df degrees of freedom, and SE the standard error of the point estimate.

In case of a parallel design with the **Welch-Satterthwaite correction** (see [Bioequivalence settings](#)), the formula is different:

$$CI_{\text{diff}} = \beta_{\text{test}} \pm t(1 - \alpha, df) \times \sqrt{\frac{S_R^2}{n_R} + \frac{S_T^2}{n_T}}$$

with S_X the samples standard deviation of the individuals and n_X the number of individuals having received formulation X. A correction is also applied to the degrees of freedom, which are calculated as:

$$df = \frac{\left(\frac{S_T^2}{n_T} + \frac{S_R^2}{n_R} \right)^2}{\frac{(S_T^2/n_T)^2}{n_T - 1} + \frac{(S_R^2/n_R)^2}{n_R - 1}}$$

The confidence interval for the ratio is then calculated using the following formula:

- without log-transformation: $CI_{\text{ratio}} = (1 + CI_{\text{diff}} / \text{LSM}_{\text{ref}}) \times 100$ with LSM_{ref} the least square mean for the reference formulation (see below).
- with log-transformation: $CI_{\text{ratio}} = \exp(CI_{\text{diff}}) \times 100$

Adjusted means (least square means)

The model $\text{NCAparam} \sim \text{SEQ} + \text{ID} + \text{PERIOD} + \text{FORM}$ allows to calculate the expected value (mathematical expectation) for any value of SEQ, ID, PERIOD and FORM. The least square mean for the reference represents the expected value for the reference formulation leaving the value for the other factors undefined. To calculate it, we average the model-predicted values across the levels of ID, PERIOD and SEQ.

The weights assigned to each combination of levels depend if the factors are nested or not. A factor A is said nested in factor B if each level of A appears in only one level of B. By design, in crossover bioequivalence studies, ID is nested in SEQUENCE (i.e each individual belongs to only one sequence). Thus, in PKanalix, when ID and SEQUENCE are defined as factors in the linear model, **we assume that ID is nested in SEQUENCE when calculating the least square means**. No further nesting is assumed, nor can be specified. Note that the nesting definition only affects the adjusted means and not the point estimate (ratio) and its confidence interval.

Excluded individuals

In case of a parallel design, individuals for which the NCA parameter cannot be calculated are excluded.

In case of a crossover design, incomplete individuals are excluded. Incomplete individual are those that do not have a value for each period, either because these were no concentration data for this period (period missing in the data set) or because the computed NCA parameter could not be calculated for this period (e.g because the data were insufficient to calculate the terminal slope). In case of a non-repeated crossover design, excluding the individuals with missing data or not has no impact on the calculation of the point estimate and confidence interval. In case of a repeated crossover design (i.e when the individuals receive several times the same formulation), excluding the individuals for which the NCA parameter for one or more periods is missing has an impact on the point estimate calculation and confidence interval. The common practice is to exclude incomplete individuals (see also the *EMA Guideline On The Investigation Of Bioequivalence*, page 14).

Excluded individuals are calculated for each NCA parameter separately, as some individuals may have values for all periods for the Cmax but not the AUCINF_obs for instance.

In the Results > Confidence intervals, the number of individuals "N" does not count the excluded individuals. Thus, in case of a crossover design, the number of individuals contributing to ref and to test is the same. **In the BE plots, only the individuals included in the bioequivalence analysis are shown.**

ANOVA

The sum of squares presented in the ANOVA table of the Results tab are **type-I sequential sum of squares**. In case of an unbalanced design (i.e not the same number of individuals receiving RT versus TR), the type-I sum of squares depends on the order of the included factors. In PKanalix, **the enforced order is SEQUENCE + ID + PERIOD + FORMULATION (+ ADDITIONAL)**, following the SAS sample codes provided by the FDA (see *Statistical Approaches to Establishing Bioequivalence*, 2001, Appendix E) and EMA (see [Questions & Answers: positions on specific questions addressed to](#)

Coefficient of variation

In Results > Coefficients of variation, the SD corresponds to the standard deviation of the residuals, i.e to the standard deviation of the normal random variable ε in the linear model $\log(\text{NCAparam}) = \text{intercept} + \beta_{\text{test}}[\text{if FORM=test}] + \varepsilon$ (for the typical model of a parallel design for instance).

The coefficient of variation is then calculated as:

- without log-transformation: $CV = SD/LSM_{ref}$ (with LSM_{ref} the least square mean for the reference formulation, see above)
- with log-transformation: $CV = 100 \times \sqrt{\exp(SD^2) - 1}$


5.2. Bioequivalence settings

This page describes the settings of the bioequivalence calculations.

Parameters to compute

In the section “Parameters to compute”, the user can move the parameters from the “Parameters” column (list of all available parameters) to the “Computed parameters” column to select them for the NCA calculations. Among the parameters selected for the NCA calculations, the user can choose which parameters should be included in the bioequivalence analysis by ticking the box in the “BE” column. Regulatory guidelines usually recommend to perform the bioequivalence on the log-transformed NCA parameters. This can be set using the toggles in the “log-transform” column.

- **BE:** parameters on which to perform the bioequivalence analysis. *Defaults:* Cmax, AUClast, and AUCINF_obs if plasma single-dose, or AUCtau and Cmax if plasma steady-state, or AURC_last and Max_Rate if urine data.
- **log-transform:** whether the bioequivalence should be computed on the log-transformed NCA parameter (toggle on) or on the NCA parameter directly (toggle off). *Default:* log-transform toggle on, except for Tmax.




Bioequivalence design

The study design, parallel or crossover (repeated or not), is detected automatically based on the structure of the ID and OCCASION columns of the data set, see the [Calculation rules](#) for more details. The user can then select one or several fixed effect factors to be included in the linear model. By default, for a parallel design, only the formulation is included. For a crossover design, the factors ID, PERIOD, FORMULATION and SEQUENCE are included by default, if present in the data set. When both ID and SEQUENCE are given, ID is considered as nested in SEQUENCE. Additional factors can also be added, to be chosen among the (STRATIFICATION) CATEGORICAL and CONTINUOUS COVARIATES columns of the data set. **All factors are considered as fixed effects and random effects are not available. Nesting of ID in SEQUENCE is considered automatically, but further nesting cannot be specified.** Note that the nesting affects the calculation of the adjusted means but not of the bioequivalence ratio.

- **Study design:** parallel or crossover (repeated or not-repeated). Detected automatically based on the structure of the ID and OCCASION columns of the data set (see [Calculation rules](#)). Cannot be changed by the user.
- **Factors in the linear model:**
 - **id [if crossover]:** data set column representing the ID of the individuals. *Default:* the column tagged as ID in the data set. This can be set to “none” if the ID should not be used as a factor in the linear model. When a “sequence” column is defined, “id” is considered as nested in “sequence”.
 - **period [if crossover]:** data set column representing the different periods for each individual. *Default:* the first column tagged as OCCASION in the data set. This setting can be set to “none” if the period should not be used as a factor in the linear model.
 - **formulation:** data set column representing the different groups for which the bioequivalence is calculated, typically the different formulations of the drug products. *Default:* the column with header ‘trt’, ‘treatment’, ‘form’, or ‘formulation’ (case insensitive) if it exists, otherwise first CATEGORICAL COVARIATE column. It is mandatory to include the formulation in the linear model.
 - **reference:** among the categories of the categorical covariate selected as formulation, the one considered as the reference. The bioequivalence will be computed for all categories compared to the reference category. *Default:* ‘r’ or ‘ref’ or ‘reference’ (case insensitive), otherwise the category of the first occasion of the first indiv of the data set.


- **sequence [if crossover]:** data set column representing the different sequences. Default: the column with header 'seq', or 'sequence' (case insensitive) if it exist, otherwise "none". This setting can be set to "none" if the sequence should not be used as a factor in the linear model.
- **additional:** data set columns tagged as CATEGORICAL or CONTINUOUS covariates, which have not already been selected above, can be added as additional factors in the linear model.



Confidence interval

Several settings are more specific to the calculation of the confidence interval. They can be changed in the pop-up settings window of the Bioequivalence task.

- **Level:** level of the confidence interval. *Default:* 90% confidence interval. This means that there is a 90% probability that the true ratio is included in the calculated confidence interval.
- **Bioequivalence limits:** limits within which the confidence interval should be in order to conclude for bioequivalence. *Default:* [80,125].
- **Degrees of freedom [if parallel design]:** method to calculate the degrees of freedom involved in the confidence interval calculations, either using the residuals (this assumes equal variance between the groups or using the Welch-Satterthwaite formula (this assumes possibly unequal variances). *Default:* Welch-Satterthwaite



6. Plots

For the CA task, PKanalix generates automatically specific plots to help in model diagnosis. Plots are interactive for better data exploration. You can modify display options (lin/log scales, axis limits, etc), stratify (split, color, filter) them using covariates from a dataset, and change the graphical preferences (font size, colors etc).

In the PLOTS tab you will find the following plots:


- **Individual fits** This plot displays the individual predictions using the individual estimated parameters (or one set of common parameters if you enable the "pooled fit" option) w.r.t. time on a continuous grid, with the observed data overlaid.
- **Individual parameters vs covariates** This plot displays estimated individual parameters w.r.t. the covariates in the dataset.
- **Distribution of the individual parameters.** This plot displays the estimated population distributions of the individual parameters.
- **Correlation between individual parameters** This plot displays scatter plots for each pair of individual parameters.
- **Observations vs predictions** This plot displays observations w.r.t. the predictions computed using the individual parameters (or one set of common parameters if you enable the "pooled fit" option).

Stratification & Preferences


Stratification of plots in PKanalix works in the same way as in Monolix. Complete description of available options and examples is in the Monolix documentation [here](#).

Save/Export plots and charts data


You can export each single plot as an image (png or svg) using an icon on top of each plot. They are saved in the results folder > PKanalix > ChartsFigures.



To export all plots at once use the top menu option: "Export > Export Plots". You can also export plots data as tables (in .txt files), for example to re-generate them with external tools, with "Export > Export charts data".



In addition, you can export plots automatically after each CA run. To select this option enable the toggle buttons in the PKanalix Preferences: top menu > Settings > Preferences:



Next section: [CA Individual fits plot](#)


6.1. Observed data

It is always good to have a look first at the observed data before running the parameter estimation. Indeed, it is very convenient to see if all the

data is consistent, or if some outliers appear.


The purpose of this plot, also called a spaghetti plot, is to display the original data w.r.t. time.

In the example below, the concentration of warfarin from the [warfarin data set](#) is displayed. A subject is highlighted in yellow by hovering on the line.




One can plot the output in a log-scale to have a better evaluation of the


elimination part for example as in the figure below.




The data can be displayed as dots and lines (by default), or only as dots or only as lines, as in the figure below. In addition, a trend line can be overlaid on the plot, based on mean values of the observed data pooled in bins. The user can choose arithmetic or geometric mean, and can add error bars representing either the standard deviations or standard errors. The mean and error values can be displayed as fixed labels next to the bars, or as tooltips when hovering on a bar.



Bin limits can be displayed as vertical lines, and the user can change the number of bins as well as binning criteria.




An interesting feature is the possibility to display the dosing time as on the figure below. The individual dosing time of the individual can be displayed when the user hovers an individual or always visible. Starting in version 2021R1 on, dosing times corresponding to doses with null amounts are not displayed.




The units as well as additional descriptive statistics of the observations can be displayed. We propose


- The total number of subjects
- The average number of doses per subject
- The total, average, minimum and maximum number of observations per individual.



In addition, if we split the graphic with a covariate, these numbers are recomputed to display the information of the group as in the following plot.




In PKanalix 2023R1 it is now possible to display the mean curves, split by covariate, in a single plot. By switching on the "Merged splits" option, the mean curves of the two wieght groups, which appeared in previous versions in separate plots, are merged into a single plot.




6.1.1. Covariate viewer

It is possible to show the matrix of all the covariates as on the following figure:




It is possible to display one covariate vs another one, selecting in the display panel which covariate to look at.


Here we display the weight versus the height and show the correlation coefficient as an information:



Categorical covariates w.r.t. other categorical covariates are displayed as a histogram (stacked or grouped),




and continuous covariates w.r.t. categorical covariates are displayed as a boxplot as below ("+" are outlier points below $Q1 - 1.5IQR$ (interquartile range) or above $Q3 + 1.5IQR$):



Covariate statistics

Starting from the 2021 version, it is also possible to get the covariate statistics in a dedicated frame of the interface:




All the covariates (if any) are displayed and a summary of the statistics is proposed. For continuous covariates, minimum, median and maximum values are proposed along with the first and third quartile, and the standard deviation. For categorical covariates, all the modalities are displayed along with the number of each. Note the "Copy table" button that allows to copy the table in Word and Excel. The format and the display of the table will be preserved.

6.1.2. Bivariate data viewer

In case of several continuous outputs, one can plot one type of observation w.r.t. another one as on the following figure for the warfarin data set. The direction of time is indicated by the red arrow. Linear interpolation is used for observations of different types that are not at the same time.

Labels for x and y axes correspond to the observation names in the model for data types that are mapped to an output of the model, or **observation ids** used in the dataset for data types not captured by the model. It is possible to choose which observation id to display on the X or Y axis in the settings panel on the right.



6.2. NCA individual fits

Purpose


This plot can be used for two purposes:

- show the terminal slope (lambdaZ) linear regression for each occasion of each individual (default)
- show the individual concentration data with one individual per subplot and all occasions

LambdaZ terminal slope

By default, the settings let you see each occasion of each individual on a subplot, in y log-scale with the lambdaZ linear regression. **Note that the NCA parameters Lambda_z and Lambda_z intercept need to be selected in the “Computed parameters” of NCA settings in order to display the lambdaZ regression line.**


Selecting the “Information” toggle allows you to display information about the rule used to calculate the lambdaZ (“Adjusted R2”, “R2”, “Interval”, “Points” or “custom” – when points have been manually selected), the number of points included in the calculation and the R2 or adjusted R2 criteria. Again, **the NCA parameters Rsq_adjusted, Rsq, and No_points_lambda_z must be selected in “Computed parameters” in the NCA settings in order to appear in the information box.**



Individual concentration data

The plot can also be used to visualize the individual concentration data with one individual per subplot and occasions all together. For a crossover study with two periods and two formulations for instance, one can:

- display the “Lines”
- display the “unused data” as “colored”
- hide the “LambdaZ regression”
- untick the “split occasions”
- set the y axis in linear scale
- color by the Formulation in the stratify tab



Settings

- **Grid arrange (top left above the plot):** The user can define the number of subplots displayed, as well as the number of rows and the number of columns. This choice is not saved. See below the “Default number of individuals per page” for a setting that is saved.
- **General**
 - Legend: hide/show the legend. The legends adapts automatically to the elements displayed on the plot. The same legend box applies to all subplots and it is possible to drag and drop the legend at the desired place.
 - Grid : hide/show the grid in the background of the plots.
 - Information: hide/show
 - the continuous covariates: continuous covariates values are displayed for each individual
 - the goodness of fit: the rule used to calculate the lambdaZ (“Adjusted R2”, “R2”, “Interval”, “Points” or “custom” – when points have been manually selected), the number of points used to calculate the lambdaZ and the R2 value (if rule=R2) or the adjusted R2 (if rule different from R2). **Note that the NCA parameters Rsq_adjusted, Rsq, and No_points_lambda_z must be selected as “Computed parameters” in the NCA settings in order to appear in the information box.**
 - Units: hide/show the units next to the NCA parameter names. The color and surrounding characters for the units can be chosen in the plot “Preferences”.

- Dosing times: hide/show dosing times as vertical lines for each subject.
- Link zoom between plots: activate the linked zoom for all subplots. The same zooming region can be applied on all individuals only on the x-axis, only on the Y-axis or on both (option “none” in the “zoom constraint”).
- *Display*
 - Dots: hide/show the observed data as dots
 - Lines: hide/show the connecting lines between the observed data points
 - Censored data [if censored data present]: hide/show the data marked as **censored (BLQ)**, shown as a red dot according to the “BLQ method before/after Tmax” choice in the **NCA settings**.
 - Unused data: hide/show the data points not used for the lambdaZ calculation
 - greyed: display the data points not used in the lambdaZ calculation as greyed
 - colored: display the data points not used in the lambdaZ calculation with the default color (set in the plot “Preferences”)
 - Split occasions [if occasions present]: Split the individual subplots by occasions or group occasions on a single subplot for each individual.
 - Default number of individuals per page: Set the number of subplots to display. This setting can be saved, whereas the “grid arrange” on the top left corner is not.
 - LambdaZ regression: hide/show the lambdaZ linear regression
- *Sorting*: Sort the subjects by ID, goodness of fit (R2 or adjusted R2 value) or continuous covariate value in ascending or descending order.


6.2.1. NCA parameters with respect to covariates

Purpose

The figure displays the individual parameters as a function of the covariates. It allows to identify correlation effects between the individual parameters and the covariates.


Identifying correlation effects

In the example below, we can see the parameters Cl_{obs} and C_{max} with respect to the covariates: the weight WT, the age AGE and the sex category.



Visual guidelines

In order to help identifying correlations, regression lines, spline interpolations and correlation coefficients can be overlaid on the plots for continuous covariates. Here we can see a strong correlation between the parameter Cl_{obs} and the age and the weight.




Overlay of values on top of the boxplot

The NCA parameter values used to generate the boxplots for the categorical covariates can be overlaid on top of the boxplots, which enables to better visualize the distribution in case of a small number of individuals.




Highlight

Hovering on a point reveals the corresponding individual and, if multiple individual parameters have been simulated from the conditional distribution for each individual, highlights all the points points from the same individual. This is useful to identify possible outliers and subsequently check their behavior in the observed data.



Selection


It is possible to select a subset of covariates or parameters, as shown below. In the selection panel, tick the items you would like to see and click "Accept". This is useful when there are many parameters or covariates.



Stratification


Stratification can be applied by creating groups of covariate values. As can be seen below, these groups can then be split, colored or filtered, allowing to check the effect of the covariate on the correlation between two parameters. The correlation coefficient is updated according to the split or filtering.

In addition to regression lines, correlation coefficients can be added to see the correlation between random effects, as well as spline interpolations.




Selection

It is possible to select a subset of parameters, whose pairs of correlations are then displayed, as shown below. In the selection panel, a set of contiguous rows can be selected with a single extended click, or a set of non-contiguous rows can be selected with several clicks while holding the Ctrl key.




Highlight

Similarly to other plots, hovering on a point provides information on the corresponding subject id, and highlights other points corresponding to the same individual.



Stratification: coloring and filtering

Stratification can be applied by creating groups of covariate values. As can be seen below, these groups can then be split, colored and/or filtered allowing to check the effect of the covariate on the correlation between two parameters. The correlation coefficient is updated according to the stratifying action. In the following case, We split by the covariate SEX and color bay 2 categories of AGE.



Settings

- **General**
 - Legend and grid : add/remove the legend or the grid. There is only one legend for all plots.
 - Information: display/hide the correlation coefficient associated with each scatter plot.
- **Display**
 - **Selection.** The user can select some of the parameters to display only the corresponding scatter plots. A simple click selects one parameter, whereas multiple clicks while holding the Ctrl key selects a set of parameters.
 - **Visual cues.** Add/remove the regression line or the spline interpolation.

6.3. Plots

- What kind of plots can be generated by PKanalix?
- Exporting plots

What kind of plots can be generated by PKanalix?

The list of plots below corresponds to all the plots that PKanalix can generate. They are partly generated directly after importing the data, partly after running an analysis task.

Data

- **Observed data:** This plot displays the original data w.r.t. time as a spaghetti plot, along with some additional information.
- **Covariate viewer:** This plot displays the matrix of all the covariates in one plot.
- **Bivariate data viewer:** This plot displays data of different observation types (typically PK and PD data) in one.

NCA

Model for the observations

- **Individual fits:** This plot displays the individual fits: individual predictions using the individual parameters and the individual covariates w.r.t. time on a continuous grid, with the observed data overlaid.

Diagnosis plots based on individual parameters

- **Distribution of the NCA parameters:** This plot displays the estimated distributions of the individual parameters.
- **Correlation between NCA parameters:** This plot displays scatter plots for each pair of individual parameters.
- **NCA parameters vs covariates:** This plot displays the estimators of the individual parameters w.r.t. the covariates.

CA

Model for the observations

- **Individual fits:** This plot displays the individual fits: individual predictions using the individual parameters and the individual covariates w.r.t. time on a continuous grid, with the observed data overlaid.
- **Observations vs predictions:** This plot displays observations w.r.t. the predictions computed using the population parameters or the individual parameters.

Diagnosis plots based on individual parameters


- **Distribution of the CA parameters:** This plot displays the estimated distributions of the individual parameters.
- **Correlation between CA parameters:** This plot displays scatter plots for each pair of individual parameters.
- **CA parameters vs covariates:** This plot displays the estimators of the individual parameters w.r.t. the covariates.

Bioequivalence

- **Sequence by period:** This plot displays scatter plots for each pair of individual parameters.
- **Subject by formulation:** This plot displays for each parameter, the subject-by-formulation interaction and the inter-subject variability.
- **Confidence intervals:** This plot displays the confidence interval for the ratio for each parameter, as well as the BE limits.
- **Individual concentration data:** This plot displays the individual concentration data with one individual per subplot and occasions all together

Exporting plots


The user can choose to export each plot as an image with an icon on top of it, or all plots at once with the menu Export. It is also possible to export plots data as table, for example to build new plots with external tools.



Note that:

- the export starts after the display of the plots,
- the plots are exported in the result folder,
- only plot selected in Plots tasks are exported,
- legends and information frames are not exported.

Automatic exporting can be chosen in the project Preferences (in Settings), as well as the exporting format (png or svg):



Customizing time grid in plots


6.3.1. CA individual fits

Purpose

The figure displays the observed data for each subject, as well as prediction using the individual parameters.


Individual parameters

Information on individual parameters can be used in two ways, as shown below. By clicking on Information (marked in green on the figure) in the General panel, individual parameter values can be displayed on each individual plot. Moreover, the plots can be sorted according to the values for a given parameter, in ascending or descending order (Sorting panel marked in blue). By default, the individual plots are sorted by subject id, with the same order as in the data set.




Special zoom

User-defined constraints for the zoom are available. They allow to zoom in according to one axis only instead of both axes. Moreover, a link between plots can be set in order to perform a linked zoom on all individual plots at once. This is shown on the figure below with observations from the *remifentanyl* example, and individual fits from a two-compartment model. It is thus possible to focus on the same time range or observation values for all individuals. In this example it is used to zoom on time on the elimination phase for all individuals, while keeping the Y axis in log scale unchanged for each plot.



• Censored data

When a data is censored, this data is different to a “classical” observation and has thus a different representation. We represent it as a bar from the censored value specified in the data set and the associated limit.



Settings

- **Grid arrange.** The user can define the number of subjects that are displayed, as well as the number of rows and the number of columns. Moreover, a slider is present to be able to change the subjects under consideration.
- **General**
 - Legend: hide/show the legend. The legends adapts automatically to the elements displayed on the plot. The same legend box applies to all subplots and it is possible to drag and drop the legend at the desired place.
 - Grid : hide/show the grid in the background of the plots.
 - Information: hide/show the individual parameter values for each subject (conditional mode or conditional mean depending on the "Individual estimates" choice is the setting section "Display").
 - Dosing times: hide/show dosing times as vertical lines for each subject.
 - Link between plots: activate the linked zoom for all subplots. The same zooming region can be applied on all individuals only on the x-axis, only on the Y-axis or on both (option "none").
- **Display**
 - Observed data: hide/show the observed data.
 - Censored intervals [if censored data present]: hide/show the data marked as **censored (BLQ)**, shown as a rectangle representing the censoring interval (for instance [0, LOQ]).
 - Split occasions [if IOV present]: Split the individual subplots by occasions in case of **IOV**.
 - Number of points of the calculation for the prediction
- **Sorting:** Sort the subjects by ID or individual parameter values in ascending or descending order.

By default, only the observed data and the individual fits are displayed.


6.3.2. CA parameters vs covariates

Purpose

The figure displays the individual parameters as a function of the covariates. It allows to identify correlation effects between the individual parameters and the covariates.


Identifying correlation effects

In the example below, we can see the parameters Cl and $V1$ with respect to the covariates: the weight WT , the age AGE and the sex category.




Visual guidelines

In order to help identifying correlations, regression lines, spline interpolations and correlation coefficients can be overlaid on the plots for continuous covariates. Here we can see a strong correlation between the parameter CI and the age and the weight. This is the same with V1.



Highlight

Hovering on a point reveals the corresponding individual and, if multiple individual parameters have been simulated from the conditional distribution for each individual, highlights all the points from the same individual. This is useful to identify possible outliers and subsequently check their behavior in the observed data.



Selection

It is possible to select a subset of covariates or parameters, as shown below. In the selection panel, a set of contiguous rows can be selected with a single extended click, or a set of non-contiguous rows can be selected with several clicks while holding the Ctrl key. This is useful when there are many parameters or covariates.



Stratification

Stratification can be applied by creating groups of covariate values. As can be seen below, these groups can then be split, colored or filtered, allowing to check the effect of the covariate on the correlation between two parameters. The correlation coefficient is updated according to the split or filtering.



Settings

• General

- Legend and grid : add/remove the legend or the grid. There is only one legend for all plots.
- Information: display/hide the correlation coefficient associated with each scatter plot.


• Display

- **Selection.** The user can select some of the parameters or covariates to display only the corresponding plots. A simple click selects one parameter (or covariate), whereas multiple clicks while holding the Ctrl key selects a set of parameters.
- **Visual cues.** Add/remove a regression line or a spline interpolation.

6.3.3. Observations vs predictions plot


The observations vs predictions plot displays observations (y_{ij}) versus the corresponding predictions (\hat{y}_{ij}) estimated by the CA task. It is a useful tool to detect misspecifications in the structural model. The points shall scatter as close as possible along the $y = x$ line (black solid line). If a large amount of points are not scattered close to the identity line, then this might be due to the chosen structural model not capturing certain kinetics. Moreover, the distribution of the dots should be symmetrical around the $y=x$ line, meaning that the model prediction is sometimes above and sometimes below the observed data.

This plot is generated automatically when you run the compartmental analysis task. You can access it from the CA section in the plots list.



Visual Guides

In addition to the observed data, in the DISPLAY settings panel you can add: the identity line $y = x$ and spline interpolation.




Highlight

Hovering on any point of observed data shows:

- the subject id and time corresponding to this point
- highlights all the points corresponding to this subject
- segments linking all points corresponding to the same individual to visualize the time chronology.


Highlighting works across plots. An individual data highlighted in this plot, will be highlighted in every other plot to help you in the model diagnosis.



Log scale

A log scale allows to focus on low observation values. You can set it in the AXIS section of the settings panel for each axis separately or for both together.

The example below displays the predicted concentrations of **remifentanyl**, modeled via a two-compartment model with a linear elimination. In this case, the log-log scale reveals a misspecification of the model: the small observations are under-predicted. These observations correspond to late time points, meaning that the last phase of elimination is not properly captured by the two-compartment model. A three-compartment model might give better results.



Settings

- **General**
 - Legend and grid : add/remove the legend or the grid.
- **Display**
 - BLQ data : show and put in a different color the data that are BLQ (Below the Limit of Quantification)
 - Visual cues: add/remove visual guidelines such as the line $y = x$ or a spline interpolation.
- **Axis**
 - Use log scale
 - Set automatic or manual (user – defined step) tick values
 - Edit axis labels

Stratification

You can split, color or filter the plot as described for the **observed data**.


6.3.4. Correlation between CA parameters

Purpose


This plot displays scatter plots for each pair of parameters. It allows to identify correlations between parameters, which can be used to see the results of your analysis and see the coherence of the parameters for each individuals.

Example

In the following example, one can see pairs of parameters estimated for all parameters.



In addition to regression lines, correlation coefficients can be added to see the correlation between random effects, as well as spline interpolations.



Selection

It is possible to select a subset of parameters, whose pairs of correlations are then displayed, as shown below. In the selection panel, a set of contiguous rows can be selected with a single extended click, or a set of non-contiguous rows can be selected with several clicks while holding the Ctrl key.

▼
Display

All

CI


V1

Q

V2


Highlight

Similarly to other plots, hovering on a point provides information on the corresponding subject id, and highlights other points corresponding to the same individual.



Stratification: coloring and filtering

Stratification can be applied by creating groups of covariate values. As can be seen below, these groups can then be split, colored and/or filtered allowing to check the effect of the covariate on the correlation between two parameters. The correlation coefficient is updated according to the stratifying action. In the following case, We split by the covariate SEX and color by 2 categories of AGE.



Settings

- *General*
 - Legend and grid : add/remove the legend or the grid. There is only one legend for all plots.
 - Information: display/hide the correlation coefficient associated with each scatter plot.
- *Display*
 - *Selection.* The user can select some of the parameters to display only the corresponding scatter plots. A simple click selects one parameter, whereas multiple clicks while holding the Ctrl key selects a set of parameters.
 - *Visual cues.* Add/remove the regression line or the spline interpolation.

6.4. Sequence-by-period plot

Purpose


The sequence-by-period plot allows to visualize, for each parameter, the mean and standard deviation for each period, sequence and formulation.

Description

For each formulation, period and sequence (as defined in the [Bioequivalence settings](#)), the mean and standard deviation is displayed. The dots are colored according to the formulation, and dots belonging to the same sequence connected by lines. Dots are placed on the x-axis according to the occasion value.

If the parameter has a log-transform (see [Bioequivalence settings](#)), the dots represents the geometric mean and the error bars \pm the geometric standard deviation. If the parameter is not log-transformed, the dots represent the arithmetic mean and the error bars \pm the standard deviation. NCA parameter values being NaN are excluded and a warning message is displayed.


The example below shows a repeated crossover with sequences TRTR and TRTR. Three NCA parameters are displayed.



Troubleshooting


• The connecting lines to do appear:

If the plot setting "Lines" is selected and the connecting lines do not appear, it is because no categorical covariate column as been defined as "Sequence" in the [Bioequivalence settings](#).



• There are vertical connecting lines:

This means that the column defined as "Sequence" in the [Bioequivalence settings](#) is not consistent with the sequence of formulation according to the Occasion column. In the tab "Data", you can define an "additional covariate" and select the "covariate sequence (XXX)" with XXX being the column representing the formulation. Then use this column as "Sequence" when defining the factors in the model.



• There are several periods on the x-axis although no column was indicated as "period" in the Bioequivalence settings:

When no column is indicated as "period", we use the column tagged as "occasion" in the data set. If no column has been tagged as "occasion" (parallel design), the occasion is considered to be 1 for all individuals.

Settings

- **General**


- Legend: hide/show the legend. The legends adapts automatically to the elements displayed on the plot. The same legend box applies to all subplots and it is possible to drag and drop the legend at the desired place.
- Grid : hide/show the grid in the background of the plots.
- Units: hide/show the units of the NCA parameters. Note that the color and surrounding character can be chosen in the plot "Preferences".

- **Display**

- Dots: hide/show the dots representing the means
- Lines: hide/show the lines connecting the dots belonging to the same sequence
- Error bars: hide/show the error bars representing the standard deviation
- Individual parameters: select which individual parameters to show on the plot

Stratify

This plot can be split and filtered by the covariates but not colored. The categorical covariate indicated of Formulation in the **Bioequivalence settings** is automatically used to color the dots. **The colors can be changed in the plot "Preferences > Formulations"**:




6.4.1. Subject-by-formulation plot

Purpose

This plot allows to visualize, for each parameter, the subject-by-formulation interaction and the inter-subject variability.

Description

Each subplot corresponds to one NCA parameter and one non-reference (test) formulation. The individual parameter values for the test and reference formulations are displayed side by side. Values corresponding to the same individual are connected by a line. When an individual has received several times the same formulation (repeated crossover), the average for each formulation is calculated and plotted.



Settings

- **General**
 - Legend: hide/show the legend. The legends adapts automatically to the elements displayed on the plot. The same legend box applies to all subplots and it is possible to drag and drop the legend at the desired place.
 - Grid : hide/show the grid in the background of the plots.
 - Units: hide/show the units of the NCA parameters. Note that the color and surrounding character can be chosen in the plot "Preferences".
- **Display**
 - Dots: hide/show the dots representing the individual parameter values
 - Lines: hide/show the lines connecting the dots belonging to the same individual
 - Individual parameters: select which individual parameters to show on the plot


Stratify

The plot can split and filtered but not (yet) colored.

6.4.2. Confidence intervals plot

Purpose

This plot allows to visualize the confidence interval for the ratio for each parameter, as well as the BE limits. Confidence intervals within the BE limits are shown in green while the others are shown in red.



Settings

- **General**
 - Legend: hide/show the legend. The legends adapts automatically to the elements displayed on the plot. The same legend box applies to all subplots and it is possible to drag and drop the legend at the desired place.
 - Grid : hide/show the grid in the background of the plots.
 - Units: hide/show the units of the NCA parameters. Note that the color and surrounding character can be chosen in the plot "Preferences".
- **Display**
 - Median: hide/show the line at ratio=100%.
 - Bioequivalence limits: hide/show the horizontal lines at the BE limits, which can be set in the [BE settings](#)
 - Individual parameters: select which individual parameters to show on the plot

Stratify

This plot cannot be stratified by covariates.

7. Reporting

Starting from the 2023R1 version of the MonolixSuite, PKanalix includes a feature of automated report generation. This feature is available in PKanalix through Export > Generate report. After defining settings in the window that pops up and clicking on the Generate report button, a report in the form of a Microsoft Word (docx) document will be generated in the specified directory. Note that reports generated using command line and lixoftConnectors will contain tables but no plots. Reporting module is not available on CentOS 7.

Here we explain:

- [Different strategies for report generation](#)
- [How to generate a default report](#)
- [How to generate a report using a user-provided template file](#)
 - [How to use placeholders](#)
 - [How to use conditional display of report components](#)
- [How to rename words and phrases used by PKanalix in generated reports](#)
- [How to generate a report using command line arguments](#)

To see examples of a full workflow of report generation in PKanalix, you can visit [this page](#).

Report generation strategies

When selecting Export > Generate report..., there are two different strategies available for report generation:

1. **Using a default template file** – choosing this option will generate the report containing tables of NCA settings, individual NCA metrics, summary of NCA metrics, BE settings and results, CA settings, individual CA parameters, summary of CA parameters and CA cost, as well as all plots available in the interface using their current settings. This option allow to obtain a basic report with all results in one click.
2. **Using a custom template file** – if users choose this option, they will be required to provide a path to the template word file which contains placeholders that will be replaced by tables, plots and keywords when generating the report.

Several other options are available in the pop-up window after the Generate report option is clicked on:

1. **Document watermark** – users can specify a custom text that will appear as a watermark in the report. Font family, font size, color, layout and transparency can be specified, if text is entered into the Text input widget.
2. **Generated report file** – users can choose if the generated report will be saved next to the project file (in that case only the name of the generated report file needs to be given to PKanalix) or in a custom directory (absolute path that includes the generated report file name needs to be given).
3. **Table styles** (only when the custom template file option is selected) – users can choose a style for tables from the list of styles available by default in Microsoft Word. This setting can be overridden by specifying the style directly in the table placeholder. This way custom user-created or company-specific styles can be used.

Default template file

Selecting the option to use the default template file to generate a report gives users a possibility to quickly generate a document containing results of an analysis.

Such a report will consist of:


1. name of the PKnaalix project
2. data set file name,
3. a date of report generation,
4. a table containing settings used in NCA analysis (if NCA was run),
5. a table of individual NCA parameters with mean and CV% (if NCA was run),
6. a table with a summary of NCA parameters (if NCA was run),
7. a table with BE settings (if BE was run),
8. a table with results of BE analysis (if BE was run),
9. a table of compartmental analysis (CA) calculation settings (if CA was run),
10. a table of individual CA parameters and costs with mean and CV% (if CA was run),
11. a table with a summary of CA parameters (if CA was run),
12. a table with total cost, -2LL, AIC and BIC (if CA was run),
13. all plots displayed in the GUI with settings currently defined in the interface.

By clicking on View default template, a Word document containing placeholders for all the elements that will be present in the generated report can be downloaded. This document can then be used as a basis for the creation of a custom template file.

Custom template file

The option to use the custom template file can be selected if a user has previously created a Word document (a template) that typically contains placeholders to be replaced by project-specific information, such as tables of results, plots or settings used in a project. Placeholders are strings with a specific syntax that define the settings and layout with which tables and plots will be generated during the report generation process.

For example, in the picture below, a template file on the right contains a placeholder for the observed data plot. So, after report generation, the report file on the right contains a plot for the observed data instead of this placeholder.



Using placeholders


There are three types of placeholders that can be used:

1. placeholders for project settings,
2. placeholders for plots,
3. placeholders for tables.



Placeholders for project settings are of format `%keyword%` and, during the report generation step, they will be replaced by project or global PKanalix settings. For example, if a custom template file contains the placeholder `%NCA_integralMethod%`, this placeholder will be replaced by the integral method (e.g. "linear up log down") used by a PKanalix project for which the report is being generated. The list of all project settings placeholders can be found [here](#).

Placeholders for plots are of format `<lixoftPLH>settings</lixoftPLH>` and, during the report generation step, they will be replaced by plots that are available in the PKanalix user interface. Placeholders for plots contain settings that describe how the generated plots should look like. All settings that are available for plots in the interface are available in reporting as well, with addition of several others that describe plot dimensions and display of captions. The list of the main plot settings placeholders can be found [here](#).

Placeholders for tables are also of format `<lixoftPLH>settings</lixoftPLH>`. Placeholders for tables contain settings that describe how the generated table should look like and they allow more flexibility in the layout compared to the display of tables in the GUI. Indeed, table placeholders can contain settings that result in generated tables being different from the ones available in the Results tab of the PKanalix interface, such as an option to choose which rows or columns to include in the table, or to split tables in a different direction than the one present in interface. [Here](#) is the list of all settings that can be used for table placeholders.

Placeholders do not need to be written by the users. They can be generated through the interface, using the **Copy reporting placeholder**  button which can be found next to every plot and results table. Clicking on this button will display the corresponding placeholder in a pop-up window. The displayed placeholder can then be copied to a template file and will be replaced by the table or plot, exactly as they are shown in the interface, with all stratification options and display settings applied.

Let's take a look at how "Copy reporting placeholder" buttons work on a specific example. If we load the demo project `censoring.pkx`, go to the Observed data plot and click on the Copy reporting placeholder button above the plot, a pop-up window will appear containing the placeholder for the plot, which can then be used in a report template to generate the plot that looks exactly like the plot currently shown in the interface. As you can see on the screenshot below, the placeholder contains information about plot type, plot size, caption and axes labels.



If we now close the pop-up, enable legend and open the pop-up again, the placeholder will contain the enabled legend setting and, if we put the placeholder in a report template, the generated plot will have a legend present in its corner, as shown in the screenshot on the right.





Figure 1: Concentration-time profiles

If we would like our report to have the Observed data plot stratified by a categorical covariate, STUDY, we can stratify the plot in the interface and click on the Copy reporting placeholder button again to generate a placeholder for this plot, as shown on the images below. You can notice that the generated placeholder now contains settings stratification and layout, which control stratification and plot layout.





Figure 1: Concentration-time profiles

There are several options present inside the "Placeholder for Report Template" pop-up:

1. **indent** – determines number of spaces used for indentation of settings in the placeholder,
2. **flow level** – determines the compactness of placeholders, i.e. how often we start a new line (lower the number, more compact the placeholder will be),
3. **dimensions units** – gives the ability to choose if width and height will be in cm or inches,
4. **lock aspect ratio** – if on, changing width or height will change the other dimension to keep the aspect ratio present on the screen.

It is possible to edit placeholders directly inside the pop-up. Clicking on the **Preview** button will generate and open a Word document containing a plot or a table described with the placeholder. This option is especially useful when a user manually modifies placeholders, to check if the syntax and the behavior are correct.

Copy to clipboard button will copy the placeholder with the <lixoftPLH> tag, ready to be pasted to the custom Word template.

Using areTaskResultsAvailable

Besides placeholders, templates can include the `areTaskResultsAvailable` tag, which can be used to display part of the template conditionally in the report.

For example, if a template contains the text:

```
1. <? areTaskResultsAvailable(nca) ?>
2. Show this if NCA has been run
3. </?>
```

The text “Show this if NCA has been run” will be present in the generated report, only if NCA results are available, otherwise the tags and the text will not be present in the generated report. Besides text, it is possible to include placeholders inside the `areTaskResultsAvailable` tag. Keywords that can be used inside the parentheses are “nca” for non-compartmental analysis, “be” for bioequivalence and “ca” for compartmental analysis.

Renamings

In the PKanalix Preferences, there is an option to replace default names of NCA parameters by custom names (see screenshot below). Changing these options will change the names in the interface (after restarting the program) and in generated reports.

A section Reporting placeholder entry renamings in Preferences can be used to change words and phrases that will appear in generated reports by default. For example, if we give alias “LeastSquareMean” to a placeholder entry “AdjustedMean”, the word “AdjustedMean” in the BE table generated by the reporting feature will be replaced by the word “LeastSquareMean”.

NCA parameters renamings

Re-run or reload the project for changes to take effect

Parameter	Alias
AUCINF_D_obs	AUC _{0-inf} /Dose
AUCINF_D_pred	AUC _{0-inf,pred} /Dose
AUCINF_obs	AUC _{0-inf}
AUCINF_pred	AUC _{0-inf,pred}
AUC_Int	AUC _{0-t1}
AUC_Int_D	AUC _{0-t1} /Dose
AUC_PerCentBack_Ext_obs	% AUC backextrapolated
AUC_PerCentBack_Ext_pred	% AUC backextrapolated _{pred}
AUC_PerCentExtrap_obs	% extrapolated AUC

Reporting placeholder entry renamings

Use the __void__ alias to get an empty table cell

Placeholder entry	Alias
-2LL	
AIC	
AdjustedMean	LeastSquareMean
BE	
BIC	
BLQ	
Bioequivalence	
CLower	
CLRawLower	

In addition, it is possible to rename words found in generated tables, by specifying an argument renamings inside the placeholder. This allows users to change words present in the data set as well. For example, if an imported data set contains a column “Formulation” with values T (for the test drug) and R (for the reference drug), specifying the following inside a placeholder will replace all occurrences of word “T” inside the table with “test” and “R” with “reference”:

```
1. renamings:
2.   T: test
3.   R: reference
```

For categorical covariates, it is possible to specify both the covariate name and the covariate modality to be replaced, in order to avoid ambiguities if category “0” appears in several covariates. In addition, quotes can be used if spaces are needed within the springs:

```
1. renamings:
2.   SEX#0: Male
3.   SEX#1: Female
4.   RACE#0: Caucasian
5.   RACE#1: "African american"
6.   RACE#2: Hispanic
```

Respecting the indentation before the modalities to replace and the space after the “:” is crucial to have the placeholder properly recognized.

Generating reports from the command line

It is possible to generate a report directly through the command line, without using the PKanalix interface. This can be useful when running PKanalix on a server.



Note that reports generated via command line will not contain plots. This is a known limitation of the current implementation of the reporting feature. If generating reports with plots from R or from the command line is important for you, please [let us know](#) so that we take it into account for future versions.

To generate a report via command line, the executable `reportGenerator` can be called from the lib folder (typically `$HOME/Lixoft/MonolixSuite2023R1/lib/`):

```
1. reportGenerator -p pkanalix_project_path
```

Multiple options can be provided to `reportGenerator`:

`-?, -h, -help`

Displays the help.

`-p, -project <path>`

Project file to load

`-tpl, -template <path>`

Template .docx file used as reporting base (if not provided a default report file is generated)

`-onxt, -output-next-project <[true], false>`

Generate report next to project

-op, -output-path <path>	Generated report file name (if output-next-project is true) or complete path (if output-next-project is false)
-tab, -tables-style <style>	Table style sheet applied on generated tables
-wt, -watermark-text <text>	Watermark text
-wff, -watermark-font-family <[Arial], ...>	Watermark font family
-wfs, -watermark-font-size <integer>	Watermark font size
-wc, -watermark-color <[255:0:0]>	Watermark color (RGB)
-wl, -watermark-layout <[horizontal], diagonal>	Watermark layout
-wst, -watermark-semi-transparent <[true], false>	Watermark semi-transparent

7.1. Reporting examples #

NCA summary table split by a covariate

- example: demo **project_Theo_extravasc_SS.pkx**

```
1. <lixoftPLH>
2. data:
3.   task: nca
4.   metrics: [Nobs, mean, CV, min, median, max]
5.   parameters: [Cmax, AUClast, CLss_F, Vz_F, HL_Lambda_z]
6.   display:
7.     units: true
8.     inlineUnits: true
9.     metricsDirection: horizontal
10.    significantDigits: 4
11.    fitToContent: true
12.    stratification:
13.      state: {split: [FORM]}
14.      splitDirection: [h]
15. </lixoftPLH>
```

FORM	ref						test					
Metrics	N	mean	CV(%)	min	median	max	N	mean	CV(%)	min	median	max
C _{max} (mg·L ⁻¹)	12	17.82	24.05	11.74	17.27	27.23	12	13.39	21.24	10.08	12.95	18.88
AUC _{0-last} (h·mg·L ⁻¹)	12	258.35	35.86	138.39	238.73	469.12	12	234.49	25.33	150.55	230	335.41
CL _{18F} (L·h ⁻¹)	12	3.47	35.63	1.71	3.35	5.78	12	3.62	25.32	2.39	3.49	5.31
V _d /F (L)	12	34.16	30	19.47	32.63	62.1	12	41.77	20.99	28.25	41.75	61.96
T _{1/2} (h)	12	7.26	29.4	4.27	7.81	11.59	12	8.51	39.16	5.49	7.36	18.01

BE table for absolute bioavailability

- example: demo **project_parallel_absBioavailability.pkx**

```
1. <lixoftPLH>
2. data:
3.   task: be
4.   table: confidenceIntervals
5.   metrics: [Ratio, CILower, CIUpper]
6.   parameters: [AUCINF_obs, AUClast]
7.   display:
8.     units: true
9.     inlineUnits: true
10.    significantDigits: 4
11.    fitToContent: true
12.    renamings:
13.      bioequivalence: "Relative bioavailability (%) based on:"
14.      CILower: "Lower 90% CI"
15.      CIUpper: "Upper 90% CI"
16. </lixoftPLH>
```

Parameter	Relative bioavailability (%) based on:		
	Ratio	Lower 90% CI	Upper 90% CI
AUC _{0-inf} (h·ug·L ⁻¹)	87.9105	81.0865	95.3089
AUC _{0-last} (h·ug·L ⁻¹)	87.856	81.2746	94.9703

Table of NCA settings

- example: demo **project_censoring.pkx**

Data type	%data_type%
Dose type	%NCA_administrationType%
Weighting for λz slope calculation	%NCA_lambdaWeighting%
Point selection method for λz slope calculation	%NCA_lambdaRule%
AUC calculation method	%NCA_integralMethod%
BLQ before T _{max} treated as	%NCA_blqMethodBeforeTmax%
BLQ after T _{max} treated as	%NCA_blqMethodAfterTmax%

Data type	plasma
Dose type	extravascular
Weighting for λz slope calculation	uniform
Point selection method for λz slope calculation	adjusted R2
AUC calculation method	linear trapezoidal linear
BLQ before T _{max} treated as	zero
BLQ after T _{max} treated as	LOQ/2

Individual NCA parameters for a group

- example: demo **DoseAndLOQ_byCategory.pkx**

```
1. <lixoftPLH>
2. data:
3.   task: nca
4.   metrics: [ID]
5.   parameters: [AUCINF_D_obs, AUCINF_obs, AUClast, AUClast_D, Cmax, Cmax_D, Tmax]
6. display:
7.   units: true
8.   inlineUnits: true
9.   metricsDirection: vertical
10.  significantDigits: 2
11.  fitToContent: true
12. stratification:
13.   state: {split: [STUDY], filter: [[STUDY, [1]]]}
14.   splitDirection: [v]
15. </lixoftPLH>
```

STUDY	ID	AUC _{0-inf} /Dose	AUC _{0-inf}	AUC _{0-last}	AUC _{0-last} /Dose	C _{max}	C _{max} /Dose	T _{max}
SD_400mg	1	0.19	77.18	76.69	0.19	3.63	0.0091	14
	2	0.26	103.32	103.16	0.26	5.91	0.015	14
	3	0.15	60.29	60.23	0.15	6.1	0.015	8
	4	0.28	111.57	111.07	0.28	5.24	0.013	16
	5	0.23	92.78	92.54	0.23	4.83	0.012	8
	6	0.22	88.78	88.51	0.22	4.35	0.011	12
	7	0.22	89.74	89.56	0.22	7.18	0.018	12
	8	0.55	218.6	209.84	0.52	8.35	0.021	8
	9	0.34	134.01	132.98	0.33	5.77	0.014	16
	10	0.22	89.08	88.98	0.22	6.19	0.015	16
	11	0.19	74.35	74.05	0.19	6.28	0.016	10
	12	0.25	100.25	100.13	0.25	6.15	0.015	10

NCA summary table with summary metrics filtered by group

- example: demo **DoseAndLOQ_byCategory.pkx**

```
1. <lixoftPLH>
2. data:
3.   task: nca
4.   metrics: [Nobs, mean, SE, CV, min, median, max, geoMean, harmMean]
5.   parameters: [AUCINF_D_obs, AUCINF_obs, AUClast, AUClast_D, Cmax, Cmax_D, Tmax]
6. display:
7.   units: true
8.   inlineUnits: true
9.   metricsDirection: vertical
10.  significantDigits: 2
11.  fitToContent: true
12. stratification:
13.   state: {split: [STUDY], filter: [[STUDY, [1]]]}
14.   splitDirection: [v]
15. </lixoftPLH>
```

STUDY		AUC _{0-inf} /Dose	AUC _{0-inf}	AUC _{0-last}	AUC _{0-last} /Dose	C _{max}	C _{max} /Dose	T _{max}
SD_400mg	N	12	12	12	12	12	12	12
	mean	0.26	103.33	102.31	0.26	5.83	0.015	12
	SE	0.03	11.8	11.16	0.028	0.36	0.0009	0.92
	CV(%)	39.57	39.57	37.78	37.78	21.33	21.33	26.59
	min	0.15	60.29	60.23	0.15	3.63	0.0091	8
	median	0.23	91.26	91.05	0.23	6.01	0.015	12
	max	0.55	218.6	209.84	0.52	8.35	0.021	16
	GeoMean	0.24	97.86	97.24	0.24	5.71	0.014	11.6
	HarmMean	0.23	93.71	93.29	0.23	5.58	0.014	11.19

Individual covariate table

- example: demo **project_M2000_bolus_SD.pkx**

```
1. <lixoftPLH>
2. data:
3.   task: nca
4.   metrics: [ID]
5.   excludedParameters: all
6.   covariates: all
7.   covariatesAfterParameters: true
8. display:
9.   units: true
10.  inlineUnits: true
11.  metricsDirection: vertical
12.  significantDigits: 4
13.  fitToContent: true
14. </lixoftPLH>
```

ID	AGE	HT	WT	GENDER
271	69	157	63.7	F
272	66.8	165	66.5	M
273	67.7	163	57.1	F
274	73.7	167	81.3	M
275	68	149	61.4	F
276	73.6	189	88.6	M
277	67.8	169	82.6	M
278	73.2	165	66	F
279	66.5	178	80.4	M
280	65.6	173.5	80	M
281	66.8	170	57.9	M
282	73.5	174	78.8	M

Included in LambdaZ table

- example: demo **project_censoring.pkx**

```
1. <lixoftPLH>
2. data:
3.   task: nca
4.   metrics: [ID]
5.   excludedParameters: all
6.   covariates: all
7.   covariatesAfterParameters: true
8. display:
9.   units: true
10.  inlineUnits: true
11.  metricsDirection: vertical
12.  significantDigits: 4
13.  fitToContent: true
14. </lixoftPLH>
```

ID	AGE	HT	WT	GENDER
271	69	157	63.7	F
272	66.8	165	66.5	M
273	67.7	163	57.1	F
274	73.7	167	81.3	M
275	68	149	61.4	F
276	73.6	189	88.6	M
277	67.8	169	82.6	M
278	73.2	165	66	F
279	66.5	178	80.4	M
280	65.6	173.5	80	M
281	66.8	170	57.9	M
282	73.5	174	78.8	M

7.2. List of reporting keywords

This page provides a list of placeholders for project settings that can be used in report templates for PKanalix projects.

Keyword	Description	Possible values
%NCA_administrationType%	Administration type for NCA parameter calculation.	"intravenous" or "extravascular" if one adm id, otherwise "1: intravenous, 2: extravascular" with 1/2 corresponding to the administration id
%NCA_ObsidUsed%	Observation id used for the NCA calculations.	string of the observation id, e.g "PK" or "1"
%data_type%	Type of data defined on the Data page.	"plasma" or "urine" if one observation id, otherwise "PK: plasma, PD: plasma" with PK/PD corresponding to the observation ids
%NCA_integralMethod%	Method for AUC and AUMC calculation and interpolation.	as in the associated dropdown menu (linear trapezoidal linear, linear log trapezoidal, linear up log down, linear trapezoidal linear/log)
%NCA_blqMethodBeforeTmax%	Method by which the BLQ data before Tmax should be replaced.	as in the associated dropdown menu (missing, zero, LOQ, LOQ/2)
%NCA_blqMethodAfterTmax%	Method by which the BLQ data after Tmax should be replaced.	as in the associated dropdown menu (missing, zero, LOQ, LOQ/2)
%NCA_adjR2AcceptanceCriteria%	Value of the adjusted R2 acceptance criteria for the estimation of lambda_Z.	a number
%NCA_extrapAucAcceptanceCriteria%	Value of the AUC extrapolation acceptance criteria for the estimation of lambda_Z.	a number
%NCA_spanAcceptanceCriteria%	Span acceptance criteria for the estimation of lambda_Z.	a number
%NCA_lambdaRule%	General rule for the lambda_Z estimation. Includes information about the main rule, together with settings defined in the "rules" tab.	Examples: "R2 with maximum 500 points, minimum time 20" "Time interval between 10 and 50" "Fixed number of points (3 points)"
%NCA_lambdaWeighting%	Weighting used for the lambda_Z estimation.	as in the associated dropdown menu (uniform, 1/Y, 1/Y ²)
%CA_cost%	Cost function used for CA.	as in the associated dropdown menu ((Ypred - Yobs) ² , (Ypred - Yobs) ² / Yobs, (Ypred - Yobs) ² / Ypred, (Ypred - Yobs) ² / Yobs ² , (Ypred - Yobs) ² / Ypred ² , (Ypred - Yobs) ² / Ypred * Yobs)
%CA_method%	Fit with individual parameters or with the same parameters for all individuals.	pooled fit, individual fit
%CA_blqMethod%	BLQ method for CA.	as in the associated dropdown menu (missing, zero, LOQ, LOQ/2)
%BE_level%	Level of the confidence interval for bioequivalence	a number
%BE_limits%	Limits in which the confidence interval must be to conclude the bioequivalence is true	Two numbers separated by a comma. Example: 0.8, 1.25
%BE_linearModelFactors%	Factors tested in the linear model for bioequivalence.	Data set column names. Example: ID, PERIOD, FORM, SEQ
%BE_design%	Design of a bioequivalence study.	parallel, crossover
%BE_reference%	Categorical covariate modality selected as a reference	One of the modalities present in the data set. Example: REF or IV
%TotalNbSubjects%	Total number of subjects in the data set (all observation ids together).	a number
%TotalNbSubjectsOcc%	Total number of subjects-occasions in the dataset (all observation ids together).	a number

%AvgNbDosesPerSubject%	Average number of doses per subject (any administration id, any occasion, all observation ids).	a number
%TotalNbObservations%	Total number of observations in the data set.	if one obs id, number, if several obs id, PK: xx, PD: xx
%AvgNbObservationsPerID%	Average number of observations per id (all occasions together).	if one obs id, number, if several obs id, PK: xx, PD: xx
%MinNbObservationsPerID%	Minimum number of observations per id (all occasions together).	if one obs id, number, if several obs id, PK: xx, PD: xx
%MaxNbObservationsPerID%	Maximum number of observations per id (all occasions together).	if one obs id, number, if several obs id, PK: xx, PD: xx
%PercentCensoredObservations%	Percentage of censored observations.	if one obs id, number, if several obs id, PK: xx, PD: xx
%StructModelCode%	Mlxtran code for the structural model.	whole structural model file content, as shown in the CA model tab
%reportGenerationDateTime(yyyy-MM-dd HH:mm:ss)%	Report generation date and time. Part of yyyy-MM-dd HH:mm:ss can be removed or reordered.	yyyy-MM-dd HH:mm:ss
%system%	Operating system on which report was generated.	Linux, macOS, Windows
%version%	Version of PKanalix with which report was generated.	2023R1
%project_fileName%	File name of the project (without path).	Example: NCA_run.pkx
%data_fileName%	File name of the data set (without path).	Example: data.csv

7.3. List of table placeholder settings

The following page describes all the settings that can be used in table placeholders in reporting templates. Placeholders for tables can include four different groups of settings. These are:

- “**data**” (required) – containing settings about table content, such as type of the table, and rows and columns that the table will contain,
- “**display**” (optional) – containing information about table display (e.g., number of significant digits to which the values in the table should be rounded, table direction, style, ...
- “**stratification**” (optional) – containing information about table splitting and/or filtering,
- “**renamings**” (optional) – containing various words present in the data set or PKanalix that should be reworded in the table.

Settings inside these groups of settings need to be indented or put inside curly brackets. It is crucial to respect the indentation rules and the space after the “:” for the placeholder to be properly interpreted.

Data

Data settings define the type of the table that will be generated, as well as rows and columns that the table will contain. Here is the list of all settings of the data group of settings, along with their descriptions:

Setting	Table	Required	Possible values	Description
task	All	yes	be, ca, nca	Task of which the table represents the results.
table	Points included for lambdaZ (NCA) Anova table (BE) Table of CVs (BE) Table of confidence intervals (BE) Cost table (CA)	yes	pointsIncludedForLambdaZ, anova, coefficientsOfVariation, confidenceIntervals, cost	Table type. If empty, CA or NCA parameters table will be generated (depending on the task argument).
formulation	Table of confidence intervals (BE)	no (default: all)	name of one of the test formulations present in the data set, or “all”	In case of confidenceIntervals table and multiple test formulations, setting can be used to select one of the non-reference formulations.
metrics	Table of parameters (CA or NCA)	no	“all” or one or several of: ID, min, Q1, median, Q3, max, mean, SD, SE, CV, Ntot, Nobs, Nmiss, geoMean, geoSD, geoCV, harmMean	If equal to “all”, all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Points included for lambdaZ (NCA)	no	“all” or one or several of: ID, time, concentration, BLQ, includedForLambdaZ	If equal to “all”, all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.

metrics	Table of confidence intervals (BE)	no	"all" or one or several of: AdjustedMean, N, Difference, CIRawLower, CIRawUpper, Ratio, CILower, CIUpper, BE	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Coefficient of variation (BE)	no	"all" or one or several of: SD, CV	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Anova table (BE)	no	"all" or one or several of: DF, SUMSQ, MEANSQ, FVALUE, PR(>F)	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Cost table (CA)	no	"all" or one or several of: Cost, -2LL, AIC, BIC	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
excludedMetrics	All	no (default: none)	Same arguments as for metrics (except "all").	Metrics that will be excluded from the table. Handy when a user wants to specify metrics: "all" and exclude a few with excludedMetrics.
parameters	Table of parameters (CA/NCA)	no (default: all)	One or more of parameters calculated in NCA or CA, including Cost if CA, or "all".	Parameters that will be included in the table. Used when table setting isn't specified.
excludedParameters	Table of parameters (CA/NCA)	no (default: none)	Same arguments as for parameters, excluding "all".	Parameters that will be excluded from the table. Handy when a user wants to keep most of them.
factors	Anova table (BE)	no (default: all)	One or more of the headers of columns used in the BE linear model, including Residuals, or "all".	Factors included in the linear model.
excludedFactors	Anova table (BE)	no (default: none)	Same arguments as for factors, excluding "all".	Factors included in the linear model.
ids	Table of parameters (CA/NCA) Points included for lambdaZ (NCA)	no (default: all)	One or several of IDs present in the dataset, or "all".	IDs of subjects whose parameters should be shown in the table. The summary statistics will be calculated on all individuals, not just ones shown in the table. If not present, default is "all".
excludedIds	Table of parameters (CA/NCA) Points included for lambdaZ (NCA)	no (default: none)	Same arguments as for ids, excluding "all".	IDs of subjects whose parameters should be excluded from the table. The summary statistics will be calculated on all individuals, not just ones shown in the table.
nbOccDisplayed	Table of parameters (CA/NCA) Points included for lambdaZ (NCA)	no (default: -1 meaning all)	-1 or an integer	-1 means all occasion levels are displayed. 0 mean none of the occasion levels are displayed. 1 means only the first occasion level is displayed, etc.
covariates	Table of parameters (CA/NCA)	no (default: all)	One or more of covariates present in the data set, or "all".	Which covariates to include in the table.
covariatesAfterParameters	Table of parameters (CA/NCA)	no (default: true)	true, false	If covariates should be shown after parameters in the table.

Data

Data settings define the type of the table that will be generated, as well as rows and columns that the table will contain. Here is the list of all settings of the data group of settings, along with their descriptions:

Setting	Table	Required	Possible values	Description
task	All	yes	be, ca, nca	Task of which the table represents the results.
table	Points included for lambdaZ (NCA) Anova table (BE) Table of CVs (BE) Table of confidence intervals (BE) Cost table (CA)	yes	pointsIncludedForLambdaZ, anova, coefficientsOfVariation, confidenceIntervals, cost	Table type. If empty, CA or NCA parameters table will be generated (depending on the task argument).
formulation	Table of confidence intervals (BE)	no (default: all)	name of one of the test formulations present in the data set, or "all"	In case of confidenceIntervals table and multiple test formulations, setting can be used to select one of the non-reference formulations.
metrics	Table of parameters (CA or NCA)	no	"all" or one or several of: ID, min, Q1, median, Q3, max, mean, SD, SE, CV, Ntot, Nobs, Nmiss, geoMean, geoSD, geoCV, harmMean	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Points included for lambdaZ (NCA)	no	"all" or one or several of: ID, time, concentration, BLQ, includedForLambdaZ	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Table of confidence intervals (BE)	no	"all" or one or several of: AdjustedMean, N, Difference, CIRawLower, CIRawUpper, Ratio, CILower, CIUpper, BE	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Coefficient of variation (BE)	no	"all" or one or several of: SD, CV	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
metrics	Anova table (BE)	no	"all" or one or several of: DF, SUMSQ, MEANSQ, FVALUE, PR(>F)	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.

metrics	Cost table (CA)	no	"all" or one or several of: Cost, -2LL, AIC, BIC	If equal to "all", all available rows or columns will be present in the table. If specific metrics are given, only those rows/columns will be available in the table.
excludedMetrics	All	no (default: none)	Same arguments as for metrics, excluding "all".	Metrics that will be excluded from the table. Handy when a user wants to keep most of them.
parameters	Table of parameters (CA/NCA)	no (default: all)	One or more of parameters calculated in NCA or CA, including Cost if CA, or "all".	Parameters that will be included in the table. Used when table setting isn't specified.
excludedParameters	Table of parameters (CA/NCA)	no (default: none)	Same arguments as for parameters, excluding "all".	Parameters that will be excluded from the table. Handy when a user wants to keep most of them.
factors	Anova table (BE)	no (default: all)	One or more of the headers of columns used in the BE linear model, including Residuals, or "all".	Factors included in the linear model.
excludedFactors	Anova table (BE)	no (default: none)	Same arguments as for factors, excluding "all".	Factors included in the linear model.
ids	Table of parameters (CA/NCA) Points included for lambdaZ (NCA)	no (default: all)	One or several of IDs present in the dataset, or "all".	IDs of subjects whose parameters should be shown in the table. The summary statistics will be calculated on all individuals, not just ones shown in the table. If not present, default is "all".
excludedIds	Table of parameters (CA/NCA) Points included for lambdaZ (NCA)	no (default: none)	Same arguments as for ids, excluding "all".	IDs of subjects whose parameters should be excluded from the table. The summary statistics will be calculated on all individuals, not just ones shown in the table.
covariates	Table of parameters (CA/NCA)	no (default: all)	One or more of covariates present in the data set, or "all".	Which covariates to include in the table.
covariatesAfterParameters	Table of parameters (CA/NCA)	no (default: true)	true, false	If covariates should be shown after parameters in the table.

Display

Display group of settings defines how the information in the tables will be displayed. Here is the list of all settings for the display group of settings:

Setting	Table	Required	Possible values	Description	Example(s)
units	Table of parameters (CA/NCA) Anova table (BE) Table of CVs (BE) Table of confidence intervals (BE) Cost table (CA)	no (default: true)	true, false	If true, units will be output next to parameters.	units: true
CDISCNames	Table of parameters (CA/NCA) Anova table (BE) Table of CVs (BE) Table of confidence intervals (BE) Cost table (CA)	no (default: false)	true, false	If true, CDISC names of parameters will be shown in the table.	CDISCNames: true
style	All	no (default: selected in Generate report pop-up)	Name of table styles present in the template document.	Microsoft Word document template style to apply to the table. This setting overrides the setting in the Generate report pop-up.	style: "Medium Grid 1 – Accent 1"
metricsDirection	All	no (default: vertical)	vertical, horizontal	Direction of metrics.	metricsDirection: horizontal
significantDigits	All	no (default: taken from Preferences)	positive integers	Number of significant digits values in the table will be rounded to.	significantDigits: 3
trailingZeros	All	no (default: taken from Preferences)	true, false	If trailing zeros should be shown in the tables.	trailingZeros: true
inlineUnits	Table of parameters (CA/NCA) Anova table (BE) Table of CVs (BE) Table of confidence intervals (BE) Cost table (CA)	no (default: true)	true, false	If units should be inline with parameters. If false, there will be a new line between a parameter name and a unit.	inlineUnits: false
fitToContent	All	no (default: true)	true, false	If true, width of the table will be equal to the width of content, otherwise width of the table will be equal to the width of the page.	fitToContent: false
caption	All	no (default: no caption)	a phrase inside quotation marks	Caption that will appear next to the table.	caption: "Individual NCA parameters"
captionAbove	All	no (default: false)	true, false	If true, caption will be positioned above the table. If false, caption will be positioned below the table.	captionAbove: true
fontSize	All	no (default: 12)	a number	Font size of the table content (will not be applied to the caption).	fontSize: 10

Stratification

Stratification group of settings defines how the table will be split. List of stratification settings:

Setting	Table	Required	Possible values	Description	Example(s)
---------	-------	----------	-----------------	-------------	------------

state	Table of parameters (CA/NCA)	no	Two possible subsetings: – split – can contain header names of covariates and/or occasions – filter – can contain header names of covariates and/or occasions and number(s) that describe alphabetical order of modalities	Contains information of how to split the table.	state: {split: [FORM]} state: {split: [FORM, OCC]} state: {filter: [[FORM, [1]]]} state: {filter: [[SEQ, [1, 2]]]} state: {filter: [[FORM, [1]], [OCC, [1]]]} state: {split: [FORM, OCC], filter: [[OCC, [1]], [SEQ, [1]]]}
splitDirection	Table of parameters (CA/NCA)	no	One or more of v, h	In which direction the table should be split. Number of arguments corresponds to number of covariates by which the table was split.	splitDirection: [v] splitDirection: [v, v] splitDirection: [v, h]

Let's take a look at the demo project project_crossover_bioequivalence.pkx. It contains data from a 2×2 crossover bioequivalence study of two controlled release formulations of theophylline. We would like to generate a placeholder for the summary table of NCA parameters and split it across sequences and formulations. If we select that we only want to compute Cmax and AUClast, here is how a default placeholder looks like when we click on the placeholder button next to the NCA summary table and which table it produces:

1.	<lixoftPLH>
2.	data:
3.	task: nca
4.	metrics: [min, Q1, median, Q3, max, mean, SD, SE, CV, Ntot, Nobs, Nmiss, geoMean, geoSD, geoCV, harmMean]
5.	parameters: [AUClast, Cmax]
6.	display:
7.	units: true
8.	inlineUnits: true
9.	metricsDirection: horizontal
10.	significantDigits: 4
11.	fitToContent: true
12.	</lixoftPLH>

	min	Q1	median	Q3	max	mean	SD	SE	CV(%)	Ntot	N	Nmiss	GeoMean	GeoSD	GeoCV	HarmMean
AUC _{0-tlast} (h·mg·L ⁻¹)	45.6709	116.8202	142.4099	179.027	236.6814	148.5547	46.7732	7.7955	31.4855	36	36	0	140.3561	1.4369	37.4696	130.4717
C _{max} (mg·L ⁻¹)	3.56	6.12	8.76	9.94	11.53	8.1806	2.2192	0.3699	27.1282	36	36	0	7.8445	1.3597	31.4694	7.4633

If we want to split the table by the period and the formulation, we can select the appropriate checkboxes on the right and the default placeholder will change automatically. However, doing this always splits the table vertically. We can also remove some of the statistics that we want to exclude from the table by modifying the metrics argument.

1.	<lixoftPLH>
2.	data:
3.	task: nca
4.	metrics: [min, median, max, mean, CV, Nobs, geoMean, geoCV]
5.	parameters: [AUClast, Cmax]
6.	display:
7.	units: true
8.	inlineUnits: true
9.	metricsDirection: horizontal
10.	significantDigits: 4
11.	fitToContent: true
12.	stratification:
13.	state: {split: [FORM, Period]}
14.	splitDirection: [v, v]
15.	</lixoftPLH>

FORM	Period	Metrics	N	mean	CV(%)	GeoMean	GeoCV
ref	1	AUC _{0-tlast} (h·mg·L ⁻¹)	9	154.8982	26.3033	150.3623	26.1424
		C _{max} (mg·L ⁻¹)	9	9.82	11.8862	9.7571	12.1369
	2	AUC _{0-tlast} (h·mg·L ⁻¹)	9	147.6183	36.178	137.7505	43.7537
		C _{max} (mg·L ⁻¹)	9	8.3244	32.019	7.8713	39.0231
test	1	AUC _{0-tlast} (h·mg·L ⁻¹)	9	133.3989	39.8426	122.3564	49.8543
		C _{max} (mg·L ⁻¹)	9	6.4689	28.8987	6.2384	29.1102
	2	AUC _{0-tlast} (h·mg·L ⁻¹)	9	158.3035	26.9691	153.1323	28.1126
		C _{max} (mg·L ⁻¹)	9	8.1089	22.1858	7.9035	25.4489

If we want to change the direction of splitting, for example, make a split by Period horizontal, instead of vertical, we can modify the splitDirection argument.

1.	<lixoftPLH>
2.	data:
3.	task: nca
4.	metrics: [Nobs, mean, CV, geoMean, geoCV]
5.	parameters: [AUClast, Cmax]
6.	display:
7.	units: true
8.	inlineUnits: true
9.	metricsDirection: horizontal
10.	significantDigits: 4
11.	fitToContent: true
12.	stratification:
13.	state: {split: [FORM, Period]}
14.	splitDirection: [v, h]
15.	</lixoftPLH>

Period		1					2				
FORM	Metrics	N	mean	CV(%)	GeoMean	GeoCV	N	mean	CV(%)	GeoMean	GeoCV
ref	AUC _{0-tlast} (h·mg·L ⁻¹)	9	154.8982	26.3033	150.3623	26.1424	9	147.6183	36.178	137.7505	43.7537
	C _{max} (mg·L ⁻¹)	9	9.82	11.8862	9.7571	12.1369	9	8.3244	32.019	7.8713	39.0231
test	AUC _{0-tlast} (h·mg·L ⁻¹)	9	133.3989	39.8426	122.3564	49.8543	9	158.3035	26.9691	153.1323	28.1126
	C _{max} (mg·L ⁻¹)	9	6.4689	28.8987	6.2384	29.1102	9	8.1089	22.1858	7.9035	25.4489

Renamings

Renamings setting can be provided to the placeholder to replace certain words or expressions that will appear in the table with a user-defined word or expression. We will explain the usage of renamings setting on a concrete example.

Let's take a look at the demo project `project_parallel_absBioavailability.pkx`. The project contains data from 18 subjects, of which 9 received an IV formulation of a certain drug and other 9 received an oral formulation. Let's say we want to use the Bioequivalence feature of PKanalix to calculate bioavailability of the oral formulation and generate a report. It would be suitable to put a bioequivalence confidence intervals table in the report, with certain renamings.

A generated placeholder for a BE confidence intervals table, with its default settings, is shown below, including the table it generates.

1.	<lixoftPLH>
2.	data:
3.	task: be
4.	table: confidenceIntervals
5.	metrics: [AdjustedMean, N, Difference, CIRawLower, CIRawUpper, Ratio, CILower, CIUpper, BE]
6.	parameters: [AUCINF_obs]
7.	display:
8.	units: true
9.	inlineUnits: true
10.	significantDigits: 4
11.	fitToContent: true
12.	</lixoftPLH>

ORAL

Parameter	ORAL	IV		Bioequivalence							
	AdjustedMean	N	AdjustedMean	N	Difference	CIRawLower	CIRawUpper	Ratio	CILower	CIUpper	BE
AUC _{0-inf} (h·ug·L ⁻¹)	223.9609	9	254.76	9	-0.1289	-0.2097	-0.04805	87.9105	81.0865	95.3089	yes

Let's say we would like to hide some of the columns and rename some of the expressions in the table. For example, it would be convient to rename "Bioequivalence" to "Absolute bioavailability", "IV" to "intravenous", "ORAL" to "oral formulation", and "CILower" to "90% CI (lower)". We can then modify the placeholder to look like the one on the left and it will generate the table on the right.

1.	<lixoftPLH>
2.	data:
3.	task: be
4.	table: confidenceIntervals
5.	metrics: [AdjustedMean, N, Ratio, CILower, CIUpper]
6.	parameters: [AUCINF_obs]
7.	renamings:
8.	IV: "intravenous"
9.	ORAL: "oral formulation"
10.	Bioequivalence: "Absolute bioavailability"
11.	CILower: "90% CI (lower)"
12.	CIUpper: "90% CI (upper)"
13.	AdjustedMean: "Adjusted mean"
14.	</lixoftPLH>

Oral capsule

Parameter	Oral capsule		IV solution		Absolute bioavailability		
	Adjusted mean	N	Adjusted mean	N	Ratio	90% CI (lower)	90% CI (upper)
AUC _{0-inf}	223.9609	9	254.76	9	87.9105	81.0865	95.3089

7.4. List of plot settings

This page describes all the settings that can be used in plot placeholders and are not editable by changing settings in the interface. All of them are optional (default values are described in the Default column).

Setting	Group	Default	Possible values	Description	Example
widthCm	-	16	positive numbers	Width of the plot in the report document in centimetres.	widthCm: 16
widthInches	-	6.3	positive numbers	Width of the plot in the report document in inches.	widthInches: 6.3
heightCm	-	11.1	positive numbers	Height of the plot in the report document in centimetres.	heightCm: 11.09
heightInches	-	4.37	positive numbers	Height of the plot in the report document in inches.	heightInches: 4.37

zoom	–	100	positive numbers	Larger values will increase the label font sizes and margins.	zoom: 80
caption	–	no caption	a phrase	Caption that will appear next to the plot.	caption: Individual fits
captionAbove	–	false	true, false	If true, caption will be positioned above the plot. If false, caption will be positioned below the plot.	captionAbove: true
legendPosition	settings	ne	n, ne, e, se, s, sw, w, nw	Position of the legend, based on abbreviations of compass points.	legendPosition: sw

For all other settings, select the settings of your plot in the GUI and then click the “[copy placeholder](#)” button to see the corresponding placeholder.

8. FAQ

This page summarizes the frequent questions about PKanalix.

- [Download, Installation, Run and Display issues](#)
- [Regulatory](#)
- [Running PKanalix](#)
- [Input data](#)
- [Settings \(options\)](#)
- [Settings \(output results\)](#)
- [Results](#)
- [Reporting](#)
- [Share your feedback](#)

Regulatory

- **Are NCA and CA analyses done with PKanalix accepted by the regulatory agencies like the FDA and EMA?** Yes.
- **How to cite Monolix, Simulx or PKanalix?** Please reference each application as here (adjusting to the version you used):

Monolix 2023R1, Lixoft SAS, a Simulations Plus company
 Simulx 2023R1, Lixoft SAS, a Simulations Plus company
 PKanalix 2023R1, Lixoft SAS, a Simulations Plus company

Running PKanalix

- **On what operating systems does PKanalix run?** PKanalix runs on Windows, Linux and MacOS platform. The requirements for each platform can be found on our [download](#) page.
- **Is it possible to run PKanalix in command line?** It is possible to run PKanalix from the R command line. This [complete R API](#) provides the full flexibility on running and modifying PKanalix projects with R functions.

Input data

- **Does PKanalix support sparse data?** No. If supporting sparse data is important for you, please [let us know](#) so that we take it into account for future versions.
- **Does PKanalix support drug-effect or PD models?** It is possible to use custom models in CA, and PD models are available via built-in model libraries. However, a time column is required in the dataset.
- **What type of data can PKanalix handle?** Extravascular, intravascular infusion, intravascular bolus for single-dose or steady-state plasma concentration and single-dose urine data can be used. See [here](#).
- **Can I give the concentration data and dosing data as separate files?** Yes, this is possible with the [data formatting](#) module.
- **Can I give the dosing information directly via the interface?** Yes, this is possible with the [data formatting](#) module.
- **Can I have BLQ data?** Yes, see [here](#).
- **Can I define units?** Yes, see [here](#).
- **Can I define variables such as “Sort” and “Carry”?** Yes, check [here](#).
- **It is possible to define several dosing routes within a single data set?** Yes, check the [ADMINISTRATION ID](#) column.
- **Can I use dose normalization to scale the dose by a factor?** Yes. In PKanalix 2023R1, the amount can now be specified in mass of the administered dose per body mass or body surface area.

Settings (options)

- **How do I indicate the type of model/data?** Extravascular versus intravascular is set in the [Settings window](#). Single versus steady-state and infusion versus bolus are imputed based on the data set column-types.
- **Can I exclude data with insufficient data?** The “[Acceptance criteria](#)” settings allow the user to define acceptance thresholds. The result tables can then be filtered in the interface according to these acceptance criteria.
- **Can I set the data as being sparse?** No, sparse data calculations are not supported. If supporting sparse data is important for you, please [let us know](#) so that we take it into account for future versions.
- **Which options are available for the lambda_z calculation?** Points to be included in the lambda_z calculation can be defined using the adjusted R2 criteria (called best fit in Winonlin), the R2 criteria, a time range or a number of points. In addition, the a specific range per individual can be defined, as well as points to include or exclude. Check the lambda_z options [here](#) and [here](#).
- **Can I define several partial areas?** Yes.
- **Can I save settings as template for future projects?** Not as such. However, the user can set the data set headers that should be recognized automatically in “Settings > Preferences”. Moreover, browsing a new dataset in a given project keeps the calculation settings that are

compatible with the new dataset.

- **Do you have a “slope selector”?** Yes, check here for information on the “[Check lambda_z](#)”.
- **Can I define a therapeutic response range?** No.
- **Can I set a user-defined weighting for the lambda_z?** No. But most common options are available as a setting.
- **Can I disable the calculation of lambda_z (curve stripping in Winnonlin)?** Yes, Lambda_z can be removed from the [list of parameters to compute](#).

Settings (output results)

- **Can I change the name of the parameters?** Yes, in the preferences of PKanalix (since 2023R1 version), you can rename NCA parameters as you wish. The output files will contain both your custom names and CDISC names.
- **Can I export the result table?** The result tables are automatically saved as text files in the result folder. In addition, result tables can be copy-pasted to Excel or Word using the copy button on the top right of the tables.
- **Can I choose which parameters to show in the result table and to export to the result folder file?** Yes, see [here](#).
- **Are the calculation rules the same as in Winonlin?** Yes.
- **Can I define myself the result folder?** By default, the result folder corresponds to the project name. However, you can define it by yourself. See [here](#) to see how to define it on the user interface.

Results

- **What result files are generated by PKanalix?** When running a task like NCA, CA or BE, [results](#) are automatically generated (one output file per table) in a result folder, and displayed in the RESULTS tab. Plots are also automatically generated and displayed in the PLOTS tab.
- **Can I replot the plots using another plotting software?** Yes, if you go to the menu Export and click on “Export charts data”, all the data needed to reproduce the plots are stored in text files.
- **When I open a project, my results are not loaded (message “Results have not been loaded due to an old inconsistent project”). Why?** When loading a project, PKanalix checks that the project (i.e all the information saved in the .pkx file) being loaded and the project that has been used to generate the results are the same. If not, the error message is shown. In that case the results will not be loaded because they are inconsistent with the loaded project.

Reporting

- **Can I generate a report?** Yes, starting from PKanalix 2023R1, [automated report generation](#) is possible.
- **How can I include plots and table from different PKanalix projects into the same Word document?** It is currently not possible to include plots and tables coming from different PKanalix runs into the same Word document in one step. You can add placeholders to an already generated report and use it as a template for another project, however we are aware that this is not a very convenient solution. If generating a single report for several projects is important for you, please [let us know](#) so that we take it into account for future versions.
- **The report I generated via command line or via R do not contain plots.** This is a known limitation of the current implementation of the reporting feature. If generating reports with plots from R or from the command line is important for you, please [let us know](#) so that we take it into account for future versions.

8.1. List of Known Bugs

The list of known bugs in every version of MonolixSuite can be found in the release notes of the next version.

- [Release Notes for MonolixSuite2023R1](#)
- [Release Notes for MonolixSuite2021R2](#)
- [Release Notes for MonolixSuite2021R1](#)
- [Release Notes for MonolixSuite2020R1](#)

We list below the known bugs in PKanalix 2023R1 (they will be fixed in the next version). The 2023R1 version of PKanalix can be downloaded [here](#).
Installation

- The lixoft folder inside temp is not deleted after installation.

Documentation

- In the Documentation tab, the tutorials section can be empty (due to issues with the host server).

Plots & Reporting

- If the axes in bivariate data viewer have been modified (to select which observation goes on which axis), default report generation can get frozen.
- [Observed data] numbers in Information box are based on all ids from the dataset, and not all ids in the Observed data plot (ie if there are ids without data in this plot they are still counted).
- Windows only: if a project file path contains more than 168 characters, generating a report might fail with a following error:
“System.IO.DirectoryNotFoundException: Path not found ‘pathToProject\PKanalix\ChartsFigures\covariatemodeldiagnosis_nca_0_a3e5281c97564a1fb22c63feef4f4d4c_0.svg”

Data set

- For specific projects, the warning message about ignored columns in external treatment file of data formatting appears each time we save the project.
- If an ID has only censored data, when navigating through IDs in the selection panel of stratify and reaching this ID, PKanalix freezes then shuts down.

CA

- After exporting a project with a scaling factor in the units to Monolix, the scaling is wrongly applied on the CA predictions in the Check init. tab.

Connectors

- Mapping arguments are not taken into account in `newProject()`. You need to use `setMapping()` instead.
- The documentation of `setNCASettings` is not correct. It is now **corrected online** (and will be fixed in the package for the next version).
Correction: parameters calculated in case of plasma data if partial AUC calculation intervals are provided through the *partialAucTime*: "AUC_int", "AUC_int_D", "CAVG_int".
- Filtering NCA results statistics by acceptance criteria is incorrect. It selects lines for which the value of Flag = 0, but should take values when Flag=1 as in the GUI.

8.2. Share your feedback

Help us improve PKanalix!

We value your opinion. We would love to hear from you to align our vision and concepts with your needs and expectations. Let's make sure you get the most of our software.

NB: This is not a support form. Commercial users can contact us via the support email.

Name *

Company/Institute *

Work Email *

Is there something we can improve? *

Upload

Express Consent *

☐ In compliance with applicable laws and regulations, by checking this box and sending this contact form, you are giving your express consent to Lixoft to process, and to transfer to Lixoft's affiliate located in the United States (Simulations Plus, Inc.) to process, your personal data for contacting you for only the business purpose included in the consent form. For the exercise of your rights or any other information related to processing of your personal data by us, please refer to our Privacy Policy at <https://www.simulations-plus.com/privacy-policy/>.

☐ I agree to receive the Lixoft newsletter.

[Edit this form](#)

9. Case studies

Different cases studies are presented here to show how to use Monolix and the MonolixSuite for Modeling and Simulation.

Warfarin case study

This video case study shows a simple PK modeling workflow in Monolix2018, with the example of warfarin. It explains the main features and algorithms of Monolix, that guide the iterative process of model building: from validating the structural model to adjusting the statistical model step-by-step.

It includes picking a model from the libraries, choosing initial estimates with the help of population predictions, estimating parameters and uncertainty, and diagnosing the model with interactive plots and statistical tests.

Tobramycin case study

This case study presents the modeling of the **tobramycin pharmacokinetics**, and the determination of a priori dosing regimens in patients with various degrees of renal function impairment. It takes advantage of the integrated use of [Datxplorer](#) for data visualization, [Mlxplore](#) for model exploration, [Monolix](#) for parameter estimation and [Simulx](#) for simulations and best dosing regimen determination.

The case study is presented in 5 sequential parts, that we recommend to read in order: [Part 1: Introduction](#), [Part 2: Data visualization with Datxplorer](#), [Part 3: Model development with Monolix](#), [Part 4: Model exploration with Mlxplore](#), and [Part 5: Dosing regimen simulations with Simulx](#).

Remifentanil case study

Remifentanil is an opioid analgesic drug with a rapid onset and rapid recovery time. It is used for sedation as well as combined with other medications for use in general anesthesia. It is given in adults via continuous IV infusion, with doses that may be adjusted to age and weight of patients.

This case-study shows how to use Monolix to build a **population pharmacokinetic model for remifentanil** in order to determine the influence of subject covariates on the individual parameters.

[Link to Remifentanil case study](#)

Longitudinal Model-Based Meta-Analysis (MBMA) with Monolix Suite

Longitudinal model-based meta-analysis (MBMA) models can be implemented using the MonolixSuite. These models use study-level aggregate data from the literature and can usually be formulated as non-linear mixed-effects models in which the inter-arm variability and residual error are weighted by the number of individuals per arm. We exemplify the model development and analysis workflow of MBMA models in [Monolix](#) using a real data set for rheumatoid arthritis, following publication from [Demin et al \(2012\)](#). In the case study, the efficacy of a drug in development (Canakinumab) is compared to the efficacy of two drugs already on the market (Adalimumab and Abatacept). Simulations using [Simulx](#) were used for decision support to see if the new drug has a chance to be a better drug.

[Link to MBMA case study](#)

Analysis of time-to-event data

Within the MonolixSuite, the mlxtran language allows to describe and model time-to-event data using a parametric approach. [This page](#) provides an introduction on time-to-event data, the different ways to model this kind of data, and typical parametric models. A library of common TTE models is also provided.

- [Introduction to time-to-event data](#)

Two modeling and simulation workflows illustrate this approach, using two TTE data sets:

- [Case study on veteran lung cancer data set](#)
- [Case study on NCCTG lung cancer data set <](#)

Veralipride case study

Multiple peaking in plasma concentration-time curves is not uncommon, and can create difficulties in the determination of pharmacokinetic parameters.

For example, double peaks have been observed in plasma concentrations of veralipride after oral absorption. While multiple peaking can be explained by different physiological processes, in this case site-specific absorption has been suggested to be the major mechanism. In this webinar we explore this hypothesis by setting up a population PK modeling workflow with the MonolixSuite 2018.

The step-by-step workflow includes visualizing the data set to characterize the double peaks, setting up and estimating a double absorption model, assessing the uncertainty of the parameter estimates to avoid over-parameterization, and simulations of the model.

[Link to veralipride case study](#)